

Coloplast Product Guide

Ostomy, Continence, Wound & Skin Care

Prices effective January 1st, 2013





Making life easier for people with intimate healthcare needs.

Our values are:

- Closeness... to better understand
- Passion... to make a difference
- Respect and responsibility... to guide us

Closeness to all customers makes this possible. We listen to better understand needs, and respond by finding new ways to do things better together. We lead the way by bringing the best ideas first and fast to market in the form of medical devices and service solutions.

Deeply private and personal medical conditions are our focus. Our passion to make a real difference to people's lives is what drives and unites us. Our culture supports high ambitions, and releases the full potential of our own people to achieve them.

We welcome the broader responsibility that comes with leadership - a responsibility to the environment, to society, to our shareholders, and to act with integrity in all we do. Our values define the way we think and act, both as individuals and as a company.

Ostomy Care Table of Contents

Product Category	Description	Page	
Ostomy Care		1	
SenSura® Two-Piece Ostomy System Stoma Sizing Chart		3	
	Skin Barriers	SenSura® Click Skin Barrier with Flange	4
		SenSura® Click Xpro Skin Barrier with Flange	5
		SenSura® Flex Skin Barrier with Flange	9
		SenSura® Flex Xpro Barrier with Flange	10
SenSura® Two-Piece Systems	Pouches	SenSura® Click Easiclose™ WIDE-Outlet Drainable Pouch	6
		SenSura® Click Closed Pouch	7
		SenSura® Click Urostomy Pouch	8
		SenSura® Flex Easiclose™ WIDE-Outlet Drainable Pouch	11
		SenSura® Flex Closed Pouch	11
		SenSura® Flex Urostomy Pouch	12
SenSura® One-Piece Systems		SenSura® Drainable Pouch	12
		SenSura® Xpro Drainable Pouch	13
		SenSura® Mio Closed Pouch	15
		SenSura® Closed Pouch	16
		SenSura® Urostomy Pouch	17
Assura® Two Piece Ostomy System Stoma Sizing Chart		18	
	Skin Barriers	Assura® Skin Barrier With Flange	19
		Assura® Extra-Extended Wear Skin Barrier With Flange	20
Assura® Two-Piece Systems	Pouches	Assura® Easiclose™ WIDE-Outlet Drainable Pouch	21
		Assura® Drainable Pouch	23
		Assura® Multi-chamber Urostomy Pouch	23
		Assura® Closed Pouch	24
Easiflex® Two-Piece Systems	Skin Barriers	Easiflex® Skin Barrier With Flange	25
	Pouches	Easiflex® Easiclose™ Drainable Pouch	26
		Easiflex® Closed Pouch	27
		Easiflex® Multi-chamber Urostomy Pouch	27
Assura® One-Piece Systems		Assura® Extra-Extended Wear Drainable Pouch	28
		Assura® Extra-Extended Wear Urostomy Pouch	28
		Assura® Extra-Extended Wear Multi-chamber Urostomy Pouch	29
		Assura® Easiclose™ Drainable Pouch	29
		Assura® Drainable Pouch	31
		Assura® Urostomy Pouch	32
		Assura® Closed Pouch	33

Ostomy Care Table of Contents

Product Category	Description	Page	
Pediatric Two-Piece Ostomy Systems	Skin Barrier	Assura® Pediatric Skin Barrier With Flange	34
		Easiflex® Pediatric Skin Barrier With Flange	35
	Pouch	Assura® Pediatric Pouch	34
		Easiflex® Easiclose™ Pediatric Drainable Pouch	35
Pediatric One-Piece Ostomy System	Assura® Pediatric Pouch	34	
Ostomy Accessories	Brava Protective Sheets	36	
	Brava Mouldable Rings	36	
	Brava Skin Barrier Rings	36	
	Brava Skin Barrier Spray and Wipe	36	
	BravaAdhesive Remover Spray and Wipe	36	
	Brava Ostomy Strip Paste	36	
	Brava Elastic Tape	37	
	Brava Ostomy Powder	37	
	Brava Belt	37	
	Brava Lubricating Deodorant	37	
	Ostomy Paste	38	
	SKin Barrier Rings	38	
	Coloplast Prep® Protective Liquid Skin Barrier	38	
	Coloplast Filtrodor® Pouch Filter	38	
	Coloplast Pouch Clamp	38	
	Coloplast Scissors	38	
Coloplast Stoma Caps	39		
Irrigation Products	Economy version	40	
	Stoma Cone	40	
	Irrigation Faceplate	40	
	Irrigation Sleeves	40	
SenSura® Hospital Assortment	SenSura® One-piece Drainable Pouches	41	
	SenSura® Post Op	41	
	SenSura® High Output	41	
Post-Op and Wound Management Systems	Two-piece High Output Drainable System	42	
	One piece Post Op Drainable Pouch	42	
	Urostomy Night Bag	42	
	Coloplast Fistula & Wound Management System	43	
	Flexible Lids Coloplast Fistula & Wound Management System	43	
	Accessories Coloplast Fistula & Wound Management System	43	
Coloplast Support Service		44	

Continence Care Table of Contents

Product Category	Description	Page
Continence care		45
Intermittent Catheters	SpeediCath™ Intermittent Catheters	46
	SpeediCath™ Compact Intermittent Catheters	47
	SpeediCath™ Complete Intermittent Catheter Sets	46
	Self-Cath® Intermittent Catheters	48
	Self-Cath® Plus Intermittent Catheters	49
	Self-Cath® Closed Systems	50
Male External Catheters	Conveen® Optima Self Sealing Male External Catheter	51
	Conveen® Security + Self Sealing Male External Catheter	51
	Conveen® Security + Male External Catheter and Liner	51
	Clear Advantage® Self Sealing Male External Catheter	52
	Freedom Clear™ Self Sealing Male External Catheter	52
	Freedom Clear™ 'LS' Long Sheath Self Sealing Male External Catheter	52
	Freedom Clear™ 'SS' Sport Sheath Self Sealing Male External Catheter	53
	Conveen® Self Sealing Male External Catheter	53
	Conveen® Male External Catheter and Liner	53
	Active Cath®	54
	Freedom Cath®	54
Urine Bags	Conveen® Security + Leg Bags	55
	Conveen® Security+ Contour Leg Bags	55
	Conveen® Security+ Extra Large Leg Bag/Bedside Drainage Bags	55
	Freedom® Leg Bags	56
	Freedom® Accessories	56
	Moveen® Bedside Night Bags	56
Continence Care Accessories	Moveen® Bag Hanger	57
	Conveen® Security+ Fabric Leg Bag Straps	57
Peristeen™ Anal Plug		57
Peristeen™ Anal Irrigation		58

Wound and Skin Care Table of Contents

Product Category	Description	Page
Wound Care		59
Foam Dressings	Biatain® Foam Dressing	61
	Biatain® Soft-Hold Foam Dressing	62
	Biatain® Silicone Foam Dressing	63
	Biatain® Silicone Lite Foam Dressing	64
	Biatain® Ag Foam Antimicrobial Dressing with Silver	65
	Biatain® Silicone Ag Foam Antimicrobial Dressing with Silver	66
	Biatain® Ibu Foam Dressing with Ibuprofen	67
Calcium Alginate Dressings	Biatain® Alginate Ag Calcium Alginate Antimicrobial Dressing with Silver	68
	Biatain® Alginate Calcium Alginate Dressing (SeaSorb® Soft)	68
Hydrocolloid Dressings	Comfeel® Plus Ulcer Dressing	69
	Comfeel® Plus Sacral Dressing	69
	Comfeel® Plus Transparent Dressing	69
Hydrogels	Purilon® Gel	70
	Woun'Dres® Collagen Hydrogel	70
Hydrophilic Paste	Triad™ Hydrophilic Wound Dressing Paste	70
Wound Cleanser	Sea-Clens® Wound Cleanser	70
Cutting Guide	Application and Cutting Guide	71
Skin Care		73
Bathing and Hand Wash Products	Gentle Rain® Extra Mild Cleanser	75
	Gentle Rain® Antibacterial Cleanser	75
Antiseptic Hand Cleansers	Isagel® Hand Sanitizing Gel	75
No-Rinse Full Body and Incontinence Cleanser	Bedside-Care® Sensitive Skin (Bedside-Care® Unscented Foam)	76
	Bedside-Care® Scented Foam	76
	Sproam® Cleanser	76
	Bedside-Care® Cleanser	76
	Peri Wash® II Cleanser	76
Moisturizers	Sween® 24 Cream	77
	Atrac-Tain® Cream	77
	Sween® Unscented Lotion	77
	Sween Cream®	77
Incontinence Cleanup and Protection (Barriers)	Baza Cleanse & Protect® Lotion	78
	Baza® Protect II Cream	78
	Critic-Aid® Clear Ointment	78
Skin Fold Management	InterDry® Ag Textile	78
Other	Sween® Body Powder	79
	Hex-On®	79
Dispensers	Wall Dispenser for Bag	79
	Wall Dispenser Drip Tray	79
	Dispenser Bracket	79
Discontinued Products		80

Ostomy Care



Coloplast Ostomy Care



Body profiles vary. Find the solution that fits yours with the online BodyCheck.

Bodies not only come in all shapes and sizes, they also vary over time.

That's why it's important to check if your existing stoma solution is still an ideal fit, and why we developed online BodyCheck. It's the quick, simple and confidential way to find the solution that's best for you.

Identify your body profile with the online BodyCheck and discover the way to optimal security and skin care.

Find the solution that's absolutely right at www.bodycheck.coloplast.ca









All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

SenSura® Two-Piece Ostomy System - Stoma Sizing Chart

Coloplast ostomy products are packaged in colour-coded boxes for easy product identification.

Two-piece products are packaged in boxes with a coloured stripe on the label. This colour code makes it easy to match the coupling sizes of flanges and pouches.

SenSura Click Two-Piece		
Colour code	Skin barrier size Diameter or "Foot print"	Coupling size
Green 	110 x 90mm - 4 3/8" x 3.5"	40mm - 1 9/16"
Red 	120 x 100mm - 4 3/4" x 3 7/8"	50mm - 2"
Blue 	130 x 110mm - 5 1/16" x 4 3/8"	60mm - 2 3/8"
Yellow 	140 x 120mm - 5 1/2" x 4 3/4"	70mm - 2 3/4"
SenSura Flex Two-Piece		
Colour code	Skin barrier size Diameter or "Foot print"	Coupling size
Green 	120 x 100mm - 4 3/4" x 3 7/8"	40mm - 1 9/16"
Red 	130 x 110mm - 5 1/16" x 4 3/8"	50mm - 2"
Yellow 	150 x 130mm - 5 7/8" x 5 1/16"	70mm - 2 3/4"
Silver 	170 x 150mm - 6 6/8" x 5 7/8"	90mm - 3 1/2"

Quick Metric to Imperial Conversion Charts for SenSura


Cut-to-fit	
Metric	Imperial
10-33mm	3/8" - 1 1/4"
10-48mm	3/8" - 2"
10-68mm	3/8" - 2 3/4"
10-88mm	3/8" - 3 1/2"
10-35mm	3/8" - 1 3/8"
10-45mm	3/8" - 1 3/4"
10-55mm	3/8" - 2 1/8"
10-65mm	3/8" - 2 1/2"
10-66mm	3/8" - 2 5/8"
10-76mm	3/8" - 3"
15-23mm	5/8" - 7/8"
15-33mm	5/8" - 1 1/4"
15-43mm	5/8" - 1 3/4"
15-53mm	5/8" - 2 1/16"
15-56mm	5/8" - 2 1/8"
40-55mm	1 9/16" - 2 1/8"

Pre-Cut	
Metric	Imperial
20mm	3/4"
21mm	7/8"
25mm	1"
28mm	1 1/8"
30mm	1 1/8"
31mm	1 1/4"
35mm	1 3/8"
40mm	1 9/16"

Coloplast Ostomy Care

SenSura® Click Two-Piece System

SenSura Click is the two-piece standard in ostomy care, offering a unique audible click-lock security ring. The patented click-lock technology provides the ultimate in security, convenience and simplicity with just a click to lock and a click to unlock. When the click system is locked, the pouch is securely fastened to the skin barrier.



SenSura Click combines the revolutionary double-layer adhesive with the added security of the click-lock security ring.

SenSura Click Two-Piece Ostomy Systems

Non-convex=Flat Convex Light= Shallow Convexity



SenSura® Click Skin Barrier With Flange

SenSura offers the flexible skin barrier with patented double layer technology for secure, safe, long-lasting wear. With belt loops.

	Convexity	Coupling Size	Stoma Size	Product Code	Units	Suggested Retail Price
Cut-to-Fit	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	10011	5	\$45.57
		50mm - 2" ●	10-45mm - 3/8"- 1 3/4"	10021		
		60mm - 2 3/8" ●	10-55mm - 3/8"- 2 1/8"	10031		
		70mm - 2 3/4" ●	10-65mm - 3/8"- 2 1/2"	10041		
	Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	11011		\$61.99
		50mm - 2" ●	15-33mm - 5/8"- 1 1/4"	11021		
		60mm - 2 3/8" ●	15-43mm - 5/8"- 1 3/4"	11031		
		70mm - 2 3/4" ●	15-53mm - 5/8"-2 1/16"	11041		
Pre-Cut	Non-Convex	40mm - 1 9/16" ●	20mm - 3/4"	10012	5	\$47.39
			25mm - 1"	10013		
		50mm - 2" ●	30mm - 1 1/8"	10022		
			35mm - 1 3/8"	10023		
	60mm - 2 3/8" ●	40mm - 1 9/16"	10032	\$64.44		
	Convex Light	40mm - 1 9/16" ●	21mm - 7/8"			11012
			25mm - 1"			11022
		50mm - 2" ●	28mm - 1 1/8"			11023
			31mm - 1 1/4"			11024
	60mm - 2 3/8" ●	35mm - 1 3/8"	11032			

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® Click Xpro Two-Piece Ostomy Systems

Non-convex= Flat
 Convex Light= Shallow Convexity



SenSura Click Xpro Skin Barrier with Flange

SenSura Click Xpro has the same double layer adhesive as SenSura Click, however two extra ingredients have been added to protect against aggressive output from an ileostomy or urostomy.
 With belt loops.

	Convexity	Coupling Size	Stoma Size	Product Code	Units	Suggested Retail Price
Cut-to-Fit	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	10015	5	\$57.13
		50mm - 2" ●	10-45mm - 3/8"-1 3/4"	10025		
		60mm - 2 3/8" ●	10-55mm - 3/8"-2 1/8"	10035		
		70mm - 2 3/4" ●	10-65mm - 3/8"-2 1/2"	10045		
	Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	11015		\$82.85
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	11025		
		60mm - 2 3/8" ●	15-43mm - 5/8"-1 3/4"	11035		
		70mm - 2 3/4" ●	15-53mm - 5/8"-2 1/16"	11045		
Pre-Cut	Non-Convex	40mm - 1 9/16" ●	20mm - 3/4"	10016	5	\$59.41
			25mm - 1"	10017		
		50mm - 2" ●	30mm - 1 1/8"	10026		
			35mm - 1 3/8"	10027		
			60mm - 2 3/8" ●	40mm - 1 9/16"		
	Convex Light	40mm - 1 9/16" ●	21mm - 7/8"	11016		\$86.17
			25mm - 1"	11026		
		50mm - 2" ●	28mm - 1 1/8"	11027		
			31mm - 1 1/4"	11028		
			60mm - 2 3/8" ●	35mm - 1 3/8"		

Coloplast Ostomy Care

SenSura® Click Two-Piece Ostomy Systems

SenSura Click Easiclose™ WIDE-Outlet Drainable Pouch

An extremely hygienic and secure pouch for people with colostomies and ileostomies that features an easy-to-use patented WIDE-Outlet Velcro closure system. No clamps or clips required. Soft, non-woven comfort backing. The one and only "Click" lock secure coupling system. Can be rotated/repositioned for comfort without removing pouch. Low-profile for maximum discretion. Developed in cooperation with enterostomal therapy nurses worldwide.

SenSura Click Easiclose WIDE-Outlet Drainable Pouches must be used with the corresponding colour coded SenSura Click double layer skin barrier.



Midi Drainable Opaque Pouch – 27cm (10 1/2")





Mid-size, opaque, drainable pouch with unique Easiclose WIDE-Outlet. Only available with filter.

Coupling Size	Product Code	Units
40mm - 1 9/16" 	11114	10
50mm - 2" 	11115	



Maxi Drainable Opaque Pouch - 30cm (12")

Large-size, opaque, drainable pouch with unique Easiclose WIDE-Outlet. Available with or without filter.

Coupling Size	Product Code		Units
	With filter	Without filter	
40mm - 1 9/16" 	11124	N/A	10
50mm - 2" 	11125	11135	
60mm - 2 3/8" 	11126	11136	
70mm - 2 3/4" 	11127	N/A	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com





Coloplast Ostomy Care

SenSura® Click Two-Piece Ostomy Systems



Maxi Drainable Transparent Pouch - 30cm (12")

Large-size, transparent, drainable pouch with unique Easyclose WIDE-Outlet feature. Available with or without filter.

Coupling Size	Product Code		Units
	With filter	Without filter	
40mm - 1 9/16" 	11194	N/A	10
50mm - 2" 	11195	11185	
60mm - 2 3/8" 	11196	11186	
70mm - 2 3/4" 	11197	N/A	

SenSura Click Closed Pouch








With soft, non-woven comfort backing. The one and only "Click" lock coupling system. Can be rotated and repositioned for comfort without removing pouch. Low-profile for maximum discretion. Developed in cooperation with enterostomal therapy nurses worldwide.

SenSura Click Closed Pouches must be used with the corresponding colour coded SenSura Click Skin Barrier.



Closed Opaque Pouch

Only available with filter.

Pouch Size	Coupling Size	Product Code	Units
Midi 18cm (7")	40mm - 1 9/16" 	10154	30
	50mm - 2" 	10155	
	60mm - 2 3/8" 	10156	
Maxi 22cm (8 1/2")	40mm - 1 9/16" 	10164	
	50mm - 2" 	10165	
	60mm - 2 3/8" 	10166	
	70mm - 2 3/4" 	10167	

Coloplast Ostomy Care

SenSura® Click Two-Piece Ostomy Systems

Opaque Urostomy Pouch



The unique multichamber pouch introduced in the Assura product line has been enhanced to ensure the ultimate comfort and discretion. Beige, soft, non-woven, water repellent backing and anti-reflux valve. Audible "Click" locking system. Can be rotated/ repositioned for comfort without removing pouch. Secure, tuck away outlet and simple connector for night drainage bag attachment.

Pouch Size	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	40mm - 1 9/16" ●	11841	10
	50mm - 2" ●	11842	
Maxi 26cm (10 1/4")	40mm - 1 9/16" ●	11844	
	50mm - 2" ●	11845	
	60mm - 2 3/8" ●	11846	

Transparent Urostomy Pouch



Pouch Size	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	40mm - 1 9/16" ●	11851	10
	50mm - 2" ●	11852	
Maxi 26cm (10 1/4")	40mm - 1 9/16" ●	11854	
	50mm - 2" ●	11855	
	60mm - 2 3/8" ●	11856	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® Flex Two-Piece Ostomy System with Adhesive Coupling

Non-convex= Flat
Convex Light= Shallow Convexity



SenSura Flex Skin Barrier With Flange

The extremely low profile adhesive coupling system provides ultimate flexibility and freedom of movement. Patented security channel prevents leakage. Coloured tab and ample coupling surface allow for easy application and removal of pouch. SenSura Flex has the SenSura double layer adhesive that makes it the ideal solution for post-operative situations.

Convexity has belt loops.

	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm - 1 3/8" ●	10-33mm - 3/8"-1 1/4"	10101	10
		50mm - 2" ●	10-48mm - 3/8"-2"	10102	
		70mm - 2 3/4" ●	10-68mm - 3/8"-2 3/4"	10103	
		90mm - 3 1/2" ●	10-88mm - 3/8"-3 1/2"	10108	5
	Convex Light	35mm - 1 3/8" ●	15-23mm - 5/8"-7/8"	11301	5
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	11302	
			15-43mm - 5/8"-1 3/4"	11303	
70mm - 2 3/4" ●	15-56mm - 5/8"-2 1/8"	11307			
Pre-Cut	Non-Convex	35mm - 1 3/8" ●	21mm - 7/8"	10111	10
			28mm - 1 1/8"	10113	
		50mm - 2" ●	25mm - 1"	10114	
			31mm - 1 1/4"	10116	
			35mm - 1 3/8"	10117	
	Convex Light	35mm - 1 3/8" ●	21mm - 7/8"	11312	5
			25mm - 1"	11313	
		50mm - 2" ●	28mm - 1 1/8"	11314	
			31mm - 1 1/4"	11315	
			35mm - 1 3/8"	11316	

Coloplast Ostomy Care

SenSura® Flex Xpro Two-Piece Ostomy System with Adhesive Coupling

Non-convex= Flat
Convex Light=Shallow Convexity

SenSura Flex Xpro Skin Barrier with Flange



SenSura Flex Xpro has the same double layer adhesive as SenSura Flex, however two extra ingredients have been added to protect against aggressive output from an ileostomy or urostomy. The extremely low profile adhesive coupling system provides the ultimate in flexibility and freedom of movement. Patented security channel prevents leakage. Colour tab and ample coupling surface allow for easy application and removal of pouch.

Convexity has belt loops.

	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	35mm - 1 3/8" ●	10-33mm - 3/8"-1 1/4"	10104	10
		50mm - 2" ●	10-48mm - 3/8"-2"	10105	
		70mm - 2 3/4" ●	10-68mm - 3/8"-2 3/4"	10106	
	Convex Light	35mm - 1 3/8" ●	15-23mm - 5/8"-7/8"	11304	5
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	11305	
			15-43mm - 5/8"-1 3/4"	11306	
70mm - 2 3/4" ●	15-56mm - 5/8"-2 1/8"	11308			



SenSura – Double-layer adhesive

SenSura features a unique double-layer adhesive, which sets new standards within ostomy care. It has all the features required for an effective and secure adhesive.

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® Flex Two-Piece Ostomy Systems with Adhesive Coupling



Drainable Pouch with Easyclose WIDE-Outlet

An extremely hygienic and secure pouch for people with colostomies and ileostomies that features an easy-to-use patented WIDE-Outlet Velcro closure system. No clamps or clips required. Soft, non-woven comfort backing. Low-profile for maximum discretion. Developed in cooperation with enterostomal therapy nurses worldwide.

Pouch Size	Pouch Color	Coupling Size	Product Code		Units
			With filter	Without filter	
Midi 27cm (10 1/2")	Opaque	50mm - 2" ●	11505	N/A	20
Maxi 30cm (12")	Opaque	35mm - 1 3/8" ●	11514	N/A	
		50mm - 2" ●	11515	11615	
		70mm - 2 3/4" ●	11518	NA	
	Transparent	50mm - 2" ●	N/A	11612	
		90mm - 3 1/2" ●	N/A	11690	



Closed Opaque Pouch

SenSura Flex pouches must be used with the corresponding colour coded SenSura Flex Skin Barrier.

Pouch Size	Coupling Size	Product Code	Units
		With filter	
Midi 18cm (7")	35mm - 1 3/8" ●	10904	30
	50mm - 2" ●	10905	
Maxi 22cm (8 1/2")	35mm - 1 3/8" ●	10914	
	50mm - 2" ●	10915	
	70mm - 2 3/4" ●	10918	

Coloplast Ostomy Care

SenSura® Flex Two-Piece Ostomy Systems with Adhesive Coupling



Urostomy Pouch

SenSura Flex Urostomy Pouches feature the same multichamber design found in Easiflex Multi-Chamber Urostomy Pouches. The innovative design provides a more discreet and quiet urostomy pouch. SenSura Flex Urostomy Pouches feature soft, moisture absorbent backing and anti-reflux valve

Pouch Size	Pouch Type	Coupling Size	Product Code	Units
Midi 24cm (9 1/2")	Opaque	35mm - 1 3/8" ●	11821	20
	Opaque	50mm - 2" ●	11822	
Maxi 26cm (10 1/4")	Opaque	35mm - 1 3/8" ●	11824	
	Opaque	50mm - 2" ●	11825	
	Transparent	35mm - 1 3/8" ●	11827	
	Transparent	50mm - 2" ●	11828	

SenSura® One-Piece Ostomy Systems

SenSura offers a unique double layered adhesive. The lower layer provides absorption and moisture balance while the upper layer is erosion resistant. The new adhesive formulation also offers optional adhesion, flexibility and ease of removal. The SenSura one-piece pouch features an improved WIDE Outlet.



Maxi Drainable Opaque Pouch – 30cm (12")

Large, opaque, drainable pouch with soft fabric on both sides and unique Easiclose WIDE-Outlet. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-76mm - 3/8"-3"	15531	10
	Convex Light	15-33mm - 5/8"-1 1/4"	15607	
		15-43mm - 5/8"-1 3/4"	15608	
Pre-Cut	Non-Convex	25mm - 1"	15532	
		30mm - 1 1/8"	15533	
		35mm - 1 3/8"	15534	
	Convex Light	25mm - 1"	15633	
		28mm - 1 1/8"	15634	
		31mm - 1 1/4"	15635	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® One-Piece Ostomy Systems



Maxi Drainable Transparent Pouch – 30cm (12")

Large, transparent, drainable pouch with soft fabric on back and unique Easyclose WIDE-Outlet. Available with and without filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-76mm - 3/8"-3"	15521	10
		10-76mm - 3/8"-3"	15821*	
	Convex Light	15-33mm - 5/8"-1 1/4"	15605	
		15-43mm - 5/8"-1 3/4"	15606	
Pre-Cut	Non-Convex	25mm - 1"	15522	
		30mm - 1 1/8"	15523	
		35mm - 1 3/8"	15524	

* Without filter



Xpro Maxi Drainable Transparent Pouch - 30cm (12")

Large, transparent, drainable pouch with soft fabric on back and unique Easyclose WIDE-Outlet. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-76mm - 3/8"-3"	15981	20
	Convex Light	15-33mm - 5/8"-1 1/4"	15695	10
		15-43mm - 5/8"-1 3/4"	15696	

Coloplast Ostomy Care



Fits individual body shapes
Feels secure

SenSura® Mio responds to natural body movements

Introducing SenSura Mio, a new colostomy appliance designed with an adhesive that is formulated to enhance elasticity and fitting performance. The elastic adhesive is based on BodyFit™ Technology. It is pliable to individual body contours and maintains a secure fit as it responds to natural body movements.

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care


SenSura® One-Piece Ostomy Systems

SenSura® Mio

SenSura Mio features the new inspection window covering. It combines the discretion of an opaque pouch with the functionality of a transparent pouch.

SenSura Mio is a new colostomy appliance designed with individual body shapes in mind. It features a new elastic adhesive:

- pliable to individual body contours, creating a secure fit
- responsive to natural body movements, maintaining a secure fit



SenSura Mio has a new elastic adhesive.



The elastic adhesive is responsive to natural body movements.



SenSura Mio

Opaque, closed pouch with inspection window.
Only available with filter.



Pouch Size	Stoma Size		Product Code	Units
Mini 14cm (5 1/2")	Cut-to-Fit	15-45mm -5/8"-1 3/4"	10811	30
Midi 18cm (7')	Cut-to-Fit	15-45mm -5/8"-1 3/4"	10821	
	Pre-Cut	25mm - 1"	10822	
		30mm - 1 1/8"	10823	
		35mm - 1 3/8"	10824	
Maxi 22cm (8 1/2')	Cut-to-Fit	15-45mm -5/8"-1 3/4"	10831	
		40-55mm -1 9/16" - 2 1/8"	10835	
	Pre-Cut	30mm - 1 1/8"	10833	
		35mm - 1 3/8"	10834	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® One-Piece Ostomy Systems



Midi Closed Opaque Pouch – 18cm (7")

Medium, opaque, closed pouch with soft fabric on both sides. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-66mm - 3/8"-2 5/8"	15450	30
	Convex Light	15-33mm - 5/8"-1 1/4"	15403	10
		15-43mm - 5/8"-1 3/4"	15404	
Pre Cut	Convex Light	25mm - 1"	15413	
		28mm - 1 1/8"	15414	
		31mm - 1 1/4"	15415	

Maxi Closed Opaque Pouch – 22cm (8 1/2")

Large, opaque, closed pouch with soft fabric on both sides. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-76mm - 3/8"-3"	15480	30
	Convex Light	15-33mm - 5/8"-1 1/4"	15407	10
		15-43mm -5/8"-1 3/4"	15408	
Pre-Cut	Non-Convex	25mm - 1"	15482	30
		30mm - 1 1/8"	15483	
		35mm - 1 3/8"	15484	

Midi Closed Transparent Pouch – 18cm (7")

Medium, transparent, closed pouch with soft fabric on back. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-66mm - 3/8"-2 5/8"	15441	30

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura® One-Piece Ostomy Systems



Maxi Closed Transparent Pouch – 22cm (8 1/2")

Large, transparent, closed pouch with soft fabric on back. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-76mm - 3/8"-3"	15470	30
	Convex Light	15-33mm - 5/8"-1 1/4"	15405	10
		15-43mm - 5/8"-1 3/4"	15406	
Pre-Cut	Non-Convex	25mm - 1"	15472	30
		30mm - 1 1/8"	15473	
		35mm - 1 3/8"	15474	



Urostomy Pouch Midi - 24 cm (9 1/2")

SenSura Urostomy Pouches have the same double layer adhesive as SenSura and come with soft backing, odour-proof film, an anti-reflux valve and simple closure with adaptor for night drainage bag.

	Convexity	Pouch Colour	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	Transparent	10-66mm - 3/8"-2 5/8"	11803	10








Urostomy Pouch Maxi - 26 cm (10 1/4")

	Convexity	Pouch Colour	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	Opaque	10-76mm - 3/8"-3"	11802	10
		Transparent	10-76mm - 3/8"-3"	11804	
	Convex Light	Opaque	15-33mm - 5/8"-1 1/4"	11811	
			15-43mm - 5/8"-1 3/4"	11812	
		Transparent	15-33mm - 5/8"-1 1/4"	11814	
			15-43mm - 5/8"-1 3/4"	11815	

Coloplast Ostomy Care

Assura® and Easiflex® Two-Piece System Stoma Sizing Chart

Coloplast ostomy products are packaged in colour-coded boxes for easy product identification. Two-piece products are packaged in boxes with a coloured stripe on the label. This colour code makes it easy to match the coupling sizes of flanges and pouches.

Assura® Two-Piece skin barriers with flange		
Colour code	Skin barrier size Diameter or "Foot print"	Coupling size
Green 	100mm - 3 15/16"	40mm - 1 9/16"
Red 	100mm - 3 15/16"	50mm - 2"
Blue 	110mm - 4 5/16"	60mm - 2 3/8"
Easiflex® Two-Piece skin barriers with flange		
Colour code	Skin barrier size Diameter or "Foot print"	Coupling size
Green 	100mm - 3 15/16"	40mm - 1 9/16"
Red 	110mm - 4 5/16"	50mm - 2"
Yellow 	130mm - 5 1/8"	70mm - 2 3/4"
Silver 	150mm - 5 7/8"	90mm - 3 1/2"

Quick Metric to Imperial Conversion Charts for Assura and Easiflex

Cut to fit	
Metric	Imperial
10-48mm	3/8" - 2"
10-68mm	3/8" - 2 3/8"
15-88mm	5/8" - 3 3/8"
10-35mm	3/8" - 1 3/8"
10-45mm	3/8" - 1 3/4"
10-55mm	3/8" - 2 1/8"
15-23mm	5/8" - 7/8"
15-33mm	5/8" - 1 1/4"
15-43mm	5/8" - 1 3/4"
15-61mm	5/8" - 2 3/8"

Pre-Cut	
Metric	Imperial
15mm	5/8"
18mm	3/4"
21mm	7/8"
25mm	1"
28mm	1 1/8"
31mm	1 1/4"
35mm	1 3/8"
38mm	1 1/2"
41mm	1 5/8"

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-820-8206 email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems

Non-convex=Flat Convex Light=Shallow Convexity Convex=Deep Convexity



Skin Barrier With Flange

Flexible, skin barrier with patented Assura spiral technology for secure, long-lasting wear. Belt loops on all skin barriers.

	Convexity	Colour Code	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	2881	5
		50mm - 2" ●	10-45mm - 3/8"-1 3/4"	2882	
		60mm - 2 3/8" ●	10-55mm - 3/8"-2 1/8"	2883	
	Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	14261	
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	14262	
		60mm - 2 3/8" ●	15-43mm - 5/8"-1 3/4"	14263	
	Convex	50mm - 2" ●	15-33mm - 5/8"-1 1/4"	12716	
		60mm - 2 3/8" ●	15-43mm - 5/8"-1 3/4"	12719	
Pre-Cut	Non-Convex	40mm - 1 9/16" ●	20mm - 3/4"	12841	
			25mm - 1"	12842	
		50mm - 2" ●	30mm - 1 1/8"	12844	
			35mm - 1 3/8"	12845	
		60mm - 2 3/8" ●	40mm - 1 9/16"	12847	
			Convex Light	40mm - 1 9/16" ●	18mm - 3/4"
	21mm - 7/8"	14272			
	50mm - 2" ●	25mm - 1"		14273	
		28mm - 1 1/8"		14274	
		31mm - 1 1/4"		14275	
	60mm - 2 3/8" ●	35mm - 1 3/8"		14276	
		38mm - 1 1/2"	14277		
	41mm - 1 5/8"	14278			
		Convex	40mm - 1 9/16" ●	15mm - 5/8"	12701
	18mm - 3/4"			12702	
	21mm - 7/8"			12703	
	50mm - 2" ●		25mm - 1"	12704	
			28mm - 1 1/8"	12705	
			31mm - 1 1/4"	12706	
	60mm - 2 3/8" ●		35mm - 1 3/8"	12707	
38mm - 1 1/2"			12708		
41mm - 1 5/8"	12709				

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems



Non-convex=Flat Convex Light=Shallow Convexity Convex=Deep Convexity

Assura Extra-Extended Wear Skin Barrier With Flange

Assura Extra-Extended Wear offers superior wear-time compared to traditional systems. Assura Extra-Extended Wear barrier comes with a unique transparent adhesive that provides security and flexibility while preserving the integrity of the surrounding skin. Patented Assura spiral technology.

Belt loops available on all skin barriers.

	Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	2831	5
		50mm - 2" ●	10-45mm - 3/8"-1 3/4"	2832	
		60mm - 2 3/8" ●	10-55mm - 3/8"-2 1/8"	2833	
	Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	14281	
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	14282	
		60mm - 2 3/8" ●	15-43mm - 5/8"-1 3/4"	14283	
	Convex	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	14243	
		50mm - 2" ●	15-33mm - 5/8"-1 1/4"	14246	
		60mm - 2 3/8" ●	15-43mm - 5/8"-1 3/4"	14249	
Pre-Cut	Non-Convex	40mm - 1 9/16" ●	20mm - 3/4"	14251	
			25mm - 1"	14252	
		50mm - 2" ●	30mm - 1 1/8"	14254	
			35mm - 1 3/8"	14255	
			60mm - 2 3/8" ●	45mm - 1 3/4"	14258
				Convex Light	40mm - 1 9/16" ●
	21mm - 7/8"	14292			
	50mm - 2" ●	25mm - 1"	14293		
		28mm - 1 1/8"	14294		
		31mm - 1 1/4"	14295		
		60mm - 2 3/8" ●	35mm - 1 3/8"		14296
	Convex	40mm - 1 9/16" ●	18mm - 3/4"	14232	
			21mm - 7/8"	14233	
		50mm - 2" ●	25mm - 1"	14234	
			28mm - 1 1/8"	14235	
			31mm - 1 1/4"	14236	
			60mm - 2 3/8" ●	35mm - 1 3/8"	14237
	38mm - 1 1/2"	14238			
		41mm - 1 5/8"	14239		

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-820-8206 email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems

Easiclose™ WIDE Outlet Drainable Pouch

Assura Easiclose WIDE Outlet drainable pouches are extremely hygienic and secure for people with colostomies and ileostomies. The WIDE Outlet makes emptying easier for people with thicker output. Slim line filter enhances discretion while controlling odor. Slick interior pouch film also makes emptying easier and faster. Velcro fixation on outlet creates a hygienic pouching solution.

Easiclose Drainable Pouches must be used with the corresponding colour coded skin barrier.



Mini Soft Front Classic Drainable Pouch - 23 cm (9")

Small, opaque, drainable pouch with soft cloth on front and back, and unique Easiclose outlet with integrated filter.

Coupling Size	Product Code	Units
40mm - 1 9/16" ●	13924	10
50mm - 2" ●	13925	



Midi Graphic Drainable Pouch - 26 cm (10 1/4")

Mid-sized, drainable pouch with printed soft cloth front and back, and unique Easiclose outlet with integrated filter.

Coupling Size	Product Code	Units
40mm - 1 9/16" ●	13964	10
50mm - 2" ●	13965	
60mm - 2 3/8" ●	13966	

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems



Maxi Drainable Transparent Pouch - 28 cm (11")

Large-size, transparent, drainable pouch with unique Easyclose outlet. Only available with filter.

Coupling Size	Product Code		Units
	With filter		
40mm - 1 9/16" ●	13974		10
50mm - 2" ●	13975		
60mm - 2 3/8" ●	13976		



Maxi Drainable Soft Front Classic Pouch - 28 cm (11")

Large-size, opaque, drainable pouch with soft cloth front and back and unique Easyclose outlet.

Coupling Size	Product Code		Units
	With filter	Without filter	
40mm - 1 9/16" ●	13984	14494	10
50mm - 2" ●	13985	N/A	
60mm - 2 3/8" ●	13986	14498	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems



Maxi Drainable Transparent Pouch - 30cm (12")

Can be rotated/repositioned for comfort without removing pouch. Includes 5 pouch clamps and 10 twist-ties. Only available without filter.

Coupling Size	Product Code	Units
50mm - 2" ●	2792	10
60mm - 2 3/8" ●	2793	



Maxi Multi-Chamber Urostomy Pouch - 25cm (10")

A unique and innovative design featuring multiple chambers to provide a more discreet and quiet urostomy pouch. Soft moisture-absorbent backing with anti-reflux valve and an audible click locking system. Can be rotated/repositioned for comfort without removing pouch and features a secure outlet closure and simple locking connector for night a drainage bag attachment. Includes 10 adaptors.

Pouch Colour	Coupling Size	Product Code	Units
Opaque	40mm - 1 9/16" ●	14224	10
	50mm - 2" ●	14225	
	60mm - 2 3/8" ●	14226	
Transparent	40mm - 1 9/16" ●	14227	
	50mm - 2" ●	14228	
	60mm - 2 3/8" ●	14229	

Coloplast Ostomy Care

Assura® Two-Piece Ostomy Systems



Closed Pouch

Soft, moisture-absorbent backing, odour-control filter and an audible click locking system. Can be rotated/repositioned for comfort without removing pouch.

Pouch Colour	Pouch Size	Coupling Size	Product Code	Units
Opaque	18cm (7")	40mm - 1 9/16" ●	12354	30
		50mm - 2" ●	12355	
		60mm - 2 3/8" ●	12356	
Transparent		40mm - 1 9/16" ●	12344	
		60mm - 2 3/8" ●	12346	
		Opaque	22cm (8 1/2")	
50mm - 2" ●	12385			
60mm - 2 3/8" ●	12386			
Transparent	40mm - 1 9/16" ●			12374
	50mm - 2" ●			12375
	60mm - 2 3/8" ●			12376



All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Easiflex® Two-Piece Ostomy Systems with Adhesive Coupling



Skin Barrier With Flange

The two-piece adhesive coupling system is a great option for a flexible lifestyle. This ostomy pouching system is based on adhesives that make it easy to apply the pouch to the skin barrier- and easy to remove. This pouching system is hygienic, easy to clean, soft and comfortable, safe, and best of all... simple to handle. Easiflex is a great solution for post-operative situations.

Convex Light skin barriers include belt loops.

		Convexity	Coupling Size	Stoma Size	Product Code	Units
Cut-to-Fit	Easiflex Extra-Extended Wear	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	14304	10
			50mm - 2" ●	10-48mm - 3/8"-2"	14305	
			70mm - 2 3/4" ●	10-68mm - 3/8"-2 3/8"	14306	
		Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	14601	5
			50mm - 2" ●	15-33mm - 5/8"-1 1/4"	14602	
			50mm - 2" ●	15-43mm - 5/8"-1 3/4"	14603	
	70mm - 2 3/4" ●		15-61mm - 5/8"-2 3/8"	14408		
	Easiflex	Non-Convex	40mm - 1 9/16" ●	10-35mm - 3/8"-1 3/8"	14301	10
			50mm - 2" ●	10-50mm - 3/8"-2"	14302	
			70mm - 2 3/4" ●	10-68mm - 3/8"-2 3/8"	14303	
			90mm - 3 1/2" ●	15-88mm - 5/8"-3 1/2"	14309	5
		Convex Light	40mm - 1 9/16" ●	15-23mm - 5/8"-7/8"	14401	5
50mm - 2" ●			15-33mm - 5/8"-1 1/4"	14402		
50mm - 2" ●			15-43mm - 5/8"-1 3/4"	14403		
70mm - 2 3/4" ●			15-61mm - 5/8"-2 3/8"	14406		
Pre-Cut	Easiflex	Convex Light	40mm - 1 9/16" ●	21mm - 7/8"	14642	5
			50mm - 2" ●	25mm - 1"	14643	
			50mm - 2" ●	28mm - 1 1/8"	14644	

Coloplast Ostomy Care

Easiflex® Two-Piece Ostomy Systems with Adhesive Coupling



Easicloose™ Drainable Pouch

An extremely hygienic and secure pouch for people with ileostomies or colostomies featuring an easy-to-use patented Velcro closure system. No clamps or clips required. Soft, comfort backing and a low-profile for maximum discretion. Easiflex Drainable Pouches must be used with the corresponding colour coded Easiflex Skin Barrier.

Maxi pouches available with or without charcoal filter.

Pouch Size	Pouch Colour	Coupling Size	Product Code		Units
			With filter	Without filter	
Midi 26cm (10 1/4")	Opaque	40mm - 1 9/16" ●	14346	N/A	20
		50mm - 2" ●	14347	N/A	
	Transparent	50mm - 2" ●	14342	N/A	
Maxi 28cm (11")	Opaque	40mm - 1 9/16" ●	14361	14524	
		50mm - 2" ●	14362	14525	
		70mm - 2 3/4" ●	14349	N/A	
	Transparent	40mm - 1 9/16" ●	14356	14520	
		50mm - 2" ●	14357	14521	
		70mm - 2 3/4" ●	14344	14528	
		90mm - 3 1/2" ●	NA	14527	30
XL 30cm (12")	Opaque	70mm - 2 3/4" ●	14364	N/A	20
	Transparent	70mm - 2 3/4" ●	14359	N/A	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com








Coloplast Ostomy Care

Easiflex® Two-Piece Ostomy Systems with Adhesive Coupling



Closed Pouch



Soft, comfort backing with charcoal filter for odour control. Easiflex Closed Pouches must be used with the corresponding colour coded Easiflex Skin Barrier. Only available with filter.

Pouch Size	Coupling Size	Product Code		Units
		Opaque	Transparent	
Midi 18cm (7")	40mm - 1 9/16" 	14316	NA	30
	50mm - 2" 	14317	NA	
Maxi 21cm (8 1/4")	40mm - 1 9/16" 	14331	NA	
	50mm - 2" 	14332	NA	
	70mm - 2 3/4" 	14319	NA	
XL 26cm (10 1/4")	90mm - 3 1/2" 	14355	NA	
	70mm - 2 3/4" 	NA	14329	



Maxi Urostomy Pouch - 25cm (10")

Easiflex Urostomy Pouches have a unique and innovative design featuring multiple chambers to provide a more discreet and quiet urostomy pouch. Soft moisture-absorbent backing and anti-reflux valve. Easiflex must be used with the corresponding colour coded Easiflex Skin Barrier.

Coupling Size	Product Code		Units
	Opaque	Transparent	
40mm - 1 9/16" 	14553	14552	20
50mm - 2" 	14555	14554	

Coloplast Ostomy Care

Assura® One-Piece Ostomy Systems

Extra-Extended Wear Pouch

Assura Extra-Extended Wear offers superior wear-time to traditional systems. The unique transparent adhesive offers security and flexibility while preserving the integrity of the surrounding skin.



Maxi Extra-Extended Wear Drainable Pouch - 28cm (11")

Includes 5 pouch clamps and 10 twist-tie closures. Only available without filter.

	Convexity	Pouch Colour	Stoma Size	Product Code	Units
Cut-to-Fit	Convex	Opaque	15-43mm 5/8" - 1 3/4"	12419	10
		Transparent		12416	



Extra-Extended Wear Transparent Urostomy Pouch 27 cm (10 1/2 ")

Secure outlet closure and simple locking connector for night drainage bag attachment. Includes 5 adaptors.

	Convexity	Pouch Size	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	150 mL	10-55mm 3/8" - 2 1/8"	12472	10
	Convex	375 mL	15-33mm 5/8" - 1 1/4"	12432	
			15-43mm 5/8" - 1 3/4"	12433	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® One-Piece Ostomy Systems



Maxi Extra-Extended Wear Multi-Chamber Urostomy Pouches - 25cm (10")

A unique and innovative design featuring multiple chambers to provide a more discreet and quiet urostomy pouch. Soft, moisture-absorbent backing and anti-reflux valve. Secure outlet closure and simple locking connector for night drainage bag attachment. Includes 5 adaptors.

	Convexity	Pouch Colour	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	Opaque/ 375 mL	10-55mm 3/8" - 2 1/8"	14201	10
		Transparent/ 375 mL		14202	
Cut-to-Fit	Convex Light	Opaque/ 375 mL	15-43mm 5/8" - 1 3/4"	14701	
		Transparent/ 375 mL		14702	

Assura Drainable Pouch With Easiclose™

An extremely hygienic and secure pouch for people with ileostomies or colostomies with an easy-to-use patented Velcro closure system. No clamps or clips required. Soft, comfortable backing.



Midi Soft Front Classic Pouch – 26cm (10 1/4")

Mid-size, opaque, drainable pouch with soft cloth on both sides, and unique Easiclose outlet with integrated filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-55mm - 3/8"-2 1/8"	13840	10
	Convex Light	15-33mm - 5/8"-1 1/4"	14404	
		15-43mm 5/8"-1 3/4"	14405	
	Convex	15-43mm 5/8"-1 3/4"	14196	
Pre-Cut	Non-Convex	25mm - 1"	13844	
		30mm - 1 1/8"	13845	
		35mm - 1 3/8"	13846	

Coloplast Ostomy Care

Assura® One-Piece Ostomy Systems



Maxi Soft Front Classic Pouch – 28cm (11")

Large, opaque, drainable pouch with soft cloth on both sides and a unique Easyclose outlet. Only available with filter.

	Convexity	Stoma Size	Product Code	Units
Cut-to-Fit	Non-Convex	10-70mm - 3/8"-2 3/4"	13870	10
	Convex Light	15-33mm - 5/8"-1 1/4"	14414	
		15-43mm 5/8"-1 3/4"	14415	
	Convex	15-43mm 5/8"-1 3/4"	14106	
Pre-Cut	Non-Convex	25mm - 1"	13874	
		30mm - 1 1/8"	13875	
		35mm - 1 3/8"	13876	
	Convex Light	25mm - 1"	14422	
		28mm - 1 1/8"	14423	
		31mm - 1 1/4"	14424	
	Convex	21mm - 7/8"	14173	
		25mm - 1"	14174	
		28mm - 1 1/8"	14175	
		31mm - 1 1/4"	14176	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® One-Piece Ostomy Systems



Maxi Drainable Transparent Pouch – 28cm (11")
 Large, transparent, drainable pouch with unique Easiclose outlet.

	Convexity	Stoma Size	Product Code		Units
			With filter		
Cut-to-Fit	Non-Convex	10-70mm - 3/8"-2 3/4"		13860	10
	Convex Light	15-33mm - 5/8"-1 1/4"		14411	
		15-43mm 5/8"-1 3/4"		14412	
	Convex	15-43mm 5/8"-1 3/4"		14103	
Pre-Cut	Convex	21mm - 7/8"		14163	
		25mm - 1"		14164	
		28mm - 1 1/8"		14165	
		31mm - 1 1/4"		14166	



Maxi Drainable Pouch - 30cm (12")
 Flexible, skin barrier with patented Assura spiral technology for secure, long-lasting wear. Soft backing for comfort. Odour-proof film. Convex skin barriers include belt loops. Includes 5 pouch clamps and 10 twist-tie closures. Without filter.

	Convexity	Stoma Size	Product Code		Units
			Opaque	Transparent	
Cut-to-Fit	Non-Convex	10-70mm - 3/8"-2 3/4"	2520	2490	10
	Convex	15-43mm 5/8"-1 3/4"	12533	13706	

Coloplast Ostomy Care

Assura® One-Piece Ostomy Systems



Urostomy Pouch - 27cm (10 1/2")

Flexible adhesive, soft backing for comfort and odour-proof film.
Anti-reflux valve and simple closure with adaptor for night drainage bag.
Convex skin barriers include belt loops. Includes 5 adaptors.

	Convexity	Pouch Capacity	Stoma Size	Product Code		Units
				Opaque	Transparent	
Cut-to-Fit	Non-Convex	375mL	10-55mm - 3/8"-2 1/8"	5580	5585	10
	Convex	375mL	15-33mm - 5/8"-1 1/4"	NA	12595	
			15-43mm - 5/8"-1 3/4"	NA	12596	
Pre-Cut	Convex	375mL	15mm - 5/8"	NA	12991	
			18mm - 3/4"	NA	12992	
			21mm - 7/8"	NA	12993	
			25mm - 1"	NA	12994	
			28mm - 1 1/8"	NA	12995	
			31mm - 1 1/4"	NA	12996	
			35mm - 1 3/8"	NA	12997	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Pediatric Ostomy System



Pediatric Skin Barrier With Flange

Flexible skin barrier for secure, long-lasting wear. Flexible flange allows pouch to be secured with minimal pressure. Pediatric skin barrier flange can be used with any green Assura pouch.

Cut-to-Fit Stoma Size	Product Code	Units
10-35mm - 3/8"-1 3/8"	2181	5



Pediatric Pouch

Soft backing comfort and odour proof film. Locking system reduces pressure to abdomen.

Pouch Type	Product Code		Units
	Transparent	Opaque	
Drainable Pouch 13cm - 5"	N/A	2151	10
Closed Pouch 13cm - 5"	2166	N/A	30

Assura® One-Piece Pediatric Ostomy System



Pediatric Pouch

Flexible and skin friendly barrier with filter for odour control.

The urostomy pouch has an odour-proof film with anti-reflux valve and soft backing for comfort.

Pouch Type	Product Code		Units
	Opaque	Transparent	
Drainable 15cm - 6" 10-35mm - 3/8"-1 3/8"	NA	2115	10
Closed 13cm - 5" 10-35mm - 3/8"-1 3/8"	NA	2125	30
Urostomy 15cm - 6" 150 mL 10-35mm - 3/8"-1 3/8"	NA	8009	10

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspect@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Pediatric Ostomy System



Pediatric Skin Barrier With Flange

Flexible skin barrier for secure, long-lasting wear. Flexible flange allows pouch to be secured with minimal pressure. Pediatric skin barrier flange can be used with any green Assura pouch.

Cut-to-Fit Stoma Size	Product Code	Units
10-35mm - 3/8"-1 3/8"	2181	5



Pediatric Pouch

Soft backing comfort and odour proof film. Locking system reduces pressure to abdomen.

Pouch Type	Product Code		Units
	Transparent	Opaque	
Closed Pouch 13cm - 5"	N/A	2151	10
Closed Pouch 13cm - 5"	2166	N/A	30

Assura® One-Piece Pediatric Ostomy System



Pediatric Pouch

Flexible and skin friendly barrier with filter for odour control.

The urostomy pouch has an odour-proof film with anti-reflux valve and soft backing for comfort.

Pouch Type	Product Code		Units
	Opaque	Transparent	
Drainable 15cm - 6" 10-35mm - 3/8"-1 3/8"	NA	2115	10
Closed 13cm - 5" 10-35mm - 3/8"-1 3/8"	NA	2125	30
Urostomy 15cm - 6" 150 mL 10-35mm - 3/8"-1 3/8"	NA	8009	10

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Assura® Two-Piece Pediatric Ostomy System



Pediatric Skin Barrier With Flange

Flexible skin barrier for secure, long-lasting wear. Flexible flange allows pouch to be secured with minimal pressure. Pediatric skin barrier flange can be used with any green Assura pouch.

Cut-to-Fit Stoma Size	Product Code	Units
10-35mm - 3/8"-1 3/8"	2181	5



Pediatric Pouch

Soft backing comfort and odour proof film. Locking system reduces pressure to abdomen.

Pouch Type	Product Code		Units
	Transparent	Opaque	
Drainable Pouch 13cm - 5"	N/A	2151	10
Closed Pouch 13cm - 5"	2166	N/A	30

Assura® One-Piece Pediatric Ostomy System



Pediatric Pouch

Flexible and skin friendly barrier with filter for odour control.

The urostomy pouch has an odour-proof film with anti-reflux valve and soft backing for comfort.

Pouch Type	Product Code		Units
	Opaque	Transparent	
Drainable 15cm - 6" 10-35mm - 3/8"-1 3/8"	NA	2115	10
Closed 13cm - 5" 10-35mm - 3/8"-1 3/8"	NA	2125	30
Urostomy 15cm - 6" 150 mL 10-35mm - 3/8"-1 3/8"	NA	8009	10

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspect@coloplast.com

Coloplast Ostomy Care

New Brava™ Ostomy Accessories

		Size	Product Code	Units
	Brava Protective Sheets Effectively protects peristomal skin by providing a ready surface for good adhesion.	10 x 10 cm 4" x 4"	3210	10
		15 x 15 cm 6" x 6"	3215	5
		20 x 20 cm 8" x 8"	3220	
	Brava Mouldable Rings Easy to handle, mould and apply for a tight seal between the stoma and the appliance.	2.0 mm	12030	30
		4.2 mm	12042	
	Brava Skin Barrier Spray and Wipe Protects against body waste and adhesives Sting free	Spray 50 ml	12020	1
		Single Application Wipe	12021	30
	Brava Adhesive Remover Spray and Wipe For easy removal of adhesive Sting free	Spray 50 ml	12010	1
		Single Application Wipe	12011	30
	Brava Ostomy Strip Paste Fills out uneven skin surfaces for a better skin barrier fit. Alcohol-Free	57g - 2oz	2655	10

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

New Brava™ Ostomy Accessories

		Size	Product Code	Units
	Brava Elastic Tape Increases skin barrier wear time. Skin friendly and elastic to follow natural body movements	Single strip	12070	20
	Brava Ostomy Powder Absorbs moisture to keep skin dry and reduce irritation	25g - 1oz	1907	16
	Brava Belt Ensures position of skin barrier for extra security. Discreet, soft and comfortable to wear. Latex free	Regular 100 cm - 40"	0421	15
		XL 135cm - 54"	0422	
	Brava Lubricating Deodorant Neutralizes odour with a multi-lubricating effect to ensure that the content is always at the bottom of the pouch and not around the stoma.	7.5mL	12060	20
		237mL 8oz	12061	1

Everybody can need a little extra security and comfort

We all have different needs.

The shape of your body and the characteristics of your skin are as unique as you personality. That's why we offer a little extra help to make your ostomy bag fit securely and to care for the skin around your stoma.



Coloplast Ostomy Care

Coloplast Accessories

		Size	Product Code	Units
	<p>Ostomy Paste Provides a tight seal around the stoma to protect against output.</p> <p>Low Alcohol Content</p>	60 g	2650	12
	<p>Skin Barrier Rings Soft and flexible rings provide skin protection for stomas, fistulas and draining wounds.</p>	20mm - 3/4"	2320	30
		25mm - 1"	2325	
		30mm - 1 1/8"	2330	
		40mm - 1 9/16"	2340	
		50mm - 2"	2350	
	<p>Coloplast Prep™ Protective Liquid Skin Barrier Clear film forms a breathable barrier protecting against skin irritants such as adhesives, urine and feces. No need to remove between applications.</p> <p>Not to be used with Extra Extended Wear Skin Barriers</p>	Single Application Wipe	2041	12 boxes / 54
		59mL - 2oz	0925	12
	<p>Coloplast Filtrodor® Pouch Filter Charcoal-activated deodorizing filter with self-sealing foam center. Adheres to pouch.</p>		0509	5 pouches (10 units per pouch)
	<p>Coloplast Pouch Clamp Simple closure clamp.</p>		9500	20
	<p>Ostomy Scissors Stainless steel medical grade ostomy scissors with rounded ends.</p>		9505	30

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

Coloplast Accessories



Coloplast Uro Minicap

The Coloplast Uro Minicap was designed for discreet situations. Soft front and super-absorbent material holds up to 100mL of urine, lasting between 30 minutes and an hour. Features an audible click locking system and odour barrier. Compatible with Sensura Click and Assura 2 piece skin barriers.

Type	Coupling Size	Product Code	Units
Opaque	40mm - 1 9/16" 	2807	30
	50mm - 2" 	2808	



Coloplast Stoma Cap

The Coloplast Stoma Cap was designed for discreet situations. Mini security pouch with filter features an audible click locking system and can be used with Sensura Click and Assura 2 piece skin barriers.

Colour	Coupling Size	Product Code	Units
Opaque	40mm - 1 9/16" 	2801	30
	50mm - 2" 	2802	
	60mm - 2 3/8" 	2803	



Coloplast Stoma Cap - 1 piece


The Coloplast Stoma Cap was designed for discreet situations. The 1-piece mini security pouch system is only available with filter.

Cut-to-Fit Stoma Size	Opaque	Transparent	Units
20-55mm - 3/4"-2 1/8"	2501	NA	30



Coloplast Ostomy Care

Assura® Irrigation Products


Assura® Irrigation Set

Assura® Irrigation Set		Contents	Product Code	Units
	Economy Version	Includes: 1 Water Bag/Regulator 1 Cone	1500	1



Assura® Irrigation Accessories

Assura® Irrigation Accessories		Description	Product Code	Units
	Stoma Cone	Soft, flexible stoma cone with tubing	1110	4
	Irrigation Faceplate	For use with ostomy Belt. (Use with #12836 Irrigation Sleeve).	12820	1

Assura® Irrigation Sleeves

Assura® Irrigation Sleeves		Description	Product Code	Units
	Transparent, disposable irrigation sleeve with adhesive back.	1-pc sleeve	1003	10

The irrigation sleeves below must be used with their corresponding colour coded flange or faceplate. Consult this chart to order the correct transparent sleeves.

Description	Coupling Size	Use with Assura® Product Code	Product Code	Units
Irrigation sleeves	50mm - 2" 	2832,2882,12704-12706, 12844-12845, 12716, 14246	12835	5
	60mm - 2 3/8" 	2833, 2883, 12707-12709, 12847, 12719, 12820, 14249	12836	

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspeg@coloplast.com

Coloplast Ostomy Care

SenSura Hospital Assortment



SenSura 1-Piece Drainable Pouch

SenSura 1-piece drainable pouch has the features of a post-op appliance, but offers the discretion of a standard pouch. The SenSura 1-piece drainable pouches are connectable with different night bags, come without filter and are available in sterile and non-sterile.

	Sterility	Stoma Size	Volume	Product Code	Units
Cut-to-Fit	Sterile*	10 - 76mm - 3/8"-3"	630 ml	19001	5
	Non-sterile	10 - 76mm - 3/8"-3"	630 ml	19003	6



SenSura Post-Op Pouch

SenSura Post-Op pouch offers easy access to the stoma through a window, which is fixed to the pouch and easy to handle. The SenSura Post-Op pouches are available in sterile and non-sterile, with 70 and 100 mm windows, as well as without window. The soft outlet connects easily to various night bags.

	Sterility	Stoma Size	Window Size	Volume	Product Code	Units
Cut-to-Fit	Sterile*	10 - 76mm - 3/8"-3"	70 mm	680 ml	19010	5
		10 - 76mm - 3/8"-3"	100 mm	725 ml	19011	
	Non-sterile	10 - 76mm - 3/8"-3"	70 mm	680 ml	19020	6
		10 - 76mm - 3/8"-3"	100 mm	725 ml	19021	
		10 - 76mm - 3/8"-3"	No window	680 ml	19030	



SenSura High Output Pouch

SenSura High Output pouch is the solution for people with a need for extra capacity. The SenSura High Output Pouches are available in transparent and opaque for each coupling size. They ensure effective drainage through a wide, funnel shaped opening and the option of outside connection to night bags. Only available with filter.

Coupling Size	Volume	Product Code		Units
		Opaque	Transparent	
40mm - 1 9/16" ●	1160 ml	19040	19050	10
50mm - 2" ●	1120 ml	19041	19051	
60mm - 2 3/8" ●	1080 ml	19042	19052	


* Sterile products available March 2013

Coloplast Ostomy Care

Assura Post Op


Two-piece High Output Drainable System

Two-piece drainable system for the first days following ostomy surgery when effluent output is high. Longer wear-time for fewer night time disturbances. Available in 2L Ileo Night Pouch for initial use. Skin barriers sold separately.


	Coupling Size	Stoma Size	Product Code	Units
 <p>ILEO NIGHT POUCH 2L (transparent)</p>	50mm - 2" ●	13-50mm - 9/16" - 2"	2836	5

One-piece Post-op Drainable Pouch Cut-to-Fit

Large 12" (30cm), transparent, one-piece pouch for use in post-surgical ostomy and wound management. Flexible skin barrier can be cut to fit up to 70 or 100 mm.

	Pouch Type	Stoma Size	Product Code	Units	
	Non-Sterile	With Window	13-70mm - 9/16"-2 3/4"	5	
		Without Window	13-100mm - 9/16"-3 15/16"		
	Sterile	With Window	13-70mm - 9/16"-2 3/4"		12805
		Without Window	13-100mm - 9/16"-3 15/16"		12808

Coloplast Urostomy Night Bag

	Product Code	Units
 <p>Includes a fixed connector that fits on to the urostomy pouch. The bag holds up to 2 litres and has a 120cm long flexible tube with anti-reflux valve.</p>	21365	10

All Coloplast Ostomy products are LATEX-FREE

Customer Service Department Tel: 1-888-880-8605 - Fax: 1-877-820-8206 - email: cacustomerservice@coloplast.com
 Consumer Support Department Tel: 1-866-293-6349 - Fax: 1-877-820-8206 - email: ca_conspect@coloplast.com

Coloplast Ostomy Care

Coloplast Fistula & Wound Management System



The Coloplast Fistula and Wound Management System is designed with a flexible adhesive. The pouch conforms to unique abdominal contours and molds to the skin with ease. This revolutionary "all in one" solution provides improved patient comfort and saves nurses' time.

Pouch Type	Cutting surface	Product Code	Units
Mini	104-159mm - 4 1/8"-6 1/4"	14050	6
Midi	156-228mm - 6 3/16"-9"	14060	3
Maxi	208-297mm - 8 1/4"-11 3/4"	14070	

Each retail box unit comes with 1 drain port, 1 flexible lid, 1 sheet, 1 clamp. There is also 1 pump per retail box (not sold separately).



Flexible Lids For The Coloplast Fistula & Wound Management System

The system lid is easy to apply and remove allowing nurses access to dress fistulas/wounds.

Pouch Type	For	Product Code	Units
Mini	Mini Fistula & Wound Management system	14021	10
Midi	Midi Fistula & Wound Management system	14031	
Maxi	Maxi Fistula & Wound Management system	14041	

If a filter is needed please use the Filtrodor®, sold separately.



Accessories For Coloplast Fistula & Wound Management System

The Drain Port can be easily applied to the lid to allow for catheter insertion.

Pouch Type	For	Product Code	Units
Bed Drainage Bag	All Fistula & Wound Management systems (Capacity 2000 mL)	14010	10
Drain Port	All Fistula & Wound Management systems	14015	12

Coloplast Ostomy Care

Coloplast Support Services



“Every individual is unique and may need to try a number of products to find the best solution”

Change is never easy, but it’s important to know you are not alone. For over 50 years it’s been our mission to make life easier for people with intimate healthcare needs. That’s why we’ve developed the following programs and services:

Sample Program:

You have unique needs that may require a unique solution. We want to help you find that solution. The Coloplast Sample program provides you with the opportunity to trial our products before purchasing them.

Consumer Support:

You may have a lot of questions. Some might be as simple as finding out where you can purchase your ostomy products, or some might be more intimate. Whatever your questions, we have a knowledgeable, courteous staff to answer your question and provide you with the information and literature you need.

Call **1-866-293-6349** or visit us online at **www.coloplast.ca** to obtain samples or to reach consumer support.

Hours of Operation:

Monday – Friday: 9am – 6pm (EST)

Continence Care



Continence Care

SpeediCath®

Intermittent Catheters

The only ready-to-use hydrophilic catheter. Clinical research indicated that twice as many individuals using SpeediCath were free of UTI's compared to catheters with gel¹. SpeediCath catheters have polished eyelets and are all latex, phthalate and PVC Free. Now with new and improved packaging that makes it easier to open and with less risk of spillage.

		Size (colour indicator)	Product Code	Units
 <p>SpeediCath Male Nelaton Intermittent Catheter 16" long</p>		8 FR (blue)	27408	30
		10 FR (black)	27410	
		12 FR (white)	27412	
		14 FR (green)	27414	
		16 FR (orange)	27416	
		18 FR (red)	27418	
<p>SpeediCath Male Coude Intermittent Catheter 16" long</p>		10 FR (black)	27490	
		12 FR (white)	27492	
		14 FR (green)	27494	
<p>SpeediCath Female Nelaton Intermittent Catheter 8" long</p>		8 FR (blue)	27508	
		10 FR (black)	27510	
		12 FR (white)	27512	
		14 FR (green)	27514	
		16 FR (orange)	27516	
<p>SpeediCath Pediatric Nelaton Intermittent Catheter 8" long</p>		6 FR (lt. green)	27706	
		8 FR (blue)	27708	
		10 FR (black)	27710	


1. Intermittent Catheterisation with Hydrophilic-Coated Catheters (SpeediCath) Reduces the Risk of Clinical Urinary Tract Infection in Spinal Cord Injured Patients: A Prospective Randomised Parallel Comparative Trial. D.J.M.K. De Ridder et al. European Urology 48 (2005) pp 991-995.

SpeediCath® Complete

Intermittent Catheters

The only ready-to-use hydrophilic catheter now with a non-touch handle and insertion guide, complete with an integrated collection bag.

Latex, phthalate and PVC Free.

		Size (colour indicator)	Product Code	Units
 <p>SpeediCath Complete Male Intermittent Catheter with Non-Touch Handle and Collection bag</p>		12 FR (white)	28462	20
		14 FR (green)	28464	


Continence Care

SpeediCath® Compact Male

SpeediCath® Compact Male is designed for discretion and practicality, both in and out of the home. Sleek and compact, it's less than half the size of standard catheters. This makes it easier to store, carry, use and dispose.

Latex, phthalate and PVC free.




	Size (incl. connector)	Product Code	Units
<p>SpeediCath Compact Male Male intermittent catheter 33 cm long (including connector)</p> 	12/18 FR (one size)	28692	30

SpeediCath® Compact Female

Intermittent Catheters

SpeediCath Compact is Coloplast's ready-to-use intermittent hydrophilic catheter designed specifically for women.

Latex, phthalate and PVC free.

	Size	Product Code	Units
<p>SpeediCath Compact Female Female intermittent catheter Under 6" / 15 cm long</p> 	8 FR	28578	30
	10 FR	28580	30
	12 FR	28582	30
	14 FR	28584	30

Continence Care

Self-Cath®

Intermittent Catheters

Siliconized catheters with fire-polished eyelets for intermittent catheterization. 100% Latex-free. Sterile.

Self-Cath		Size	Product Code	Reference Code	Units
Male Straight Tipped 16" (15cm) long 		8 FR	408	504500	30
		10 FR	410	504510	
		12 FR	412	504520	
		14 FR	414	504530	
		16 FR	416	504540	
		18 FR	418	504550	
Male Straight Tipped Curved Packaging 16" (15cm) long 		12 FR	460	502200	50
		14 FR	450	502210	
Soft-Straight Tipped 16" (15cm) long 		10 FR	110	501800	30
		12 FR	112	501810	
		14 FR	114	501820	
		16 FR	116	501830	
Coude Tapered Tip with guide stripe 16"(15cm)long 		8 FR	608	504650	30
		10 FR	610	504660	
		12 FR	612	504670	
		14 FR	614	504680	
Coude Olive Tip with Guide Stripe 16" (15cm) long 		6 FR	806	504700	30
		8 FR	808	504710	
		10 FR	810	504720	
		12 FR	812	504730	
		14 FR	814	504740	
		16 FR	816	504750	
		18 FR	818	504760	

Continence Care

Self-Cath®

Intermittent Catheters

Siliconized catheters with fire-polished eyelets for intermittent catheterization. 100% Latex-free. Sterile.

	Size	Product Code	Reference Code	Units	
Self-Cath Female Straight Tipped 6" (15cm) long 	8 FR	208	504600	30	
	10 FR	210	504610		
	12 FR	212	504620		
	14 FR	214	504630		
Self-Cath Female Straight Tipped 6" (15cm) long	14 FR	240*	501960		
Self-Cath Pediatric Catheters 10" (15cm) long 	5 FR	305	504400		
	6 FR	306	504410		
	8 FR	308*	504420		
	10 FR	310*	504430		
Self cath Extension Tube	26"	475	507850		12

* without connectors

Self-Cath® Plus

Intermittent Catheters

Siliconized catheters with fire-polished eyelets for intermittent catheterization. The hydrophilic coating activates with exposure to water for low friction catheterization. The Self-Cath Plus has an uncoated grip zone for better control and handling. 100% Latex-free. Sterile

	Size	Product Code	Reference Code	Units
Self-Cath Plus Hydrophilic Catheters Male 	10 FR	4410	504910	30
	12 FR	4412	504920	
	14 FR	4414	504930	

Continence Care

Self-Cath®

Closed System

The Self-Cath Closed System is a system that provides a safe, simple catheterization designed to reduce the risk of urinary tract infections. It includes a 100% latex-free Self-Cath catheter.

	Size	Product Code	Reference Code	Units
Self-Cath Closed System Straight tip	8 FR	1108	502700	50
	12 FR	1112	502720	50
	14 FR	1114	502730	50
	16 FR	1116	502740	50
Self-Cath Closed Coude Tapered Tip with Guide Stripe	14 FR	2614	502850	50



Contenance Care

Conveen® Optima

Self-Sealing Male External Catheter

Latex Free

Self-sealing male external catheter features push ring for secure connection to a leg bag and an anti-kink bulb to prevent backflow of urine. Integrated skin-friendly adhesive is gentle, yet stays in place. Pull-tab enables easy application. Transparent to allow easy monitoring of the skin.

	Size	Product Code	Units
Conveen Optima Self-Sealing Male External Catheter Standard 	25 mm	22025	30
	30 mm	22030	30
	35 mm	22035	30
	40 mm	22040	30
Conveen Optima Self-Sealing Male External Catheter Short	21 mm	22121	30
	25 mm	22125	30
	30 mm	22130	30
	35 mm	22135	30

Conveen®

Male External Catheter

Latex Free

	Size	Product Code	Units
Conveen Security+ Self Sealing Male External Catheter 	21 mm (red)	5221	35
	25 mm (yellow)	5225	35
	30 mm (green)	5230	35
	35 mm (blue)	5235	35
	40 mm (grey)	5240	35
	25 mm (yellow)	5225H	100
	30 mm (green)	5230H	100
	35 mm (blue)	5235H	100
	40 mm (grey)	5240H	100
Conveen Security+ Male External Catheter and Adhesive Liner	30 mm	5030	35
	35 mm	5035	35
	25 mm	5025H	100
	30 mm	5030H	100
	35 mm	5035H	100

Continence Care

Clear Advantage®

Male External Catheters

		Size	Product Code	Reference Code	Units	Suggested Retail Price
Clear Advantage Silicone Self-Adhering Male External Catheter		Small - 23 mm	C6100	504311	100	
		Medium - 28 mm	6200	504331		
		Intermediate - 31 mm	6300	504341		
		Large - 35 mm	6400	504351		
		X-Large - 40 mm	6500	504361		

Freedom Clear™

Male External Catheters

		Size	Product Code	Reference Code	Units
Freedom Clear Non-Latex Self-Adhering Male External Catheter		Small - 23 mm	5100	505711	100
		Medium - 28 mm	C5200	505731	
		Intermediate - 31 mm	5300	505741	
		Large - 35 mm	5400	505751	
		X-Large - 40 mm	5500	505761	

Freedom Clear™ 'LS' Long Seal


Is designed with a longer length adhesive seal on a thin sheath providing added security and comfort.

		Size	Product Code	Reference Code	Units
Freedom Clear 'LS' Long Seal Silicone Self-Adhering Male External Catheter		Small - 23 mm	5190	505411	100
		Medium - 28mm	5290	505421	
		Intermediate - 31 mm	5390	505431	
		Large - 35 mm	5490	505341	
		X-Large - 40 mm	5590	505451	

Continence Care

Freedom Clear™ 'SS' Sport Sheath


Offers a thin sheath for less roll ring at the base of the catheter. It is a good choice for men with a retracted penis.

		Size	Product Code	Reference Code	Units
Freedom Clear 'SS' Sport Sheath Silicone Self-Adhering Male External Catheter		Small - 23 mm	5110	505311	100
		Medium - 28 mm	C5210	505321	
		Intermediate - 31 mm	C5310	505331	
		Large - 35 mm	C5410	505441	

Conveen®

Male External Catheters


Latex

		Size	Product Code	Units
Conveen Self-Sealing Male External Catheters Contains Latex		20 mm	5212	35
		25 mm	5200	
		30 mm	5205	
		35 mm	5210	
		40 mm	5215	
		40 mm	5215H	100
Conveen Male External Catheters and Adhesive Liner Contains Latex		25 mm	5125	35
		30 mm	5130	
		35 mm	5135	

Continence Care

Active Cath®

One-piece, self adhering, extended wear catheter combines a comfortable latex sheath with a wide, water-tight seal.

		Size	Product Code	Reference Code	Units
Active Cath Self-Adhering Male External Catheter Contains Latex		Small - 23 mm	8100	506220	100
		Medium - 28 mm	8300	506200	
		Intermediate - 31 mm	8305	506240	
		Large - 35 mm	8500	506260	

Freedom Cath®

One-piece self adhering male external catheter for standard every day use.

		Size	Product Code	Reference Code	Units
Freedom Cath Self-Adhering Male External Catheter Contains Latex		Small - 23 mm	8000	506110	100
		Medium - 28 mm	8200	506130	
		Intermediate - 31 mm	8205	506140	
		Large - 35 mm	8400	506150	



Continence Care

Conveen® Urine Bags

	Size	Product Code	Units
<p>Conveen Security+ Leg Bag® with Non-Latex Fabric Leg Bag Straps</p> <p>Lightweight leg bag includes pre-attached, adjustable, anti-kink tubing, soft cloth backing, anti-reflux valve, and centre weld for even distribution of urine.</p> <p>Latex and PVC - Free</p> 	17 oz. (500 ml) Medium with 20" tubing	5161	10
	26 oz. (750 ml) Large with 20" tubing (no center weld)	5167	
<p>Conveen® Security+ Contour Leg Bag with Non-Latex Fabric Leg Bag Straps</p> <p>Deluxe leg bag with pre-attached adjustable, standard tubing. With soft cloth backing and anti-reflux valve.</p> <p>Latex - Free</p> 	26 oz. (750 ml)	5170	
	28 oz. (800 ml) Extra Large with 18" tubing	5174	
<p>Conveen® Security+ Extra Large Leg Bag/Bedside Drainage Bag</p> <p>Extra large leg bag/drainage bag with anti-reflux valve, volume markings, pre-attached, standard tubing, and soft cloth backing. Secure, user-friendly closure. (Leg bag straps sold separately)</p> <p>Latex - Free</p> 	51 oz. (1500 ml) Extra Large with 36" tubing	5062	
<p>Conveen® Standard Sterile Bag</p> <p>Sterile, extra large leg bag/drainage bag with anti-reflux valve, volume markings, pre-attached, standard tubing, and soft cloth backing. Secure, user-friendly closure. (Leg bag straps sold separately)</p> <p>Sterile and Latex - Free</p>	51 oz. (1500 ml) Extra Large with 36" tubing	5063	

Continence Care

Freedom® Leg Bags

A secure, comfortable, soft, sterile, vinyl bag. Available in 500ml and 1000 ml sizes with 3 easy drain ports: t-tap, fas-tap and twist. Kits include bag with 18" latex free tubing, connectors, and velcro straps.

	Size	Product Code	Reference Code	Units
Freedom Flocked Back Fas-Tap Drainage Port Leg Bags & Kit 	500ml	9005	508750	12 Kits/Box
	1000ml	9007	508760	
Freedom T-Tap Drainage Port Leg Bags & Kit 	500ml	7078	508130	
	1000ml	7088	508240	
Freedom Twist Drainage Port Leg Bags & Kit 	500ml	C7094	508420	
	1000ml	7096	508430	

Freedom® Accessories

	Size	Product Code	Reference Code	Units
Clear Non-Latex Tube & Connector Kit	18"	7092	508350	12 Kits/Box
Comfort Leg Bag Straps - Adjustable / Stretchable, Latex-Free	1"	7090	508300	1 Pair

Moveen® Bedside Night Bag


	Size	Product Code	Units
Moveen® Bedside Night Bag Large volume (2 L) bedside drainage bag. Non-Sterile package includes 30 bags with anti-kink tubing and 1 hanger. Sterile package includes 10 bags with anti-kink tubing and 2 hangers. Latex-Free 	(Non-Sterile) 2L	21346	30
	(Sterile) 2L	21356	10

Continence Care

Accessories

		Size	Product Code	Units
<p>Moveen® Bag Hanger</p> <p>Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability.</p>		One-size	5070	1
<p>Conveen® Security+ Fabric Leg Bag Straps</p> <p>Soft, elastic fabric straps that are reusable. Latex-Free</p>		25" pair	50501	10






Peristeen™ Anal Plug

		Size	Product Code	Units
<p>Peristeen Anal Plug</p> <p>A simple and effective aid to stop faecal incontinence. The Peristeen anal plug is made from a comfortable foam that retains its softness in use. Now with a less absorbent string for greater discretion. Each box comes with a tube of lubricant Latex Free</p>		Small	1450	20
		Large	1451	

Continence Care

Peristeen® Anal Irrigation

Peristeen Anal Irrigation is a unique way of emptying the bowel and is used to prevent faecal incontinence, constipation and to reduce the amount of time spent on bowel management. The system has been specially designed to make it portable and easy to use.

	Product Code	Description	Units
	29121	System * 90 irrigations = 6 months	1 control unit 2 rectal catheters 1 bag (with lid) 2 straps 1 toiletries bag
	29122	Accessory unit * 15 irrigations = 1 month	15 rectal catheters 1 bag (no lid)
	29123	Rectal Catheter Single use only	10 rectal catheters
	29124	Leg Strap	1 set of 2 straps
	29125	Tube	2 tubes with blue connectors

* if used every other day



Peristeen is only available for direct purchase from Coloplast Canada by customers. To order Peristeen, contact our Consumer Support Team at 1-866-293-6349 or email us at ca_conspeg@coloplast.com

Continence Care

Folysil

Indwelling Foley Catheters

The Folysil 100% silicone indwelling catheter features a significantly larger lumen than latex versions of the same FR size. This helps to reduce obstruction and inhibit encrustation. The larger diameter eyelets allow for better drainage while an atraumatic tip adds to ease of insertion. The Folysil indwelling catheter also features a radiopaque stripe in a transparent shaft for easy detection under x-ray.



	Size	Product Code	Units
Folysil (indwelling catheters) 2-way, straight, 5-15cc, 41 cm. Latex free	12 FR	AA6112	5
	14 FR	AA6114	5
	16 FR	AA6116	5
	18 FR	AA6118	5
	20 FR	AA6120	5
	22 FR	AA6122	5
	24 FR	AA6124	5
Folysil (indwelling catheters) 2-way, straight, 30cc, 41 cm. Latex free	18 FR	AA6C18	5
	20 FR	AA6C20	5
	22 FR	AA6C22	5
	24 FR	AA6C24	5
Folysil (indwelling catheters) Paediatric, 2-way, straight, 1.5-3cc, , 30-35 cm. Latex free	6 FR, 1.5cc, 30cm	AA6106	5
	8 FR, 3cc, 30cm	AA6108	5
	10 FR, 3cc, 35cm	AA6110	5
Folysil (indwelling catheters) 2 way, Tiemann (coudé) tip, 15cc, 2-way, 41cm. Latex free	12 FR	AA6312	5
	14 FR	AA6314	5
	16 FR	AA6316	5
	18 FR	AA6318	5
	20 FR	AA6320	5
	22 FR	AA6322	5
	24 FR	AA6324	5
Folysil (indwelling catheters) Paediatric, 2 way, Tiemann (coudé) tip, 3cc, 30-35cm. Latex free	8 FR, 30cm	AA6308	5
	10 FR, 35cm	AA6310	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

Continence Care

X-Flow

Prostatic Catheters

Performance and safety

Wide irrigation and drainage channels to prevent blocking. Large bevels and eyes to ensure excellent drainage capacity. Stiff walls to facilitate evacuation.

All silicone, to avoid any risk of allergy or toxicity. Clear catheter body to monitor drainage quality.

Closed round tip for ease of placement of an introducer if needed

	Size	Product Code	Units
X-Flow (Prostatic Catheter) 3-way, straight, 50cc, 42 cm (Latex free)	18 FR	AB6018	5
	20 FR	AB6020	5
	22 FR	AB6022	5
	24 FR	AB6024	5
X-Flow (Prostatic Catheter) 3-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6118	5
	20 FR	AB6120	5
	22 FR	AB6122	5
	24 FR	AB6124	5
X-Flow (Prostatic Catheter) 3-way, Delinotte (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6218	5
	20 FR	AB6220	5
	22 FR	AB6222	5
	24 FR	AB6224	5
X-Flow (Prostatic Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6318	5
	20 FR	AB6320	5
	22 FR	AB6322	5
	24 FR	AB6324	5
X-Flow (Prostatic Catheter) 2-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6518	5
	20 FR	AB6520	5
	22 FR	AB6522	5
	24 FR	AB6524	5
X-Flow (Haematuria Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm) Curved tip allows easy introduction without stylet. Same characteristics as X-Flow prostatic catheter (AB63XX) with 30 cc balloon	18 FR	AB6A18	5
	20 FR	AB6A20	5
	22 FR	AB6A22	5
	24 FR	AB6A24	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Continence Care

Folysil

Indwelling Foley Catheters

The Folysil 100% silicone indwelling catheter features a significantly larger lumen than latex versions of the same FR size. This helps to reduce obstruction and inhibit encrustation. The larger diameter eyelets allow for better drainage while an atraumatic tip adds to ease of insertion. The Folysil indwelling catheter also features a radiopaque stripe in a transparent shaft for easy detection under x-ray.



	Size	Product Code	Units
Folysil (indwelling catheters) 2-way, straight, 5-15cc, 41 cm. Latex free	12 FR	AA6112	5
	14 FR	AA6114	5
	16 FR	AA6116	5
	18 FR	AA6118	5
	20 FR	AA6120	5
	22 FR	AA6122	5
	24 FR	AA6124	5
Folysil (indwelling catheters) 2-way, straight, 30cc, 41 cm. Latex free	18 FR	AA6C18	5
	20 FR	AA6C20	5
	22 FR	AA6C22	5
	24 FR	AA6C24	5
Folysil (indwelling catheters) Paediatric, 2-way, straight, 1.5-3cc, , 30-35 cm. Latex free	6 FR, 1.5cc, 30cm	AA6106	5
	8 FR, 3cc, 30cm	AA6108	5
	10 FR, 3cc, 35cm	AA6110	5
Folysil (indwelling catheters) 2 way, Tiemann (coudé) tip, 15cc, 2-way, 41cm. Latex free	12 FR	AA6312	5
	14 FR	AA6314	5
	16 FR	AA6316	5
	18 FR	AA6318	5
	20 FR	AA6320	5
	22 FR	AA6322	5
	24 FR	AA6324	5
Folysil (indwelling catheters) Paediatric, 2 way, Tiemann (coudé) tip, 3cc, 30-35cm. Latex free	8 FR, 30cm	AA6308	5
	10 FR, 35cm	AA6310	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

Continence Care

X-Flow

Prostatic Catheters

Performance and safety

Wide irrigation and drainage channels to prevent blocking. Large bevels and eyes to ensure excellent drainage capacity. Stiff walls to facilitate evacuation.

All silicone, to avoid any risk of allergy or toxicity. Clear catheter body to monitor drainage quality.

Closed round tip for ease of placement of an introducer if needed

	Size	Product Code	Units
X-Flow (Prostatic Catheter) 3-way, straight, 50cc, 42 cm (Latex free)	18 FR	AB6018	5
	20 FR	AB6020	5
	22 FR	AB6022	5
	24 FR	AB6024	5
X-Flow (Prostatic Catheter) 3-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6118	5
	20 FR	AB6120	5
	22 FR	AB6122	5
	24 FR	AB6124	5
X-Flow (Prostatic Catheter) 3-way, Delinotte (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6218	5
	20 FR	AB6220	5
	22 FR	AB6222	5
	24 FR	AB6224	5
X-Flow (Prostatic Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6318	5
	20 FR	AB6320	5
	22 FR	AB6322	5
	24 FR	AB6324	5
X-Flow (Prostatic Catheter) 2-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6518	5
	20 FR	AB6520	5
	22 FR	AB6522	5
	24 FR	AB6524	5
X-Flow (Haematuria Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm) Curved tip allows easy introduction without stylet. Same characteristics as X-Flow prostatic catheter (AB63XX) with 30 cc balloon	18 FR	AB6A18	5
	20 FR	AB6A20	5
	22 FR	AB6A22	5
	24 FR	AB6A24	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Continence Care

Folysil

Indwelling Foley Catheters

The Folysil 100% silicone indwelling catheter features a significantly larger lumen than latex versions of the same FR size. This helps to reduce obstruction and inhibit encrustation. The larger diameter eyelets allow for better drainage while an atraumatic tip adds to ease of insertion. The Folysil indwelling catheter also features a radiopaque stripe in a transparent shaft for easy detection under x-ray.



	Size	Product Code	Units
Folysil (indwelling catheters) 2-way, straight, 5-15cc, 41 cm. Latex free	12 FR	AA6112	5
	14 FR	AA6114	5
	16 FR	AA6116	5
	18 FR	AA6118	5
	20 FR	AA6120	5
	22 FR	AA6122	5
	24 FR	AA6124	5
Folysil (indwelling catheters) 2-way, straight, 30cc, 41 cm. Latex free	18 FR	AA6C18	5
	20 FR	AA6C20	5
	22 FR	AA6C22	5
	24 FR	AA6C24	5
Folysil (indwelling catheters) Paediatric, 2-way, straight, 1.5-3cc, , 30-35 cm. Latex free	6 FR, 1.5cc, 30cm	AA6106	5
	8 FR, 3cc, 30cm	AA6108	5
	10 FR, 3cc, 35cm	AA6110	5
Folysil (indwelling catheters) 2 way, Tiemann (coudé) tip, 15cc, 2-way, 41cm. Latex free	12 FR	AA6312	5
	14 FR	AA6314	5
	16 FR	AA6316	5
	18 FR	AA6318	5
	20 FR	AA6320	5
	22 FR	AA6322	5
	24 FR	AA6324	5
Folysil (indwelling catheters) Paediatric, 2 way, Tiemann (coudé) tip, 3cc, 30-35cm. Latex free	8 FR, 30cm	AA6308	5
	10 FR, 35cm	AA6310	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

Continence Care

X-Flow

Prostatic Catheters

Performance and safety

Wide irrigation and drainage channels to prevent blocking. Large bevels and eyes to ensure excellent drainage capacity. Stiff walls to facilitate evacuation.

All silicone, to avoid any risk of allergy or toxicity. Clear catheter body to monitor drainage quality.


Closed round tip for ease of placement of an introducer if needed

	Size	Product Code	Units
X-Flow (Prostatic Catheter) 3-way, straight, 50cc, 42 cm (Latex free)	18 FR	AB6018	5
	20 FR	AB6020	5
	22 FR	AB6022	5
	24 FR	AB6024	5
X-Flow (Prostatic Catheter) 3-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6118	5
	20 FR	AB6120	5
	22 FR	AB6122	5
	24 FR	AB6124	5
X-Flow (Prostatic Catheter) 3-way, Delinotte (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6218	5
	20 FR	AB6220	5
	22 FR	AB6222	5
	24 FR	AB6224	5
X-Flow (Prostatic Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm (Latex free) Curved tip allows easy introduction without stylet.	18 FR	AB6318	5
	20 FR	AB6320	5
	22 FR	AB6322	5
	24 FR	AB6324	5
X-Flow (Prostatic Catheter) 2-way, Couvelaire (whistle tip), 50cc, 42 cm (Latex free)	18 FR	AB6518	5
	20 FR	AB6520	5
	22 FR	AB6522	5
	24 FR	AB6524	5
X-Flow (Haematuria Catheter) 3-way, Dufour (coudé tip), 50cc, 42 cm) Curved tip allows easy introduction without stylet. Same characteristics as X-Flow prostatic catheter (AB63XX) with 30 cc balloon	18 FR	AB6A18	5
	20 FR	AB6A20	5
	22 FR	AB6A22	5
	24 FR	AB6A24	5

For information on other surgical urology products, please visit www.urologycare.coloplast.ca

TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Wound Care




Product Selection Guide

For Medium-High Exuding Wounds

Managing Exuding Wounds?


Exuding Wounds



Biatain®

- Superior absorption and exudate management - also under compression
- Minimal risk of maceration and leakage
- Soft & conformable foam
- Wear-time is up to 7 days

Exuding Wounds with Slough and/or Cavity



Biatain® Alginate (formerly SeaSorb® Soft)


- Highly absorbent alginate dressing
- Highly cohesive gel
- Easy one-piece removal
- Wear-time is up to 7 days

+

Biatain Silicone or **Biatain Silicone Lite**
as cover dressing

Managing Exuding Wounds and Infection?


Infected Wounds



Biatain® Ag

- Superior absorption for infected wounds
- Sustained release of silver during entire wear time (up to 7 days)
- Effective on broad range of bacteria

Infected Wounds with Slough and/or Cavity




Biatain® Alginate Ag (formerly SeaSorb® Ag)

- Highly absorbent antimicrobial alginate dressing
- Effective on broad range of bacteria
- Easy one-piece removal
- Wear-time is up to 7 days

+

Biatain Silicone or **Biatain Silicone Lite**
as cover dressing

Managing Exuding Wounds and Pain?



Biatain® Ibu

- Superior absorption and exudate management - also under compression
- Combines moist wound healing with local release of ibuprofen
- Wear-time is up to 7 days

All Coloplast Wound Care products are LATEX-FREE
TELEPHONE 1-888-880-8605 / 905-820-2731
FAX 1-877-820-8206


Wound Care

Biatain®

Superior absorption for faster healing

Biatain foam dressings are proven to provide faster healing of exuding wounds:

- Superior exudate management
- Minimal risk of leakage and maceration
- Long wear time, up to 7 days
- A cost-effective solution
- Packaged in sterile single use packs






Areas of Use:

- Leg/Vascular ulcers
- Pressure ulcers
- Diabetic foot ulcers
- Traumatic wounds
- Wounds with compromised peri-wound skin

The unique 3-D structure of Biatain serves as a partial fluid lock, retaining the exudate - even under compression therapy. Wound exudate is absorbed locally with minimal lateral dispersal providing an optimal moist wound healing environment.

Foam Dressing

		Size	Product Code	Units
	Biatain® Non-Adhesive Foam Dressing <ul style="list-style-type: none"> • Superior exudate management • Suitable for use on fragile skin due to absence of adhesive 	2" x 2.8" (5 x 7 cm)	6105	10
		4" x 4" (10 x 10 cm)	3410	10
		4" x 8" (10 x 20 cm)	3412	5
		6" x 6" (15 x 15 cm)	3413	5
		8" x 8" (20 x 20 cm)	3416	5
	Biatain® Adhesive Foam Dressing <ul style="list-style-type: none"> • Adhesive border based on hydrocolloid technology 	3" x 3" (7.5 x 7.5 cm)	3462	10
		4" x 4" (10 x 10 cm)	3430	10
		5" x 5" (12.5 x 12.5 cm)	3420	10
		7" x 7" (18 x 18 cm)	3423	5
	Biatain® Adhesive Foam Dressing/Shapes <ul style="list-style-type: none"> • Specifically designed to conform to difficult to dress areas of the body • The adhesive border is based on hydrocolloid technology 	Sacral 7" x 7" (17 x 17 cm)	3483	5
		Sacral 9" x 9" (23 x 23 cm)	3485	5
		Heel 7 ½" x 8" (19 x 20 cm)	3488	5

All Coloplast Wound Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Wound Care

Biatain® Soft-Hold


Superior absorption for faster healing

Biatain foam dressings are proven to provide faster healing of exuding wounds:

- Superior exudate management
- Minimal risk of leakage and maceration
- Long wear time, up to 7 days
- A cost-effective solution
- Packaged in sterile single use packs

Areas of Use:

- Leg/Vascular ulcers
- Pressure ulcers
- Diabetic foot ulcers
- Traumatic wounds
- Wounds with compromised peri-wound skin



Unique 3-D alveolar structure of Biatain foam serves as a partial fluid lock, retaining exudate even under compression. The soft adherent layer minimizes trauma to the wound and pain to the patient during dressing changes.

Foam Dressing

Foam Dressing		Size	Product Code	Units
	Biatain® Soft-Hold Foam Dressing <ul style="list-style-type: none"> • Soft adherent layer which allows for: <ul style="list-style-type: none"> - Easy Application - Easy Removal - Skin Protection 	4" x 4" (10 x 10 cm)	3470	5
		4" x 8" (10 x 20 cm)	3472	5
		6" x 6" (15 x 15 cm)	3475	5

Wound Care

Biatain® Silicone

Superior absorption for faster healing

Biatain foam dressings are proven to provide faster healing of exuding wounds:

- Superior exudate management
- Minimal risk of leakage and maceration
- Long wear time, up to 7 days
- A cost-effective solution
- Packaged in sterile single use packs


A versatile solution for a wide variety of non-infected moist wounds:

- Leg ulcers
- Second degree burns
- Post operative wounds
- Non-infected diabetic foot ulcers
- Donor sites
- Pressure ulcers*
- Skin abrasions



Superior absorption under pressure makes Biatain Silicone the obvious choice in combination with compression therapy.

Foam Dressing with Silicone Border

		Size	Product Code	Units
 <p>Biatain® Silicone Foam Dressing</p> <ul style="list-style-type: none"> • Superior absorption and exudate management allows for reduced maceration and leakage • Unique foam structure ensures high retention of exudate - also under compression • Soft & conformable foam allows for high patient comfort • Long wear time - up to 7 days • Cost effective solution for exuding wounds 		3" x 3" (7.5 x 7.5 cm)	33434	10
		4" x 4" (10 x 10 cm)	33435	10
		5" x 5" (12.5 x 12.5 cm)	33436	10
		6" x 6" (15 x 15 cm)	33437	5
		7" x 7" (17.5 x 17.5 cm)	33438	5

*Except pressure ulcers with a high degree of friction and shear force. In these cases, Coloplast recommends the use of Biatain Adhesive dressings.

Biatain® Silicone Lite

Same effective healing, more mobility


- Thinner foam for increased mobility
- Ideal for low exudate, non-infected wounds
- Easy, painless removal and repositioning
- Silicone adhesive border

Areas of Use:

- Leg ulcers
- Second degree burns
- Post operative wounds
- Non-infected diabetic foot ulcers
- Donor sites
- Pressure ulcers*
- Skin abrasions



Foam Dressing with Silicone Border

		Size	Product Code	Units
	<p>Biatain® Silicone Lite – same effective healing, more mobility</p> <ul style="list-style-type: none"> • Thinner foam for increased mobility • Ideal for low exudate, non-infected wounds • Easy, painless removal and repositioning • Silicone adhesive border 	3" x 3" (7.5 x 7.5 cm)	33444	10
		4" x 4" (10 x 10 cm)	33445	10
		5" x 5" (12.5 x 12.5 cm)	33446	10


*Except pressure ulcers with a high degree of friction and shear force. In these cases, Coloplast recommends the use of Biatain Adhesive dressings.

Wound Care

Biatain® Ag

Superior absorption for infected wounds

Biatain Ag foam dressings combine superior exudate management with sustained antimicrobial activity for improved patient outcomes. Sustained release of silver assures fast and effective kill against bacterial strains known to be detrimental to wound healing such as *Pseudomonas aeruginosa*, *Staphylococcus aureus* and *B-haemolytic Streptococcus*.




Areas of Use:

- Exudating wounds with delayed healing due to bacteria
- Exudating wounds where there is a risk of infection
- Leg/Vascular ulcers
- Pressure ulcers
- Diabetic foot ulcers
- Second degree burns
- Donor sites
- Postoperative wounds
- Skin abrasions with delayed healing due to bacteria
- Suitable for use in combination with compression therapy

Ionic silver is evenly distributed in the Biatain Ag foam. In the presence of exudate, there is a sustained release of silver into the wound bed throughout the wear time of the dressing. Silver is continuously released for up to 7 days in proportion to the level of exudate within the wound. This ensures an effective control of the local infection which is maintained during the entire wear time.

Antimicrobial Foam Dressing

		Size	Product Code	Units
	Biatain® Ag Non-Adhesive Foam Dressing <ul style="list-style-type: none"> • Superior exudate management • Suitable for use on fragile skin due to absence of adhesive 	2" x 2.8" (5 x 7 cm)	5105	5
		4" x 4" (10 x 10 cm)	9622	5
		4" x 8" (10 x 20 cm)	9623	5
		6" x 6" (15 x 15 cm)	9625	5
		8" x 8" (20 x 20 cm)	9626	5
	Biatain® Ag Adhesive Foam Dressing <ul style="list-style-type: none"> • Adhesive border based on hydrocolloid technology 	3" x 3" (7.5 x 7.5 cm)	9631	5
		5" x 5" (12.5 x 12.5 cm)	9632	5
		7" x 7" (18 x 18 cm)	9635	5
	Biatain® Ag Adhesive Foam Dressing/Shapes <ul style="list-style-type: none"> • Specifically designed to conform to difficult to dress areas of the body • The adhesive border is based on hydrocolloid technology 	Sacral 9" x 9" (23 x 23 cm)	9641	5
		Heel 7 1/2" x 8" (19 x 20 cm)	9643	5

All Coloplast Wound Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Biatain® Silicone Ag

Superior absorption for infected wounds

Biatain Silicone Ag foam dressings combine superior exudate management with sustained antimicrobial activity for improved patient outcomes.

Sustained release of silver assures fast and effective kill against bacterial strains known to be detrimental to wound healing such as *Pseudomonas aeruginosa*, *Staphylococcus aureus* and *B-haemolytic Streptococcus*.


Areas of Use:

- Exudating wounds with delayed healing due to bacteria
- Exudating wounds where there is a risk of infection
- Leg/Vascular ulcers
- Pressure ulcers*
- Diabetic foot ulcers
- Second degree burns
- Donor sites
- Postoperative wounds
- Skin abrasions with delayed healing due to bacteria
- Suitable for use in combination with compression therapy



Ionic silver is evenly distributed in the Biatain Ag foam. In the presence of exudate, there is a sustained release of silver into the wound bed throughout the wear time of the dressing. Silver is continuously released for up to 7 days in proportion to the level of exudate within the wound. This ensures an effective control of the local infection which is maintained during the entire wear time.

Antimicrobial Foam Dressing

	Size	Product Code	Units
 <p>Biatain® Silicone Ag</p> <ul style="list-style-type: none"> • Superior exudate management • Soft silicone adhesive 	3" x 3" (7.5 x 7.5 cm)	39636	5
	4" x 4" (10 x 10 cm)	39637	5
	5" x 5" (12.5 x 12.5)	39638	5


*Except pressure ulcers with a high degree of friction and shear force. In these cases, Coloplast recommends the use of Biatain Adhesive dressings.

Wound Care

Biatain® Ibu

Superior absorption for painful wounds

Biatain Ibu combines superior exudates management with the local (nonsystemic) continuous release of ibuprofen for improved patient outcomes.





Areas of Use:

- Painful exuding wounds
- Leg/Vascular ulcers
- Pressure ulcers
- Non-infected diabetic foot ulcers
- Second degree burns
- Donor sites (not applicable for Soft-Hold)
- Postoperative wounds
- Skin abrasions
- Suitable for use in combination with compression therapy

Biatain Ibu combines superior exudate management and continuous release of ibuprofen. Biatain Ibu releases ibuprofen locally with no observed systemic effect. In the presence of exudate, there is a sustained release of a low dose of ibuprofen into the wound bed throughout the wear time of the dressing. One dressing (10 x 10 cm) contains 50 mg of ibuprofen which is continuously released for up to 7 days in proportion to the level of exudate.

Foam Dressing with Ibu


		Size	Product Code	Units
	Biatain® Ibu Non-Adhesive Foam Dressing <ul style="list-style-type: none"> • Superior Exudate Management • Suitable for use on fragile skin due to absence of adhesive • Well tolerated with no observed systemic effect 	4" x 4" (10 x 10 cm)	4110	5
		4" x 8" (10 x 20 cm)	4112	5
		6" x 6" (15 x 15 cm)	4115	5
		8" x 8" (20 x 20 cm)	4120	5
	Biatain® Ibu Soft-Hold Foam Dressing <ul style="list-style-type: none"> • Gentle adherence minimizes trauma and pain at removal • Sustained release of ibuprofen reduces persistent pain during entire wear time 	4" x 4" (10 x 10 cm)	4140	5
		4" x 8" (10 x 20 cm)	4142	5

Wound Care

Biatain® Alginate (Formerly SeaSorb®)

Superior absorption for wound with slough and/or cavity



Allows for intact removal, will not stick to healing tissue.
 Absorbs, collects and contains exudate while providing a moist healing environment.
 Forms a very soft, cohesive and hydrophilic gel. Has haemostatic properties



Areas of Use:

- Painful exuding wounds
- Moderate to heavily exuding wounds
- Cavity wounds where it should be gently layered into the wound, as the gel will conform to the wound shape
- SeaSorb Ag is indicated for wounds with signs of local infection or critical colonization

Calcium Alginate Dressing

		Size	Product Code	Units
	Biatain® Alginate Ag (Formerly SeaSorb® Ag) <ul style="list-style-type: none"> • Provides an effective, sustained release of silver against common bacteria for up to 4 days • Reduces odour - Effectively controls bacteria which cause odour • Cost effective - High absorbency results in longer wear time and reduced nursing time 	4" x 4" (10 x 10 cm)	3760	10
		1" x 17½" (44 cm) rope	3780	10
	Biatain® Alginate (Formerly Seasorb® Soft) <ul style="list-style-type: none"> • Forms a very soft cohesive and hydrophilic gel on contact with wound exudate • Ideal for moderate to highly exuding wounds • Easy, one piece removal minimizes pain at dressing change 	2" x 2" (5 x 5 cm)	3705	30
		4" x 4" (10 x 10 cm)	3710	10
		6" x 6" (15 x 15 cm)	3715	10
		1" x 17½" (44 cm) rope	3740	6

New size

Wound Care

Comfeel® Plus

Hydrocolloid Dressing

Comfeel Plus hydrocolloid dressings change colour as they absorb exudate. This helps determine the optimal time for dressing change. As exudate is absorbed, a moist gel forms, and the dressing begins to turn white. Maximum absorbency has been reached when the dressing becomes opaque or when exudate extends to within ½" of the edges.

Polyurethane top film allows excess moisture to evaporate while maintaining a moist environment for optimal wound healing.

Areas of Use:

- Management of low to moderately exudating leg and pressure ulcers
- Superficial burns
- Post-operative wounds
- Superficial partial-thickness burns
- Skin abrasions



Hydrocolloid Dressing

		Size	Product Code	Units
	Comfeel® Plus Ulcer Dressing <ul style="list-style-type: none"> • Combines a highly absorbent sterile hydrocolloid with an alginate for superior exudate management, for use up to 7 days • Grid imprint on the polyurethane top film facilitates wound measurement • Tapered edges minimize edge roll 	1½" x 2½" (4 x 6 cm)	3146	30
		4" x 4" (10 x 10 cm)	3110	10
		6" x 6" (15 x 15 cm)	3115	5
		8" x 8" (20 x 20 cm)	3120	5
	Comfeel® Plus Sacral (Triangle) <ul style="list-style-type: none"> • Specially designed for use in the sacral area • Thin, tapered edges increase wear time by preventing edge roll • Grid imprint allows convenient, accurate measurement of the wound area 	7" x 8" (18 x 20 cm)	3285	5
	Comfeel® Plus Transparent Dressing (Formerly Comfeel Plus Clear Dressing) <ul style="list-style-type: none"> • Transparency allows for wound inspection without dressing change • Ideal for friction management • Low profile design minimizes edge roll 	2" x 2¾" (5 x 7 cm)	3530	10
		4" x 4" (10 x 10 cm)	3533	10
		3½" x 5½" (9 x 14 cm)	3536	10
		6" x 6" (15 x 15 cm)	3539	5


All Coloplast Wound Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Wound Care


Hydrogels

		Size	Product Code	Units
	Purilon® Gel (Sterile) <ul style="list-style-type: none"> • High viscosity hydrogel maintains moist wound environment, while assisting in autolytic debridement • Ideal for necrotic and sloughy wounds • Accordion applicator allows delivery of a precise amount of gel for easy one-handed application 	0.28 oz. / 8 g	3906	10
		0.5 oz. / 15 g	3900	10
	Woun'Dres® Collagen Hydrogel <ul style="list-style-type: none"> • Promotes autolytic debridement by rehydrating and softening dry wounds and necrotic tissue • Collagen is the main component of skin and connective tissue and is known to play an important role in all phases of wound healing 	1 oz. / 28 g tube	1166	36
		3 oz. / 84 g tube	7690	12

Hydrophilic Paste

		Size	Product Code	Units
	Triad™ Hydrophilic Wound Dressing (Sterile) <ul style="list-style-type: none"> • Absorbs excess wound exudate while retaining a moist wound environment • Zinc oxide-based hydrophilic paste gently adheres to moist, weeping wounds and spreads evenly over entire wound area • Protects peri-wound skin 	2.5 oz. / 71 g tube	1964	12
		6 oz. / 170 g tube	1967	12

Wound Cleanser

		Size	Product Code	Units
	Sea-Cleans® Wound Cleanser <ul style="list-style-type: none"> • Cost-effective, saline-based wound cleanser designed to loosen debris from the wound bed without interfering with the healing process • Saline-based cleanser provides safe and effective irrigation • Soften eschar to facilitate debridement • Non-toxic, non-ionic, pH-balanced formula • DIN 02230331 	6 oz. / 178 ml	1063	12
		12 oz./ 355 ml	1061	12

All Coloplast Wound Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Application and Cutting Guide

Application

Biatain, Biatain Ag & Biatain Ibu Non-Adhesive Foam Dressings

The dressing must be applied with the plain (non-printed) side towards the wound.

Biatain & Biatain Ibu Soft-Hold Foam Dressings

Remove the protective film and apply the dressing with the plain (non-printed) side towards the wound.

Biatain & Biatain Ag Adhesive Foam Dressings

Use the handles to prevent touching the dressing and to ensure aseptic application. Remove the protective film first and place the foam part of the dressing over the wound. Then remove the release films on the sides as part of the dressing application.

Biatain & Biatain Ag Adhesive Foam

Dressings/Sacral

Use the handles to ensure aseptic application. Remove the center protective film. Place the narrow end of the dressing as far down over the wound as possible and fix the dressing working up and outwards. Remove the handles.

Biatain & Biatain Ag Adhesive Foam Dressings/Heel

Use the handles to ensure aseptic application. Remove the center protective film. The dressing is shaped like an "arrow". Fold the dressing in a 90° angle between the "arrow head" and "arrow tail". Let the "arrow" point away from the heel tip and fix the "arrow tail" first. Secondly, fix the "arrow head". Remove the handles one by one and fix the sides carefully, in order that the foam parts touch or overlap.

Removal

To remove the dressing, gently lift the corners of the dressing away from the wound.

To remove the adhesive dressings, it is recommended to loosen the adhesive border before lifting the dressing away from the wound.

Biatain Adhesive Foam Dressing/Sacral should be removed from the top edge and down towards the anus to minimize the risk of transmitting infection.

General Notes:

Disclaimer: These are general guidelines, please check local treatment recommendations applicable in your country or healthcare institution.

- Always check the cause of the wound.
- Follow your local infection control guidelines to prevent cross infection.
- Please make sure to use a pair of sharp clean scissors.
- Always remember to cut across the peel back, otherwise the backing paper will be difficult to remove.
- When cutting dressings always round off the edges for better adhesion.
- For adhesive dressings keep the backing film as a template for next application on your patient. For non-adhesive dressings draw a template on the pouch.

- Most manufacturers recommend placing the dressing over the centre of the wound. Consider the effect of gravity and apply the dressing to account for the direction of fluid drainage.
- Do not apply dressings with tension as this may cause blistering.
- Although some dressings might provide padding and reduction in friction, they do not relieve pressure.
- For all joints, the principles of application are similar. Examples given here on heels, elbows and knees may be used on any of these joints.

Developed by: Pia Carlsen, RN, Denmark, Jacqui Fletcher, Principal Lecturer, MSc BSc (Hons) PG Dip (ED) RN ILT, UK, Maria Mousley, AHP, sultant, Podiatrist, UK, Henri Post, M.A., Nurse Practitioner, RN, Netherlands, Adrienne Taylor, MSc, RN, RM, DN, UK, Wilma Verbeek-Gijsbers, RN, Netherlands, Robert Zimmer, RN, Germany

All Coloplast Wound Care products are LATEX-FREE
TELEPHONE 1-888-880-8605 / 905-820-2731
FAX 1-877-820-8206

Application and Cutting Guide



Washing hands

Always remember to wash your hands before and after the application of a dressing.

Please use your recommended dressing technique according to infection control guidelines.

Disclaimer: These are general guidelines, please check local treatment recommendations applicable in your region or healthcare institution.



Ears



Use a thin hydrocolloid (5 x 7 cm), cut into a strip, fold in half lengthways, then cut a 'fringe' along one edge (the cuts should be approximately 0.5 cm apart). Apply the uncut edge of the dressing first, and use the cut edges to fold around the shape of the ear. If the wound is wet, use a small piece of alginate underneath.



Ears



For wounds caused in the area where oxygen tubing or foam rubs on the ears, cut half way into a thin hydrocolloid (the piece should be approx (2 x 3 cm), anchor to the flat of the head with the uncut half, and then put the other two sections around the ear. These can also be adapted for the nasal region.



Ears



Use a Non-adhesive or a Soft-Hold dressing for this kind of wound. Cut the dressing in order to fix the size of the ear. Choose a suitable secondary dressing to keep the dressing in place.



Elbows



Cut two small slits in the dressing and apply it, mould the slit around the elbow. Remember to apply the dressing with the patient bending the arm in a 90° angle.



Elbows



Use 3 strips of thin hydrocolloid. With the joints at an angle of 90° apply the first strip directly over the angle of the joint. Apply the second strip overlapping one edge approximately 30° to the first strip. Apply the third strip over the other edge of the first strip in the same way. Cup the dressing firmly to mould the shape of the joint. If the wound is wet, a foam dressing is more suitable.



Bunion



Use for the hallux toe and for inter-digital ulcers, secure with a secondary dressing.



Toes



Cut a dressing and use it for interdigital wounds. The dressing must be secured with a secondary dressing.



Toes



Use for the tip of the toes and for inter-digital ulcers, secure with tubular gauze.



Toes



Cut a "bow-tie" shaped thin hydrocolloid dressing. Use the slimmest dressing possible considering footwear. Often a non-adhesive dressing with separate fixation will work better than an adhesive dressing due to the pressures applied to the feet on walking.



Toes



For wet wounds use a foam dressing. Additional fixation is required.



Hands



For margin of hand, cut one slit towards the edge of the dressing. Apply over the first or fifth digit (as appropriate). To anchor the dressing mould the remaining dressing to the hand. Where there is fungal infection present, a silver-containing foam dressing is appropriate.



Fingers



For wounds between the fingers cut a dressing and slide onto the finger. Use a secondary dressing to secure the foam dressing.



Fingers



For small wounds on and between fingers a thin hydrocolloid dressing can easily be cut into small strips to fit the wound.



Fingers



For multiple digits.



Heels



There is a variety of different heel dressings available, consider which of these shapes is best for your patient.



Axilla



A heel dressing may be a good solution for dressing the axilla area. Keep the direction of fluid in mind and apply the largest part of the dressing as illustrated on the photo.



Sternum



Particularly suitable for women. Cut according to the size of breast.



Peristoma wounds



When the wound is close to an ostomy, care needs to be taken when deciding whether to apply the stoma base plate on top or underneath the dressing.



Sacral



Bottoms are considerable in size shape and skin texture. General principles for management include: Use the slimmest dressing possible. There is a variety of different sacral dressings available, consider which of these shapes is best for your patient's bottom. Please ensure that you differentiate between pressure ulcers and moisture lesions, as the supporting management differs greatly.



Sacral



If the wound is wet, use e.g. a foam and/ or an alginate dressing (the alginate dressing requires a secondary dressing).

Skin Care



Simplify Skin Care

ALL PATIENTS

Step 1: Gently Cleanse Using a Mild pH Balanced Cleanser



Sproam® Cleanser

All Body Cleansing

- » Add 2-3 pumps onto a warm moist washcloth and cleanse
- » Gently pat dry

Shampooing

- » Dampen hair, add 2-3 pumps directly on to dampened hair
- » Gently massage into hair and scalp for 30-60 seconds
- » Rub gently with a moist washcloth to remove dirt and oils

Incontinence Cleanup

- » Add 2-3 pumps on a warm moist washcloth, and cleanse
- » If feces, spray 2-3 pumps directly onto soiled area, wait 30 seconds and wipe off using a warm moist washcloth

Step 2: Moisturize All Patients/Residents



Sween® Lotion Unscented

- » Moisturizing body lotion with natural vitamin E
- » For all dry skin areas
- » Gently apply a small amount of lotion to patient's skin

INCONTINENCE

Step 3: Minimize Skin Exposure To Moisture From Incontinence



Critic-Aid® Clear

- » Gently clean the perineal area and apply a thin layer of **Critic-Aid®** Clear evenly over the entire area to be protected

MOISTURE CONTROL

Moisture Control in Skin Folds

InterDry® - Textile with Antimicrobial Silver Complex


- » Wicks away moisture; keeping skin dry
- » Helps relieve maceration, denudement, inflammation, itching, erythema and satellite lesions
- » Provides effective antimicrobial action for up to five days
- » Reduces odour commonly associated with skin folds
- » Place **InterDry®** between skin folds; leaving at least two inches of the textile outside the skin fold



Bathing and Hand Washing Products

		Size	Product Code	Units
	Gentle Rain® Extra Mild Cleanser Sensitive Skin Cleanser, Shampoo and Handwash Contains special ingredients and extra emollients to gently clean and moisturize fragile, sensitive skin. Mild enough for daily use. Dermatologist-tested, hypoallergenic and dye-free. pH balanced. Safe to use in bath tubs.	5 ml sachet	7231	300
		50 ml	7230	36
		600 ml	7233	12
		1000 ml	7234	12
		3.8 L	7236	4
	Gentle Rain® Antibacterial Cleanser Triclosan Antiseptic Lotion Skin Cleanser Mild cleanser with triclosan provides highly effective control against VRE. DIN 02308673	600 ml	7205	12

Antiseptic Hand Cleansers

		Size	Product Code	Units
	Isagel® Hand Sanitizing Gel Ethyl Alcohol Antiseptic No-Rinse Hand Gel 70% Ethyl Alcohol Effective, fast-acting handrub kills greater than 99.9% of most transient microorganisms (including MRSA and VRE) in 15 seconds or less. Enriched with emollients, Isagel is gentle on the skin. Compatible with gloves and CHG-based products.	115 ml	1644	36
		621 ml	1645	12
		800 ml	7041	12
		3.8 L	1647 Minimum order of 5 cases	4

Skin Care

No-Rinse Full Body and Incontinence Cleansers

		Size	Product Code	Units
	Bedside-Care® Sensitive Skin (Formerly Bedside-Care Unscented Foam) No-Rinse Full Body Shampoo and Perineal Cleanser Unique formulation provides a thicker, richer foaming action than traditional bathing or perineal cleansers. No more runny messes. Cleanses, moisturizes and conditions skin. No rinsing needed. Hypoallergenic formula is alcohol-free and fragrance-free. pH-balanced	125 ml	7147	36
		240 ml	7146	12
	Bedside-Care® Scented Foam No-Rinse Full Body Shampoo and Perineal Cleanser Same hypoallergenic formula, alcohol-free. pH-balanced	240 ml	7145	12
	Sproam® Cleanser Benzethonium Chloride Antiseptic No-Rinse All Body Spray/Foam Cleanser Complete head-to-toe cleansing, shampooing, bathing and incontinent skin care. Can be used as a spray or foam, for maximum versatility. Non-aerosol dense foaming action. Dye and fragrance free. pH-balanced DIN 02245903	175 ml	0897	12
		350 ml	0901	12
		3.8 L	0903	4
	Bedside-Care® Cleanser No-Rinse All-Body Wash and Incontinent Cleanser A complete bathing and incontinent cleanser in a bottle. Gently cleanses while moisturizing irritated skin. Soothing and deodorizing. pH-balanced	120 ml	1768	36
		250 ml	1772	12
		3.8 L	1774	4
	Peri Wash® II Cleanser Benzethonium Chloride Antiseptic No-Rinse Cleanser Remoisturizes the skin as it cleans for gentle continued comfort. No dyes or fragrances. pH-balanced DIN 02229501	250 ml	1485	12
		3.8 L	1487*	4

* To be discontinued

All Coloplast Skin Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Skin Care

Moisturizers


		Size	Product Code	Units
	<p>Sween® 24 Cream Dimethicone Cream Moisturizing Skin Protectant</p> <p>Once-a-day application for long-lasting 24 hour protection. Formulated for sensitive skin. Softens and soothes skin by sealing in moisture. Applies easily and absorbs quickly. DIN 02250322</p>	4 g sachet	7090	300
		90 g	7094	12
	<p>Atrac-Tain® Cream With Urea and Alpha Hydroxy Acid (AHA)</p> <p>Advanced moisturizing relief for severely afflicted dry, scaly, cracked skin.</p> <p>Ideal for heels, elbows, knees, lower leg areas. Contains 10% Urea and 4% AHA.</p>	2 g sachet	1843	300
		60 ml	1802	12
		140 ml	1814	12
	<p>Sween® Unscented Lotion (formerly Xtra Care Unscented Lotion) With Natural Vitamin E</p> <p>An everyday, all-over body lotion for normal maintenance of healthy skin.</p>	60 ml	7222	36
		120 ml	7223	36
		620 ml	7224	12
	<p>Sween Cream® With Natural Vitamins A & D</p> <p>Moisturizing cream designed to hydrate moderately dry skin.</p> <p>Non-occlusive formula absorbs quickly into the skin and allows adhesion of tapes and appliances.</p>	85 g	7067	12
		190 g	7068	12

All Coloplast Skin Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Incontinence Clean up And Protection (Barriers)

		Size	Product Code	Units
	<p>Baza Cleanse & Protect® Lotion Dimethicone Skin Protectant Lotion</p> <p>One-step, no-rinse lotion that cleanses, moisturizes and helps prevent and protect red, injured, chafed and chapped skin due to incontinence or diaper rash. Combines mild cleansing agents, a breathable skin barrier, and moisturizing ingredients. Safe to use with incontinence briefs. Fragrance free.</p> <p>Also available with a mild fragrance to reduce incontinent odours.</p> <p>pH-balanced DIN 02243966</p>	4 g (fragrance free)	7710	300
		240 ml (fragrance free)	7712	12
		240 ml (with fragrance)	7725	12
	<p>Baza® Protect II Cream (formerly Baza® Protect cream) Zinc Oxide and Dimethicone Skin Protectant Cream</p> <p>Special formulation of ingredients that aids in the protection, relief and/or treatment of diaper dermatitis (rash) from incontinence. Safe to use with incontinence briefs.</p> <p>DIN 02315793</p>	60 g	7086 (formerly 1877)	12
		140 g	7087 (formerly 1880)	12
	<p>Critic-Aid® Clear Skin Protectant Clear Moisture Barrier Ointment</p> <p>Clear moisture barrier ointment adheres to intact or denuded skin and manages moisture. Easy to apply and remove. Provides protection, relief and/or treatment of diaper dermatitis (rash) caused by light to severe or chronic diarrhea and enzymatic drainage. Safe to use with incontinence briefs.</p> <p>DIN 02277778</p>	4 g sachet	7565	300
		71 g	7573	12

Skin Fold Management

		Size	Product Code	Units
	<p>InterDry® Textile Textile with Silver Complex</p> <p>Manages inflammation between skin folds. Wicks away moisture from the skin. Reduces skin-to-skin friction.</p> <p>Provides odour control.</p> <p>Multi-use roll can be easily cut and applied to fit individual body location.</p>	Box of 1 roll: 10 X 144 in (25 X 366 cm)	7910	1
		10 X 36 in (25 X 92 cm)	7912	10




All Coloplast Skin Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Skin Care

Other

Other		Size	Product Code	Units
	Sween® Body Powder With natural corn starch powder to relieve the discomfort of chafed skin and minor skin irritations. Absorbs moisture in skin creases without caking. Long-lasting deodorant action.	85 g	0506	36
		227 g	0505	36
	Hex-On® Odour antagonist Eliminates odour. Non-aerosol.	Fresh Linen	2 oz. / 59 ml 7583	12

Dispensers

Dispensers		Product Code	Units
	Wall Dispenser for Bag Safe and compact, easy to dispense. Designed for 800 ml and 1000 ml sized bags. (7234 and 7041)	7251	12
	Wall Dispenser Drip Tray Accessory to be use with 7251 to prevent spillage.	7262	1
	Dispenser Bracket Made of very sturdy plastic. Made specifically to hold 621 ml sized pump bottles (7224, 1645, 7205 and 7233). Theft Resistant.	0089	12

DIN items cannot be returned, please contact customer service for details

All Coloplast Skin Care products are LATEX-FREE
 TELEPHONE 1-888-880-8605 / 905-820-2731
 FAX 1-877-820-8206

Discontinued Products

Ostomy Care


Item Code	Suggested Replacement
2156	2151
2161	2166
2837	2836
11801	11802
12345	12375
12471	12472
12801	NA
12803	NA
12834	12835
12848	10031
14257	2833
14297	14283
14298	14283
14314	14329/14319
14326	14319
14327	14332
14334	14329
14421	14414
14434	14445
14435	14445
14442	14445
14444	14445
14645	14402

Continence Care

Item Code	Suggested Replacement
5040H	5240H
3816	1116

Skin Care

Item Code	Suggested Replacement
283	7090/7067
1676	1644
1806	7223
7662	7067
7584	7583



Visit www.improvingcare.coloplast.ca for an updated list of discontinued products or call our Customer Service Department at 1-888-880-8605

Product Index

Product Code	Page
89	79
110	48
112	48
114	48
116	48
208	49
210	49
212	49
214	49
240	49
305	49
306	49
308	49
310	49
408	48
410	48
412	48
414	48
416	48
418	48
421	37
422	37
450	48
460	48
475	49
505	79
506	79
509	38
608	48
610	48
612	48
614	48
806	48
808	48
810	48
812	48
814	48
816	48
818	48
897	76
901	76
903	76
925	38
1003	40
1061	70

Product Code	Page
1063	70
1108	50
1110	40
1112	50
1114	50
1116	50
1166	70
1225	38
1450	57
1451	57
1485	76
1487	76
1500	40
1644	75
1645	75
1647	75
1676	80
1768	76
1772	76
1774	76
1802	77
1806	80
1814	77
1843	77
1907	37
1964	70
1967	70
2041	38
2115	34
2125	34
2151	34
2156	80
2161	80
2166	34
2181	34
2320	38
2330	38
2340	38
2350	38
2490	31
2501	39
2520	31
2614	50
2650	38
2655	36

Product Code	Page
2792	23
2793	23
2801	39
2802	39
2803	39
2807	39
2808	39
2831	20
2832	20
2833	20
2836	42
2837	80
2881	19
2882	19
2883	19
3110	69
3115	69
3120	69
3146	69
3210	36
3215	36
3220	36
3285	69
3410	61
3412	61
3413	61
3416	61
3420	61
3423	61
3430	61
3462	61
3470	62
3472	62
3475	62
3483	61
3485	61
3488	61
3530	69
3533	69
3536	69
3539	69
3705	68
3710	68
3715	68
3740	68

Product Code	Page
3760	68
3780	68
3816	80
3900	70
3906	70
4110	67
4112	67
4115	67
4120	67
4140	67
4142	67
4410	49
4412	49
4414	49
5030	51
5035	51
5062	55
5063	55
5070	57
5100	52
5105	65
5110	53
5125	53
5130	53
5135	53
5161	55
5167	55
5170	55
5174	55
5190	52
5200	53
5205	53
5210	53
5212	53
5215	53
5221	51
5225	51
5230	51
5235	51
5240	51
5290	52
5300	52
5390	52
5400	52
5490	52

Product Code	Page
5500	52
5580	32
5585	32
5590	52
6105	61
6200	52
6300	52
6400	52
6500	52
7041	75
7067	77
7068	77
7078	56
7086	78
7087	78
7088	56
7090	56
7090	77
7092	56
7094	77
7096	56
7145	76
7146	76
7147	76
7205	75
7222	77
7223	77
7224	77
7230	75
7231	75
7233	75
7234	75
7236	75
7251	79
7262	79
7565	78
7573	78
7583	79
7662	80
7690	70
7710	78
7712	78
7725	78
7910	78
7912	78

Product Index

Product Code	Page
8000	54
8009	34
8100	54
8200	54
8205	54
8300	54
8305	54
8400	54
8500	54
9005	56
9007	56
9500	38
9505	38
9622	65
9623	65
9625	65
9626	65
9631	65
9632	65
9635	65
9641	65
9643	65
10011	4
10012	4
10013	4
10015	5
10016	5
10017	5
10021	4
10022	4
10023	4
10025	5
10026	5
10027	5
10031	4
10032	4
10035	5
10036	5
10041	4
10045	5
10101	9
10102	9
10103	9
10104	10
10105	10

Product Code	Page
10106	10
10108	9
10111	9
10113	9
10114	9
10116	9
10117	9
10154	7
10155	7
10156	7
10164	7
10165	7
10166	7
10167	7
10811	15
10821	15
10822	15
10823	15
10824	15
10831	15
10833	15
10834	15
10835	15
10904	11
10905	11
10914	11
10915	11
10918	11
11011	4
11012	4
11015	5
11016	5
11021	4
11022	4
11023	4
11024	4
11025	5
11026	5
11027	5
11028	5
11031	4
11032	4
11035	5
11036	5
11041	4

Product Code	Page
11045	5
11114	6
11115	6
11124	6
11125	6
11126	6
11127	6
11135	6
11136	6
11185	7
11186	7
11194	7
11195	7
11196	7
11197	7
11301	9
11302	9
11303	9
11304	10
11305	10
11306	10
11307	9
11308	10
11312	9
11313	9
11314	9
11315	9
11316	9
11505	11
11514	11
11515	11
11518	11
11612	11
11615	11
11690	11
11801	17
11801	80
11802	17
11803	17
11804	17
11811	17
11812	17
11814	17
11815	17
11821	12

Product Code	Page
11822	12
11824	12
11825	12
11827	12
11828	12
11841	8
11842	8
11844	8
11845	8
11846	8
11851	8
11852	8
11854	8
11855	8
11856	8
12010	36
12011	36
12020	36
12021	36
12030	36
12042	36
12060	37
12061	37
12070	37
12110	33
12130	33
12140	33
12144	33
12145	33
12146	33
12147	33
12170	33
12174	33
12175	33
12176	33
12177	33
12344	24
12345	80
12346	24
12354	24
12355	24
12356	24
12374	24
12375	24
12376	24

Product Code	Page
12384	24
12385	24
12386	24
12416	28
12419	28
12432	28
12433	28
12471	80
12472	28
12533	31
12595	32
12596	32
12701	19
12702	19
12703	19
12704	19
12705	19
12706	19
12707	19
12708	19
12709	19
12716	19
12719	19
12800	42
12801	80
12802	42
12803	80
12805	42
12808	42
12810	42
12820	40
12834	80
12835	40
12836	40
12841	19
12842	19
12844	19
12845	19
12847	19
12848	80
12991	32
12992	32
12993	32
12994	32
12995	32

Product Index

Product Code	Page
12996	32
12997	32
13706	31
13840	29
13844	29
13845	29
13846	29
13860	31
13870	30
13874	30
13875	30
13876	30
13924	21
13925	21
13964	21
13965	21
13966	21
13974	22
13975	22
13976	22
13984	22
13985	22
13986	22
14010	43
14015	43
14021	43
14031	43
14041	43
14050	43
14060	43
14070	43
14103	31
14106	30
14163	31
14164	31
14165	31
14166	31
14173	30
14174	30
14175	30
14176	30
14196	29
14201	29
14202	29
14224	23

Product Code	Page
14225	23
14226	23
14227	23
14228	23
14229	23
14232	20
14233	20
14234	20
14235	20
14236	20
14237	20
14238	20
14239	20
14243	20
14246	20
14249	20
14251	20
14252	20
14254	20
14255	20
14257	80
14258	20
14261	19
14262	19
14263	19
14271	19
14272	19
14273	19
14274	19
14275	19
14276	19
14277	19
14278	19
14281	20
14282	20
14283	20
14291	20
14292	20
14293	20
14294	20
14295	20
14296	20
14297	80
14298	80
14301	25

Product Code	Page
14302	25
14303	25
14304	25
14305	25
14306	25
14307	35
14308	35
14309	25
14314	80
14316	27
14317	27
14319	27
14326	80
14327	80
14329	27
14331	27
14332	27
14334	80
14342	26
14344	26
14346	26
14347	26
14349	26
14355	27
14356	26
14357	26
14359	26
14361	26
14362	26
14364	26
14401	25
14402	25
14403	25
14404	29
14405	29
14406	25
14411	31
14412	31
14414	30
14415	30
14421	80
14422	30
14423	30
14424	30
14434	80

Product Code	Page
14435	80
14441	33
14442	80
14444	80
14445	33
14452	33
14453	33
14454	33
14491	31
14494	22
14498	22
14520	26
14521	26
14524	26
14525	26
14527	26
14528	26
14552	27
14553	27
14554	27
14555	27
14601	25
14602	25
14603	25
14642	25
14643	25
14644	25
14645	80
14681	35
14682	35
14691	35
14692	35
14701	29
14702	29
15403	16
15404	16
15405	17
15406	17
15407	16
15408	16
15413	16
15414	16
15415	16
15441	16
15450	16

Product Code	Page
15470	17
15472	17
15473	17
15474	17
15480	16
15482	16
15483	16
15484	16
15521	13
15522	13
15523	13
15524	13
15531	12
15532	12
15533	12
15534	12
15605	13
15606	13
15607	12
15608	12
15633	12
15634	12
15635	12
15695	13
15696	13
15821	13
15981	13
19001	41
19003	41
19010	41
19011	41
19020	41
19021	41
19030	41
19040	41
19041	41
19042	41
19050	41
19051	41
21346	56
21356	56
21365	42
22025	51
22030	51
22035	51

Product Index

Product Code	Page
22040	51
22121	51
22125	51
22130	51
22135	51
27408	46
27410	46
27412	46
27414	46
27416	46
27418	46
27490	46
27492	46
27494	46
27508	46
27510	46
27512	46
27514	46
27516	46
27706	46
27708	46
27710	46
28462	46
28464	46
28578	47
28580	47
28582	47
28584	47
28692	47
29121	58
29122	58
29123	58
29124	58
29125	58
33434	63
33435	63
33436	63
33437	63
33438	63
33444	64
33445	64
33446	64
50501	57
63968	66
69636	66

Product Code	Page
69637	66
2040H	51
5025H	51
5030H	51
5035H	51
5040H	80
5215H	53
5225H	51
5230H	51
5235H	51
5240H	51
C5200	52
C5210	53
C5310	53
C5410	53
C6100	52
C7094	56

Terms and Sales Policy- Prices Effective 01/01/2013

ALL SALES ARE SUBJECT TO THE FOLLOWING TERMS AND CONDITIONS:

PRICE: This price list supersedes all previous price lists. Coloplast Canada reserves the right to revise prices at any time; however, we do attempt to provide the trade with reasonable notification of any changes.

TERMS: Net 30 Days from invoice date to the date payment is post-marked (excluding C.O.D. and Visa/MC orders). No shipments will be released to accounts past due. If payment is not received when due, the invoice is subject to a monthly service charge equal to one and a half percent (1-1/2%) of the unpaid balance.

MINIMUM ORDER: Orders less than \$250 are subject to freight charges and a \$25 handling charge. Orders for Urology products are exempt from the \$250 minimum; however, a \$25 handling fee applies to all orders for Urology products.

SHIPPING: Prices are F.O.B. Coloplast Warehouse. Coloplast Canada Corporation pays for ground transportation charges on all orders of net \$250 or more to destinations in Canada. No freight will be added on backorders. Orders will be shipped by a carrier of our choice. When special handling is requested, the cost of the special service will be added to the invoice. Orders of less than \$250 will be shipped prepaid and freight added to the invoice.

SHIPMENT VERIFICATION: Any discrepancy in the quantity billed and the quantity received resulting from damage (concealed or otherwise), shortages, overages, or wrong product, must be reported to Coloplast Customer Service within 15 days of the date a shipment is received. Any discrepancy reported after 15 days requires Coloplast Management authorization for credit. We will not provide proof of shipment delivery after 90 days. Claims for non-delivery must be notified to Coloplast Customer Service and the Carrier within 15 days from the date of invoice, failing which, Coloplast Canada Corporation shall be under no liability for non-delivery.

INVOICE VERIFICATION: Any discrepancy in an invoiced price and corresponding order price must be reported to Coloplast Customer Service within 30 days of the date of the invoice.

DROP SHIPMENT: All drop shipments are subject to freight charges and a \$25 handling fee.

RETURNED GOODS POLICY: Coloplast Canada Corporation will accept returned goods from our customers, for credit or exchange, under the following conditions:

1. All returns must have written authorization from Customer Service, Coloplast Canada Corporation (1-888-880-8605). This authorization is an estimate, pending inspection by Coloplast of returned goods.
2. Credit will not be issued for products returned without proper documents.
3. Returned Merchandise Authorizations are valid for 30 days and product must be received in designated Coloplast location within 30 from date of issue.
4. Full credit will be given for items shipped in error by Coloplast Canada Corporation if returned within 30 days from date of invoice.
5. Full credit, less a 15% restocking fee, will be given for all products returned with proper documentation as long as the following conditions are met: All returned goods must be in saleable condition with at least three (3) months of shelf life remaining and have been purchased in the last one (1) year. No credit will be issued if the original package integrity is broken, damaged, or if the package is defaced with stickers or markers or if the product is within three (3) months of expiration date. Only full box quantities will be accepted. Credit will not be given for products no longer listed in current Coloplast Canada Corporation price lists or for products that have been replaced by updated versions with significant feature changes or box quantity changes.
6. All returns are subject to Coloplast Canada Corporation inspection. Coloplast Canada Corporation reserves the right to destroy merchandise deemed unfit for sale, whether or not accepted for credit or exchange.
7. Returned goods deemed saleable will be valued at the original cost.
8. Returns must be shipped, prepaid to:

COLOPLAST CANADA CORPORATION

3300 Ridgeway Drive, Unit 12,
Mississauga, ON L5L 5Z9

LIMITED WARRANTY: Coloplast Canada Corporation warrants that all products sold hereunder shall be free from defects in material and workmanship at time of delivery, subject to proper handling, transportation, storage and use of such products. Except as provided above, Coloplast Canada Corporation makes no other warranty, expressed or implied (save for its right to transfer good title to the products), in fact or by law, including any warranty or merchantability or fitness for a particular purpose. All warranties hereunder shall be null and void in the event any alteration of any product is made. The liability of Coloplast Canada Corporation for any product not conforming to such warranty shall be limited to the replacement thereof with a conforming product provided that Coloplast Canada Corporation will be so liable unless (1)

Coloplast Canada Corporation received prompt written notice of the defect within ten (10) days from the date of defect was of should have been discovered and (2) if requested by Coloplast Canada Corporation, the nonconforming product is returned to Coloplast Canada Corporation. In no event shall Coloplast Canada Corporation be liable in tort or in contract for direct, incidental or consequential damages of any kind.

CUSTOMER SERVICE: All orders are subject to our Terms and Sales Policies and to final acceptance and approval by our Customer Service Department at 1-888-880-8605. For your convenience, orders are accepted via telefax at 1-877-820-8206.

REMITTANCE ADDRESS:

COLOPLAST CANADA
Lockbox No T42355C
P.O. Box 4235 Postal Station A
Toronto, ON M5W 5P7

CUSTOMER SERVICE DEPARTMENT: Tel: 1-888-880-8605
CUSTOMER SERVICE DEPARTMENT: Telefax: 1-877-820-8206
FINANCE/ACCOUNTS RECEIVABLE: Tel: 1-800-726-6362

Visit our web site: www.coloplast.ca for additional information

The Coloplast story began back in 1954. Elise Sørensen is a nurse. Her sister Thora has just had an ostomy operation and is afraid to go out, fearing that her stoma might leak in public. Listening to her sister's problems, Elise creates the world's first adhesive ostomy bag. A bag that does not leak, giving Thora – and thousands of people like her – the chance to return to their normal life.

A simple solution with great significance.

Today, our business includes ostomy care, urology and continence care and wound and skin care. But our way of doing business still follows Elise's example: we listen, we learn and we respond with products and services that make life easier for people with intimate healthcare needs.

Ostomy Care
Urology & Continence Care
Wound & Skin Care

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

The Coloplast logo is a registered trademark of Coloplast A/S. © 2011-10
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.



Coloplast Canada Corp
3300 Ridgeway Dr Unit 12
Mississauga ON L5L 5Z9
1-877-820-7008

www.coloplast.ca

E0002N