



Patient Monitoring and Telemetry Supplies

PHILIPS

sense and simplicity

Changes to this edition

May 2010

Product removed

Mainstream CO₂ Supplies

| | | |
|--------|------------|----|
| M1614A | Cable ties | 20 |
|--------|------------|----|

Masimo® LNOP® SpO₂ Sensors

| | | |
|--------------|--|----|
| 989803140321 | LNOP DCI adult reusable finger sensor | 91 |
| 989803140331 | LNOP DCIP pediatric reusable finger sensor | 91 |
| 989803140341 | LNOP TC-I reusable adult ear sensor | 91 |
| 989803140231 | LNOP Adt adult finger adhesive sensor | 92 |
| 989803140261 | LNOP Pdt pediatric finger adhesive sensor | 92 |
| 989803140311 | LNOP Inf-L infant toe adhesive sensor | 92 |
| 989803140291 | LNOP Neo-L neonatal foot adhesive sensor | 93 |
| 989803140301 | LNOP SofTouch NeoPt-L neonatal foot sensor | 93 |
| 989803140351 | Replacement poseys for LNOP SofTouch NeoPt sensors | 94 |
| 989803140361 | Replacement tapes | 94 |

Masimo® LNCS® SpO₂ Sensors

| | | |
|--------------|--|----|
| 989803148231 | LNCS adult adhesive Adtx finger sensor | 95 |
| 989803148241 | LNCS Pdt adhesive finger sensor | 95 |
| 989803148291 | LNCS DC-IP pediatric reusable clip sensor | 95 |
| 989803148301 | LNCS TC-I adult reusable ear clip sensor | 96 |
| 989803148311 | LNCS TF-I adult forehead sensor with head band | 96 |

Dear Customer:

If you're looking for high-quality medical supplies to optimize the performance of your medical devices, you'll find them right here. Whether your focus is on quality patient care, reliability, or cost containment, you can turn to a single source – Philips Medical Supplies. Explore our complete family of supplies, from electrodes to SpO₂ sensors, for every area of care.

With Philips, you can buy your medical equipment and supplies from a single source, simply and easily. Use this guide to find the Philips supplies you need – or go to our website to purchase supplies: <http://www.philips.com/healthcarestore>

Using this guide

Latex Cautions

Philips Medical Supplies currently distributes over 600 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:

LATEX

Caution: Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber. The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.

CE Marking

Products labeled with this symbol conform to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the product is in conformity with the European

Medical Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

The CE symbol provided with some products in the Philips Medical Supplies Reference Guide and online store is intended only as an indication of availability in EU countries and is not intended as a legal product marking.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

Note: Some products identified in this catalogue may not be licensed in accordance with Canadian law.

Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your Philips Medical Supplies representative near you or visit our website: <http://shop.medical.philips.com>

Online Ordering

Online ordering is a simple process:

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

If you are new to the web, we can guide you through the process on the phone in minutes. Many new users tell us that the ordering process is so intuitive that they do not need help, even the first time.

You can search for supplies in the online catalog three easy ways: by application category, Philips instrument, or keyword.

If you currently transmit via MMS/EDI or auto fax, the interactive system provides real-time information and services not available through one-way transmissions.

Online ordering is available to U.S. customers at the following url: <http://www.philips.com/healthcarestore>

Patient Monitoring and Telemetry Supplies

Anesthesia-Gas

| | |
|------------------------------------|----|
| Anesthesia Gas Monitoring Supplies | 11 |
|------------------------------------|----|

Capnography

| | |
|--|----|
| IntelliVue Microstream® CO ₂ Solution | 13 |
| Philips IntelliVue Sidestream CO ₂ Solution | 15 |
| Sidestream CO ₂ Supplies | 18 |
| Philips IntelliVue Mainstream CO ₂ Solution | 19 |
| Mainstream CO ₂ Supplies | 20 |

Spirometry CO₂ / Flow

| | |
|---|----|
| IntelliVue Spirometry CO ₂ Sensors | 21 |
|---|----|

Cardiac Output

| | |
|-----------------------------|----|
| Invasive Cardiac Output | 22 |
| Non-Invasive Cardiac Output | 23 |

Patient Monitoring ECG Cables and Leads

| | |
|---|----|
| IntelliVue 10-Lead ECG Monitoring – AAMI/IEC | 24 |
| 5-Lead Sets for 10-Lead Monitoring – AAMI Color Coding | 24 |
| 5-Lead Sets for 10-Lead Monitoring – IEC Color Coding | 25 |
| IntelliVue 10-Lead (6+4) ECG Monitoring – AAMI/IEC | 26 |
| 6-Lead plus 4-Lead Sets for 10-Lead Monitoring - AAMI Color Coding | 26 |
| 6-Lead plus 4-Lead Sets for 10-Lead Monitoring – IEC Color Coding | 27 |
| 6-Lead Monitoring – AAMI/IEC | 28 |
| 6-Lead Sets for 6-Lead Monitoring – AAMI Color Coding | 28 |
| 6-Lead Sets for 6-Lead Monitoring – IEC Color Coding | 29 |
| 5-Lead Monitoring – AAMI/IEC | 30 |
| 5-Lead Sets for 5-Lead Monitoring – AAMI Color Coding | 30 |
| 5-Lead Sets for 5-Lead Monitoring – IEC Color Coding | 31 |
| 3-Lead Monitoring – AAMI/IEC | 32 |
| 3-Lead Sets for 3-Lead Monitoring – AAMI Color Coding | 32 |
| 3-Lead Sets for 3-Lead Monitoring – IEC Color Coding | 33 |
| Radiolucent ECG Lead Sets – AAMI/IEC Color Coding | 34 |
| SureSigns Combined ECG Cable/Leads | 34 |
| A1/A3/C1 ECG Cable/Leads | 35 |
| 3-Lead CMS Safety ECG Patient Trunk Cables – 12-Pin Connector | 35 |
| 5-Lead CMS Safety ECG Patient Trunk Cables – 12-Pin Connector | 36 |
| 3-Lead CMS ECG Safety Trunk Cables: 8-Pin Connector | 36 |
| 3-Lead Sets for CMS Safety ECG Patient Trunk Cables – AAMI Color Coding | 36 |
| 3-Lead Sets for CMS Safety ECG Patient Trunk Cables – IEC Color Coding | 37 |
| 5-Lead Sets for CMS Safety ECG Patient Trunk Cables – AAMI Color Coding | 37 |
| 5-Lead Sets for CMS Safety ECG Patient Trunk Cables – IEC Color Coding | 37 |
| EEG Cable and Leads (electroencephalogram) | 38 |
| Cable and Leads Accessories | 39 |

Telemetry ECG Cables and Leads

| | |
|---|----|
| IntelliVue Telemetry 3-Lead Sets – AAMI Color Coding | 40 |
| IntelliVue Telemetry 3-Lead Sets – IEC Color Coding | 40 |
| IntelliVue Telemetry 5-Lead Sets – AAMI Color Coding | 41 |
| IntelliVue Telemetry 5-Lead Sets for 5-Lead Monitoring – IEC Color Coding | 41 |
| IntelliVue Telemetry 6-Lead Sets – AAMI Color Coding | 42 |
| IntelliVue Telemetry 6-Lead Sets – IEC Color Coding | 42 |
| IntelliVue Telemetry Tether Cables | 43 |
| IntelliVue Shower Shields & Alignment Guides | 44 |
| Legacy HP/Agilent Telemetry Leads | 45 |
| ECG Extender Cables for TeleMon | 46 |

Disposable ECG Electrodes

| | |
|---|----|
| Disposable Adult ECG Monitoring Electrodes | 47 |
| Disposable Adult Metallic ECG Lead Wire Electrodes - AAMI Color Coding | 49 |
| Disposable Adult Metallic ECG Lead Wire Electrodes - IEC Color Coding | 49 |
| Disposable Adult Radiolucent ECG Lead Wire Electrode - AAMI Color Coding | 50 |
| Disposable Adult Radiolucent ECG Lead Wire Electrode - IEC Color Coding | 50 |
| Disposable Pediatric ECG Monitoring Electrodes | 50 |
| Disposable Infant/Pediatric Metallic Lead Wire Electrodes: AAMI Color Coding | 51 |
| Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: AAMI Color Coding | 51 |
| Disposable Infant/Pediatric Metallic Lead Wire Electrodes: IEC Color Coding | 52 |
| Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: IEC Color Coding | 52 |
| EEG Electrode (electroencephalogram) | 52 |

Invasive Blood Pressure

| | |
|---|----|
| Disposable Blood Pressure Kits | 53 |
| Reusable Invasive Blood Pressure Solution | 54 |

Non-Invasive Blood Pressure

| | |
|--|----|
| Multi-Patient Comfort Cuffs - Single Hose | 55 |
| Multi-Patient Easy Care Cuffs - 1 Hose | 57 |
| Multi-Patient Easy Care Cuffs - 2 Hose | 62 |
| Multi-Patient Traditional Cuffs - 1 Hose | 64 |
| Multi-Patient NIBP Cuff Kits - 1 Hose | 66 |
| Single Patient Disposable Cuffs | 68 |
| Soft Single Patient Disposable Cuffs - 1 Hose | 70 |
| Soft Single Patient Disposable Cuffs - 2 Hose | 73 |
| Blood Pressure Air Hoses | 75 |
| NIBP Cuff Caps | 76 |
| Bulb and Valve Kit for Manual Sphygmomanometer | 76 |
| NIBP Cuff Adapters | 77 |

Paper

| | |
|-----------------|----|
| Recording Paper | 79 |
|-----------------|----|

Pulse Oximetry

| | |
|--|----|
| Philips Reusable SpO ₂ Sensors | 80 |
| Philips Disposable SpO ₂ Sensors | 86 |
| SpO ₂ Adapter Cables for Philips Sensors | 88 |
| Adhesive SpO ₂ Sensors <i>(available from Philips in Europe and Japan only)</i> | 90 |
| Masimo® LNOP® SpO ₂ Adapter Cables | 91 |
| Masimo® LNCS® SpO ₂ Adapter Cable | 91 |

Temperature Monitoring

| | |
|---|----|
| Disposable Skin Surface Temperature Probe | 92 |
| Disposable Esophageal/Rectal Temperature Probes (<i>meets EN12470 Standards</i>) | 92 |
| Disposable Esophageal/Stethoscope Temperature Probes (<i>meets EN12470 Standards</i>) | 93 |
| Disposable Foley Catheter Temperature Probes (<i>meets EN12470 Standards</i>) | 94 |
| Disposable Temperature Probes Extension and Adapter Cables | 95 |
| Reusable Temperature Probes for Continuous Monitoring (<i>meets EN12470 Standard</i>) | 96 |
| Reusable Temperature Probes for SureSigns VS1 Monitors | 98 |
| Reusable Temperature Probes for SureSigns VM4/VS3/MP5 Monitors | 98 |
| Reusable Temperature Probes for SureSigns A1/A3/C1/VM4-8 Series Monitors | 99 |

Batteries, Cases, Pouches

| | |
|---|-----|
| Patient Monitoring Batteries | 100 |
| Telemetry Batteries | 100 |
| Patient Monitoring Cases, Bags, Pads | 101 |
| Telemetry Pouches and Protective Covers | 102 |

Patient Monitoring and Telemetry Supplies

Anesthesia-Gas

Anesthesia Gas Monitoring Supplies

| Product # | Description | Packaging | Use With: |
|--------------|---|------------------|--|
| 13901A CE | Reusable Nafion sample tube Length: 8 feet (2.4m) | 1 tube | Module: M1026A |
| 13902A CE | Disposable elbow airway adapter for single patient use. | 20 adapters | Modules: M1026B, M1013A (G1), M1019A (G5) |
| 13904A CE | Disposable bacteria filter (0.45 micron) for single patient use. | 20 filters | Module: M1026A |
| 13905A CE | Disposable hybrid Nafion polyethylene sample tube Nafion length: 6 feet (1.8m) Polyethylene length: 9 feet (2.7m) | 1 tube | Module: M1026A |
| 13907A | Calibration tube assembly | 1 assembly | Modules: M1025A, M1025B |
| M1612A CE | Straight airway adapter for single patient use. | 20 adapters | Modules: M1019A (G5), M1013A (G1), M1026A, M1026B |
| M1655A CE | Reusable gas exhaust return line. (not available in U.S.) Length: (3.6m)~ | 1 line/5 filters | M1026A |
| M1655B CE | Reusable gas exhaust return line. Length: 12 feet (3.6m) | 1 line/5 filters | Modules: M1013A (G1), M1019A (G5), M1026B (Note: M1026B not in the US) |



13901A



13902A



M1612A



M1657B

Anesthesia Gas Monitoring Supplies *(continued)*

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|--|
| M1656B CE | Disposable gas exhaust return filter. (use with M1655B) | 20 filters | Modules: M1013A (G1), M1019A (G5), M1026A, M1026B (Note: M1026B not in the US) |
| M1657B CE | Water trap (use with M1658A) | 25 traps | Modules: M1013A (G1), M1019A (G5), M1026B, M1026A with Opt. A02/A05 |
| M1658A CE | Gas sample tubing for single patient use. (use with M1657B) | 20 tubes | Modules: M1013A (G1), M1019A (G5), M1026B, M1026A with Opt. A02/A05 |
| M1659A | Calibration bag. (use with M1660A and M1662A) | 1 bag | Modules: M1013A (G1), M1019A (G5), M1026B |
| M1662A | Anesthetic CAL gas. Volume of gas = 5 liters. (use with M1659A) Shelf life: 1 year | 1 cylinder | Monitor Modules: M1013A (G1), M1019A (G5), M1026A, M1026B |

Capnography

Sidestream CO₂

IntelliVue Microstream® CO₂ Solution

| Product # | Description | Packaging | Use With: |
|--------------------|--|------------|--|
| 989803159571 CE | VitaLine™ H Set, Adult/Pediatric | 25 sets | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803159581 CE | VitaLine H Set, Infant/Neonatal | 25 sets | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160241 CE | FilterLine® Set Long, Adult/Pediatric | 25 sets | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160251 CE | FilterLine H Set Long, Adult/Pediatric | 25 sets | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160261 CE | FilterLine H Set Long, Infant/Neonatal | 25 sets | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160271 CE | Smart CapnoLine® O ₂ Long, Pediatric | 25 cannula | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160281 CE | Smart CapnoLine O ₂ Plus Long, Adult | 25 cannula | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| 989803160301 CE | Smart CapnoLine Plus Long, Adult | 25 cannula | Microstream Ext: M3015A, VM1: 863266; MP5: M8105A |
| M1920A CE | FilterLine Set, Adult/Pedi, Intubated | 25 sets | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M1921A CE | FilterLine H Set, Adult/Pedi, Intubated | 25 sets | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M1923A CE | FilterLine H Set, Infant/Neonatal, Intubated | 25 sets | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M2520A CE | Smart CapnoLine O ₂ , Pediatric, Non-intubated Dual Purpose | 25 cannula | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M2522A CE | Smart CapnoLine O ₂ plus, Adult/intermed, Non-intubated Dual Purpose | 25 cannula | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M2524A CE | Smart CapnoLine, Pediatric, Non-intubated Single Purpose | 25 cannula | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M2526A CE | Smart CapnoLine plus, Adult/intermed, Non-intubated Single Purpose | 25 cannula | Microstream Ext: M3015A, VM8: 863066, 863088, VM1: 863266, MP5: M8105A |
| M4680A CE | CapnoLine H O ₂ , Adult, nasal, (up to 24 hours) Non-intubated Single Purpose | 25 cannula | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M4681A CE | CapnoLine H O ₂ , Pediatric, nasal, (up to 24 hours) Non-intubated Single Purpose | 25 cannula | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M4686A CE | NIV Line, Adult, nasal (up to 8 hours), mask, single purpose | 25 lines | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M4687A CE | NIV Line, Pediatric, nasal (up to 8 hours), mask, single purpose | 25 lines | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |

IntelliVue Microstream® CO₂ Solution (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-------------|---|
| M4689A CE | Capnoline H, Adult, nasal (up to 24 hours), Non-intubated single purpose | 25 cannula | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M4690A CE | Capnoline H, Pediatric, nasal (up to 24 hours), Non-intubated single purpose | 25 cannula | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M4691A CE | Capnoline H, Infant/Neonatal, nasal (up to 24 hours), Non-intubated single purpose | 25 cannula | Microstream Ext: M3015A, VM1: 863266, MP5: M8105A |
| M2267A | Calibration Regulator | 1 regulator | Microstream Ext. M3015A, Multi- Measurement Server: M3002A, MP2: M8102A |



989803159571



989803159581



989803160261



M1920A



M1921A



M1923A



M2522A



M2524A



M2526A



M4681A



M4689A



M4690A



M4691A

Philips IntelliVue Sidestream CO₂ Solution

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|--|
| M2741A CE | Sidestream CO ₂ Sensor | 1 sensor | Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2744A CE | CO ₂ Nasal Cannula - Adult, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2745A CE | CO ₂ Nasal Cannula - Pediatric, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2746A CE | CO ₂ Nasal Cannula - Infant, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2750A CE | CO ₂ / O ₂ Nasal Cannula - Adult, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2751A CE | CO ₂ / O ₂ Nasal Cannula - Pediatric, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |



M2741A



M2744A



M2745A



M2746A



M2750A



M2751A

Philips IntelliVue Sidestream CO₂ Solution (continued)

| Product # | Description | Packaging | Use With: |
|--------------------|--|------------|--|
| 989803144471 CE | CO ₂ / O ₂ Nasal Cannula - Infant/Neonate, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2756A CE | CO ₂ Oral-Nasal Cannula - Adult, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2757A CE | CO ₂ Oral-Nasal Cannula - Pediatric, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2760A CE | CO ₂ / O ₂ Oral-Nasal Cannula - Adult, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2761A CE | CO ₂ / O ₂ Oral-Nasal Cannula - Pediatric, disposable, Non-intubated | 10 cannula | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |



989803144471



M2756A



M2757A



M2760A



M2761A

Philips IntelliVue Sidestream CO₂ Solution (continued)

| Product # | Description | Packaging | Use With: |
|--------------------|---|-----------|--|
| M2768A CE | Airway Adapter Set - ET greater than 4.0 mm, disposable, intubated | 10 sets | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| 989803144531 CE | Airway Adapter Set - ET equal to or less than 4.0 mm | 10 sets | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2772A CE | Airway Adapter Set H - ET greater than 4.0 mm, disposable, intubated | 10 sets | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2773A CE | Airway Adapter Set H - ET equal to or less than 4.0 mm, disposable, intubated | 10 sets | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2776A CE | Straight Sample Line, disposable, non-humidified, intubated | 10 lines | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2777A CE | Straight Sample Line H, humidified, intubated | 10 lines | Sidestream Sensor: M2741A with Sidestream Ext. Server: M3014A, MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |



M2768A



989803144531



M2772A



M2773A



M2776A



M2777A

Sidestream CO₂ Supplies

| Product # | Description | Packaging | Use With: |
|--------------|---|-------------|-----------------------------------|
| M1612A CE | Straight airway adapter for single patient use. | 20 adapters | Legacy HP/Agilent: M1015A |
| 13901A CE | Reusable Nafion sample tube Length: 8 feet (2.4m) | 1 tube | Legacy HP/Agilent: M1015A, M1016A |
| 13902A CE | Disposable elbow airway adapter for single patient use. | 20 adapters | Legacy HP/Agilent: M1015A, M1016A |
| 13904A CE | Disposable bacteria filter (0.45 micron) for single patient use. | 20 filters | Legacy HP/Agilent: M1015A |
| 13905A CE | Disposable hybrid Nafion polyethylene sample tube Nafion length: 6 feet (1.8m) Polyethylene length: 9 feet (2.7m) | 1 tube | Legacy HP/Agilent: M1015A |



M1612A



13901A



13902A

Mainstream CO₂

Philips IntelliVue Mainstream CO₂ Solution

| Product # | Description | Packaging | Use With: |
|-------------------------|---|-----------------|--|
| M2501A CE | Mainstream CO ₂ Sensor | 1 sensor | M3014A Capnography Extension Server, MP2: M8102A, MP5: M8105A, Multi-Measurement Module X2: M3002A |
| M2513A CE | Reusable airway adapter, Adult/Pediatric (color: black) ET tube size: greater than 4mm, deadspace less than 5cc | 1 adapter | M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2516A CE | Reusable CO ₂ airway adapter, Infant/Neonatal (color: red) ET tube size: less than 4mm, deadspace less than 1cc | 1 adapter | M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2533A CE | Single patient disposable CO ₂ airway adapter, Adult/Pediatric (color: clear) ET tube size: greater than 4mm, deadspace less than 5cc | 10 adapters | M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2536A CE | Single patient disposable CO ₂ airway adapter, Infant/Neonatal (color: violet) ET tube size: less than 4mm, deadspace less than 1cc | 10 adapters | M2501A Sensor with M3014A Capnography Extension, Monitors: MP5: M8105A, MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M2505A | Gas Cylinder Regulator | 1 gas regulator | M2506A verification gas |
| REVISED M2506A CE | Gas Verification Gas | 1 gas cylinder | Mainstream CO ₂ Sensor: M2501A |



M2501A

Mainstream CO₂ Supplies

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| 14363A CE | Small airway adapter, reusable, non-sterile Patient size: weighing less than 10kg, deadspace: 2cc | 1 adapter | CO ₂ Sensor M1460A, Legacy HP/ Agilent 78352C, 78834C, Modules: M1016A, M3016A |
| M1465A CE | Standard airway adapter, reusable, non-sterile Patient weight: 10kg or more, deadspace: 7.5cc | 1 adapter | CO ₂ Sensor M1460A, Legacy HP/ Agilent 78352C, 78834C, Modules: M1016A, M3016A |

Spirometry CO₂ / Flow

IntelliVue Spirometry CO₂ Sensors

| Product # | Description | Packaging | Use With: |
|--------------|--|------------|---|
| M2781A CE | Combined CO ₂ / Flow Sensor. Adult / Pediatric, disposable ET: greater than 5.5mm (color: clear) | 10 sensors | Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext. |
| M2782A CE | Combined CO ₂ / Flow Sensor Neonatal. Infant, disposable ET: 2.5mm to 4mm (color: violet) | 10 sensors | Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext. |
| M2783A CE | CO ₂ / Flow Sensor Pediatric, disposable ET: 3.5mm to 6mm (color: green) | 10 sensors | Spirometry Module: M1014A with M2501A Mainstream CO ₂ Sensor and M3014A Capnography Ext. |
| M2785A CE | Flow Sensor Adult/Pediatric, disposable ET: greater than 4mm (color: clear) | 10 sensors | M1014A Spirometry Module |
| M2786A CE | Flow Sensor Infant/Neonatal ET: less than 4mm, (color: violet) | 10 sensors | M1014A Spirometry Module |



M2781A



M2782A



M2783A



M2785A



M2786A

Cardiac Output

Invasive Cardiac Output

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 15244A CE | Remote Startup Switch and Cable | 1 cable | Module: M1012A Opt. C12, Hemodynamic Measurement Server: M3012A, Capnography Extension: M3014A |
| 23001A CE | Reusable Injectate Temp Probe to measure cardiac output. For use with Baxter CO-Set™ kit only. Length: 7'9 feet (2.4m) | 1 probe | Module: M1012A, Hemodynamic Measurement Server: M3012A, Capnography Extension: M3014A |
| 23001B CE | Reusable Injectate Temp Probe to measure cardiac output. For use with Baxter CO-Set™ kit only. Length: 1'9 feet (0.5m) | 1 probe | Module: M1012A, Hemodynamic Measurement Server: M3012A, Capnography Extension: M3014A |
| 23002A CE | Ice bath temperature probe | 1 probe | Module: M1012A, Hemodynamic Measurement Server: M3012A |
| M1642A CE | Cardiac Output Cable | 1 cable | Module: M1012A, Hemodynamic Measurement Server: M3012A |
| M1643A CE | Cardiac Output Cable for use with pulse contour analysis and transpulmonary thermodilution. Length: 15.75 feet (4.8m) | 1 cable | Module: M1012A, Hemodynamic Measurement Server: M3012A, Capnography Extension: M3014A |
| M1646A CE | PICCO temperature probe to measure cardiac output. Not available in N. America or Japan | 1 probe | Module: M1012A |



15244A



23001A



23002A



M1643A

Non-Invasive Cardiac Output

| Product # | Description | Packaging | Use With: |
|--------------|------------------------------------|----------------------------|--------------------------------|
| M4620A CE | ICG Cable (Impedance Cardiography) | 1 cable | 862146 ICG Hemodynamic Monitor |
| M4608A CE | ICG Disposable Sensor | 4/per pack 10 packs/box | |



M4620A

Patient Monitoring ECG Cables and Leads

IntelliVue 10-Lead (5+5) ECG Monitoring – AAMI/IEC

(for conventional diagnostic 12 lead ECG, use one limb 5-lead set and one chest 5-lead set)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1663A CE | 10-Lead ECG Patient Trunk (5 +5) - AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 6.5 feet (2m) | 1 cable | Multi-Measurement Modules: M3000A, M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; |
| M1949A CE | 10-Lead ECG Patient Trunk (5 +5) - AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 8.9 feet (2.7m) | 1 cable | Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG Modules: M1001B, M1002B |



M1663A


5-Lead Sets for 10-Lead Monitoring – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|------------------------------|------------|---|
| M1968A | ICU grabber, limb, AAMI | 1 lead set | Use with M1663A or M1949A 10-lead patient trunk cable |
| M1976A | ICU grabber, chest, AAMI | 1 lead set | |
| M1644A | ICU snap, limb, AAMI | 1 lead set | |
| M1602A | ICU snap, chest, AAMI | 1 lead set | |
| M1973A | OR grabber, limb, AAMI | 1 lead set | |
| M1979A | OR grabber, chest, AAMI | 1 lead set | |
| M1647A | Miniclip grabber, unshielded | 1 lead set | |



M1968A

M1976A



M1973A



M1979A

5-Lead Sets for 10-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|------------------------------|------------|---|
| M1971A CE | ICU grabber, limb, IEC | 1 lead set | Use with M1663A or M1949A 10-lead patient trunk cable |
| M1978A CE | ICU grabber, chest, IEC | 1 lead set | |
| M1645A CE | ICU snap, limb, IEC | 1 lead set | |
| M1604A CE | ICU snap, chest, IEC | 1 lead set | |
| M1974A CE | OR grabber, limb, IEC | 1 lead set | |
| M1984A CE | OR grabber, chest, IEC | 1 lead set | |
| M1648A CE | Miniclip grabber, unshielded | 1 lead set | |



M1971A



M1978A



M1974A



M1984A

IntelliVue 10-Lead (6+4) ECG Monitoring – AAMI/IEC

(for conventional diagnostic 12 lead ECG, use one limb 6-lead set and one chest 4-lead set)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1665A CE | 10-Lead ECG Patient Trunk (6 + 4) - AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 8.9 feet (2.7m) | 1 cable | M3000A M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18, M8102A Opt C12, M8105A Opt C12 |



M1665A


6-Lead plus 4-Lead Sets for 10-Lead Monitoring - AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|---------------------------------|------------|---|
| M1680A | 6-lead ICU grabber, limb, AAMI | 1 lead set | Use with M1665A 10-lead patient trunk cable |
| M1532A | 4-lead ICU grabber, chest, AAMI | 1 lead set | |
| M1682A | 6-lead ICU snap, limb, AAMI | 1 lead set | |
| M1537A | 6-lead ICU snap, chest, AAMI | 1 lead set | |
| M1684A | 6-lead OR grabber, limb, AAMI | 1 lead set | |
| M1557A | 4-lead OR grabber, chest, AAMI | 1 lead set | |



M1680A

M1532A



M1557A

M1684A


6-Lead plus 4-Lead Sets for 10-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|--------------------------------|------------|---|
| M1681A CE | 6-lead ICU grabber, limb, IEC | 1 lead set | Use with M1665A 10-lead patient trunk cable |
| M1533A CE | 4-lead ICU grabber, chest, IEC | 1 lead set | |
| M1683A CE | 6-lead ICU snap, limb, IEC | 1 lead set | |
| M1538A CE | 6-lead ICU snap, chest, IEC | 1 lead set | |
| M1685A CE | 6-lead OR grabber, limb, IEC | 1 lead set | |
| M1558A CE | 4-lead OR grabber, chest, IEC | 1 lead set | |



M1681A

M1533A



M1558A

M1685A

6-Lead Monitoring – AAMI/IEC

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1667A CE | 6-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 6.5 feet (2m) | 1 cable | M3000A M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18, M8102A Opt C12, M8105A Opt C12 |



M1667A

6-Lead Sets for 6-Lead Monitoring – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|--------------------------------|------------|-------------------------------------|
| M1680A | 6-lead ICU grabber, limb, AAMI | 1 lead set | Use with M1667A patient trunk cable |
| M1682A | 6-lead ICU snap, limb, AAMI | 1 lead set | |
| M1684A | 6-lead OR grabber, limb, AAMI | 1 lead set | |




M1680A



M1684A

6-Lead Sets for 6-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|-------------------------------|------------|-------------------------------------|
| M1681A CE | 6-lead ICU grabber, limb, IEC | 1 lead set | Use with M1667A patient trunk cable |
| M1683A CE | 6-lead ICU snap, limb, IEC | 1 lead set | |
| M1685A CE | 6-lead OR grabber, limb, IEC | 1 lead set | |



M1681A



M1685A

5-Lead Monitoring – AAMI/IEC

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1668A CE | 5-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector: AAMI and IEC labels included. Length: 8.9 feet (2.7m) | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863087, 863088 |



M1668A

5-Lead Sets for 5-Lead Monitoring – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|------------------------------------|------------|--|
| M1968A | ICU grabber, limb, AAMI | 1 lead set | Use with M1668A 5-lead patient trunk cable |
| M1644A | ICU snap, limb, AAMI | 1 lead set | |
| M1973A | OR grabber, limb, AAMI | 1 lead set | |
| M1647A | Miniclip grabber, unshielded, AAMI | 1 lead set | |



M1968A




M1647A



M1973A

5-Lead Sets for 5-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|-----------------------------------|------------|--|
| M1971A CE | ICU grabber, limb, IEC | 1 lead set | Use with M1668A 5-lead patient trunk cable |
| M1645A CE | ICU snap, limb, IEC | 1 lead set | |
| M1974A CE | OR grabber, limb, IEC | 1 lead set | |
| M1648A CE | Miniclip grabber, unshielded, IEC | 1 lead set | |



M1971A



M1974A



M1648A

3-Lead Monitoring – AAMI/IEC

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1669A CE | 3-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 8.9 feet (2.7m) | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088 |



M1669A

3-Lead Sets for 3-Lead Monitoring – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|--|------------|--|
| M1671A | ICU grabber, AAMI | 1 lead set | Use with M1669A 3-lead patient trunk cable |
| M1673A | ICU snap, AAMI | 1 lead set | |
| M1675A | OR grabber, AAMI | 1 lead set | |
| M1622A | Miniclip grabber, unshielded (18 inches/0.45m) | 1 lead set | |
| M1624A | Miniclip grabber, unshielded (28 inches/0.7m) | 1 lead set | |



M1671A



M1675A



M1624A

3-Lead Sets for 3-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|-----------------------------------|------------|--|
| M1672A CE | ICU grabber, IEC | 1 lead set | Use with M1669A 3-lead patient trunk cable |
| M1674A CE | ICU snap, IEC | 1 lead set | |
| M1678A CE | OR grabber, IEC | 1 lead set | |
| M1626A CE | Miniclip grabber, unshielded, IEC | 1 lead set | |



M1672A



M1678A



M1626A

Radiolucent ECG Lead Sets – AAMI/IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|--|--------------------|---|
| M1649A CE | Radiolucent unshielded 5-lead set with safety connectors. Grabber termination, 5 individual wires, IEC/AAMI color coding Length: 36 inches (914mm) | 5 individual leads | Telemetry: 862231 (TRx4851A), 862439 (TRx4841A), 862439 (US only), 862108: M2601B Transmitter, MMM X2: M3002A; MP2: M8102A; MP5: M8105A |



M1649A

SureSigns Combined ECG Cable/Leads

| Product # | Description | Packaging | Use With: |
|--------------------|---------------------------|-----------|---|
| 989803143171 CE | 3 Lead Set, grabber, IEC | 1 cable | SureSigns Monitors 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088; MP2: |
| 989803143181 | 3 Lead set, grabber, AAMI | 1 cable | |
| 989803143191 CE | 5 Lead set, grabber, IEC | 1 cable | M8102A, MP5: M8105A, MMM X2: M3002A |
| 989803143201 | 5 Lead set, grabber, AAMI | 1 cable | |



989803143171



989803143181



989803143191



989803143201

A1/A3/C1 ECG Cable/Leads

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|---|
| M3913A CE | 3-lead ECG patient trunk cable, AAMI or IEC Length: 10 feet (3m) | 1 cable | M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, 863053, 863054 |
| M3914A CE | 3-lead snap, IEC | 1 lead set | Use with M3913A patient trunk cable |
| M3915A CE | 3 lead snap, AAMI | 1 lead set | |

3-Lead CMS Safety ECG Patient Trunk Cables – 12-Pin Connector

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1500A | 3-lead ECG patient trunk cable, AAMI color coding Length: 8.9 feet | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, C3: 863474, SureSigns: 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088 |
| M1540C | 3-lead ECG patient trunk cable, AAMI color coding Length: 2.9 feet | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Module: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A |
| M1510A CE | 3-Lead ECG patient trunk cable, IEC color coding Length: 2.7M | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Module: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, C3: 862474, 862478, SureSigns: 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088 |
| M1550C CE | 3-lead IEC patient trunk cable, IEC color coding Length: 0.9m | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A |

5-Lead CMS Safety ECG Patient Trunk Cables – 12-Pin Connector

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M1520A | 5-lead ECG patient trunk cable, AAMI color coding Length: 8.9 feet | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Module: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, C3: 862474, 862478, SureSigns: 863063, 863064, 863065, 863066, 863068 |
| M1560C | 5-lead ECG patient trunk cable, AAMI color coding Length: 2.9 feet | 1 cable | ECG Modules: M1001B, M1002B, Multi-Measurement M3000A, M3001A, M3002A, Monitors: M8102A, M8105A |
| M1530A CE | 5-lead ECG patient trunk cable, IEC color coding Length: 2.7m | 1 cable | Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, MP2: M8102A, MP5: M8105A, C3: 862474, 862478, SureSigns: 863063, 863064, 863065, 863066, 863068, 863086, 863087, 863088 |
| M1570C CE | 5-lead ECG patient trunk cable, IEC color coding Length: 0.9m | 1 cable | ECG Modules: M1001B, M1002B, Multi-Measurement M3000A, M3001A, M3002A, Monitors: MP2: M8102A, MP5: M8105A |

3-Lead CMS ECG Safety Trunk Cables: 8-Pin Connector

| Product # | Description | Packaging | Use With: |
|--------------|--------------------------------------|-----------|------------------------|
| M1580A | 3-lead ECG patient trunk cable, AAMI | 1 cable | M2400A, M2406A, M2424A |
| M1600A | 5-lead ECG patient trunk cable, AAMI | 1 cable | |
| M1590A CE | 3-lead ECG patient trunk cable, IEC | 1 cable | |
| M1610A CE | 5-lead ECG patient trunk cable, IEC | 1 cable | |

3-Lead Sets for CMS Safety ECG Patient Trunk Cables – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|---|------------|--------------------------------|
| M1601A | 3-lead OR set, grabbers, AAMI | 1 lead set | Use with M1500A, M1580A or |
| M1603A | 3-lead ICU set, grabbers, AAMI | 1 lead set | M1540C CMS patient trunk cable |
| M1605A | 3-lead ICU set, snap, AAMI | 1 lead set | |
| M1608A | 3-lead miniclip, AAMI Unshielded (18 inches/0.45m) | 1 lead set | |
| M1609A | 3-lead miniclip, AAMI Unshielded (28 inches/0.7m) | 1 lead set | |

3-Lead Sets for CMS Safety ECG Patient Trunk Cables – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|------------------------------------|------------|---|
| M1611A CE | 3-lead OR set, grabbers, IEC | 1 lead set | Use with M1510A, M1590A or M1550C CMS patient trunk cable |
| M1613A CE | 3-lead ICU set, grabbers, IEC | 1 lead set | |
| M1615A CE | 3-lead ICU set, snap, IEC | 1 lead set | |
| M1619A CE | 3-lead miniclip, IEC Unshielded | 1 lead set | |

5-Lead Sets for CMS Safety ECG Patient Trunk Cables – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-----------|-------------------------------------|------------|---|
| M1621A | 5-lead OR set, grabbers, AAMI | 1 lead set | Use with M1520A, M1600A or M1560C CMS patient trunk cable |
| M1623A | 5-lead ICU set, grabbers, AAMI | 1 lead set | |
| M1625A | 5-lead ICU set, snap, AAMI | 1 lead set | |
| M1629A | 5-lead miniclip, AAMI Unshielded | 1 lead set | |

5-Lead Sets for CMS Safety ECG Patient Trunk Cables – IEC Color Coding


| Product # | Description | Packaging | Use With: |
|--------------|------------------------------------|------------|---|
| M1631A CE | 5-lead OR set, grabbers, IEC | 1 lead set | Use with M1530A, M1610A or M1570C CMS patient trunk cable |
| M1633A CE | 5-lead ICU set, grabbers, IEC | 1 lead set | |
| M1635A CE | 5-lead ICU set, snap, IEC | 1 lead set | |
| M1639A CE | 5-lead miniclip, IEC Unshielded | 1 lead set | |

EEG Cable and Leads (electroencephalogram)

| Product # | Description | Packaging | Use With: |
|--------------|--|---------------------|------------------------|
| M2268A CE | EEG 5-lead patient trunk cable, 2.7 meters | 1 cable | M1027A EEG Module |
| M1931A CE | Reusable EEG adult cup electrode | 1 lead set | M1027A EEG Module |
| M1932A CE | Reusable EEG pediatric cup electrodes | 1 lead set | M1027A EEG Module |
| M1934A CE | Reusable EEG miniclip lead set | 1 lead set | M1027A EEG Module |
| M1937A | EEG conductive adhesive paste | 1 tube | M1931A, M1932A, M1934A |
| M1935A CE | Disposable EEG Snap electrode | 25/pack 100/case | Adult EEG Monitoring |



M2268A



M1931A




M1935A

Cable and Leads Accessories

| Product # | Description | Packaging | Use With: |
|--------------------|--|---------------|---|
| 989803148841 CE | IntelliVue Cable Management Kit Keeps patient cables organized. | 5 organizers | Any cables |
| 989803148861 | IntelliVue Trunk Cable Covers | 10 covers | For M1663A or M1949A Trunk Cables |
| M1501A | Combiner for 3-lead sets. Use with CMS ECG individual lead sets. | 3 combiners | use with ECG, safety cable accessory, non-shielded |
| M1502A | Combiner for 5-lead sets. Use with CMS ECG individual leads sets | 3 combiners | use with ECG, safety cable accessory, non-shielded |
| M1503A | Organizer for shielded 3-lead sets | 10 organizers | use with all shielded lead sets (AAMI/IEC) |
| M1504A | Organizer for shielded 5-lead sets | 10 organizers | use with all shielded lead sets (AAMI/IEC) |
| M1509A | Bedsheet clip for trunk cables | 3 clips | Any cables |
| M2281A | Patient cable organizer (hook) | 8 hooks | Any cables |
| M1636A | Organizer for unshielded 3-lead sets | 10 organizers | use with all unshielded lead sets (AAMI/IEC) |
| M1638A | Organizer for unshielded 5-lead sets | 10 organizers | use with all unshielded lead sets (AAMI/IEC) |
| M1664A | Organizer for shielded 4-lead sets | 10 organizers | (AAMI) M1532A, M1537A, M1557A (IEC) M1533A, M1538A, M1558A |
| M1679A | Organizer for shielded 6-lead sets | 10 organizers | (AAMI) M1680A, M1682A, M1684A, M1532A, M1537A, M1557A (IEC) M1681A, M1683A, M1685A, M1558A, M1533A, M1538A |



989803148841



989803148861



M1501A



M1504A



M1509A



M2281A



M1636A



M1638A



M1664A



M1679A

Telemetry ECG Cables and Leads

IntelliVue Telemetry 3-Lead Sets – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-------------------------------|--|------------|---|
| REVISED 989803151971 CE | 3-lead grabber and shower shield, AAMI | 1 lead set | Telemetry: 862108: M2601B; 862439: U.S. only TRx4841A; |
| REVISED 989803151991 CE | 3-lead snap and shower shield, AAMI | 1 lead set | 862231: Non-U.S. TRx4851A |



989803151971



989803151991

IntelliVue Telemetry 3-Lead Sets – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|---------------------------------------|------------|--|
| 989803151981 CE | 3-lead grabber and shower shield, IEC | 1 lead set | Telemetry: 862108: M2601B; 862231: TRx4851A |
| 989803152001 CE | 3-lead snap and shower shield, IEC | 1 lead set | |



989803151981



989803152001

IntelliVue Telemetry 5-Lead Sets – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|----------------------------|--|------------|---|
| REVISÉD 989803152051 CE | 5-lead grabber and shower shield, multi-color wire | 1 lead set | Telemetry: 862108: M2601B; 862439: U.S. only TRx4841A; 862231: Non-U.S. |
| REVISÉD 989803152071 CE | 5-lead snap and shower shield, multi-color wire | 1 lead set | TRx4851A |



989803152051



989803152071

IntelliVue Telemetry 5-Lead Sets for 5-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|--|------------|---|
| 989803152061 CE | 5-lead grabber and shower shield, multi-color wire | 1 lead set | Telemetry: 862108: M2601B; 862231: TRx4851A |
| 989803152081 CE | 5-lead snap and shower shield, multi-color wire | 1 lead set | |



989803152061



989803152081

IntelliVue Telemetry 6-Lead Sets – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|-------------------------------|--|------------|--|
| REVISÉD 989803152131 CE | 6-lead grabber and shower shield, multi-color wire | 1 lead set | Telemetry: 862439: U.S. only TRx4841A; 862231: Non-U.S. |
| REVISÉD 989803152151 CE | 6-lead snap and shower shield, multicolor wire | 1 lead set | TRx4851A |



989803152131



989803152151

IntelliVue Telemetry 6-Lead Sets – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|--|------------|---|
| 989803152141 CE | 6-lead grabber and shower shield, multi-color wire | 1 lead set | Telemetry: 862231: Non-U.S. TRx4851A |
| 989803152161 CE | 6-lead snap and shower shield, multicolor wire | 1 lead set | |



989803152141



989803152161

IntelliVue Telemetry Tether Cables

| Product # | Description | Packaging | Use With: |
|--------------------|--|-----------|--|
| 989803143481 CE | PWD tether cable connects either the TRx4841A, TRx4851A or M2601B to TeleMon tether cable or MP5 tether cable. | 1 cable | Telemetry: 862108: M2601B; 862439: U.S. only TRx4841A; 862231: Non-U.S. TRx4851A |
| 989803143491 CE | TeleMon tether cable connects the M2636C to the PWD tether cable. | 1 cable | |
| 989803146911 CE | MP5 tether cable connects the MP5 to the PWD tether cable. | 1 cable | |



989803143481



989803143491



989803146911

IntelliVue Shower Shields & Alignment Guides

| Product # | Description | Packaging | Use With: |
|--------------|--|------------|---|
| 989803153031 | 3 Lead Detachable Shield Replacements | 10 shields | IntelliVue Telemetry lead sets |
| 989803153041 | 5 Lead Detachable Shield Replacements | 10 shields | |
| 989803153051 | 6 Lead Detachable Shield Replacements | 10 shields | |
| 989803140401 | Single ECG alignment guide protective insert, untethered | 10 guides | IntelliVue Telemetry: 862231:TRx4851A Transceiver; 862439: TRx4841A |
| 989803140411 | Single ECG alignment guide protective insert, tethered | 10 guides | Transceiver; 862108: M2601B |
| 989803140421 | Double ECG alignment guide protective insert, untethered | 10 guides | Transmitter |



989803153031



989803153041



989803153051



989803140401 / 989803140411 / 989803140421

Legacy HP/Agilent Telemetry Leads

| Product # | Description | Packaging | Use With: |
|---------------------|--|-----------|--|
| M1421A | Telemetry 3-lead grabber, AAMI | 1 set | M1403A Transmitter |
| M1422A | Telemetry 4-lead snap, AAMI | 1 set | M1403A Transmitter |
| M1423A | Telemetry 4-lead grabber, AAMI | 1 set | M1403A Transmitter |
| M1425A | Telemetry 5-lead grabber, AAMI | 1 set | M1403A Transmitter |
| M1430A | Telemetry 3 lead snap, IEC | 1 set | M1403A Transmitter |
| CE M1431A | Telemetry 3 lead grabber, IEC | 1 set | M1403A Transmitter |
| CE M1433A | Telemetry 4 lead grabber, IEC | 1 set | M1403A Transmitter |
| CE M2590A | Telemetry 3-lead snap, AAMI | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| M2591A | Telemetry 3-lead grabber, AAMI | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| M2592A | Telemetry 5-lead snap, AAMI | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| M2593A | Telemetry 5-lead grabber, AAMI | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| CE M2594A | Telemetry 3-lead snap, IEC | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| CE M2595A | Telemetry 3-lead grabber, IEC | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| CE M2596A | Telemetry 5-lead snap, IEC | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| CE M2597A | Telemetry 5-lead grabber, IEC | 1 set | Telemetry Transmitters: M2601A, M2600A, M2636A |
| M2598A | Telemetry 3-lead combiner shield (U.S. only) | 1 shield | M2601A Transmitter; M2636A TeleMon |
| M2599A | Telemetry 3-lead combiner shield (U.S. only) | 1 shield | M2601A Transmitter; M2636A TeleMon |



M2593A



M2597A

ECG Extender Cables for TeleMon

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|-----------------------------------|
| M4793A | ECG extender trunk cable, 5-lead, AAMI Recommended when TeleMon is to be used at patient's bedside. | 1 cable | M2636A, M2636B, M2636C ECG cables |
| M4794A CE | ECG extender trunk cable, 5-lead, IEC Recommended when TeleMon is to be used at patient's bedside. | 1 cable | |
| M4795A | ECG extender trunk cable, 3-lead, AAMI Recommended when TeleMon is to be used at patient's bedside. | 1 cable | |
| M4796A CE | ECG extender trunk cable, 3-lead, IEC Recommended when TeleMon is to be used at patient's bedside. | 1 cable | |



M4794A


Disposable ECG Electrodes

Disposable Adult ECG Monitoring Electrodes

| Product # | Description | Packaging | Use With: |
|--------------------|---|--|----------------------|
| 989803148801 CE | Small Adult, foam, solid gel snap ECG electrode Size: 1.6" x 1.26". Available in U.S. and Japan. | 30/pack 1500/case Shelf life: 6 months | Adult ECG Monitoring |
| 989803148811 CE | Small Adult, cloth, solid gel snap ECG electrode Size: 1.3" x 1.3" Available in U.S. and Japan. | 30/pack 450/case Shelf life: 6 months | Adult ECG Monitoring |
| 989803148821 CE | Small Adult, radiolucent foam, solid gel snap ECG electrode Size: 1.4" x 1.6" Available in U.S. and Japan. | 30/pack 600/case Shelf life: 6 months | Adult ECG Monitoring |
| 40489E CE | Adult paper tape ECG electrode Recommended for diaphoretic patients Dimensions: 2-3/8" (59mm) diameter | 30/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| 40493D CE | Adult foam, wet gel ECG electrode Dimensions: 2-1/8" (54mm) diameter | 5/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| 40493E CE | Adult foam, wet gel ECG electrode Dimensions: 2-1/8" (54mm) diameter | 30/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |



989803148801



989803148811



989803148821



40489E



40493D

Disposable Adult ECG Monitoring Electrodes (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|---|---------------------------------|
| M4612A CE | Adult foam, solid gel ECG Electrode Teardrop shape, 2" (50mm) diameter | 5/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| M4613A CE | Adult foam, solid gel ECG Electrode Teardrop shape, 2" (50mm) diameter | 30/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| 13941E CE | Adult cloth, wet gel ECG electrode. Contains non-allergenic skin balance gel. Dimensions: 2-1/8" (54mm) diameter | 30/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| 13942E CE | Adult cloth, wet gel ECG electrode. Contains non-allergenic skin balance gel. Dimensions: 2-1/8" (54mm) diameter | 30/pack 300/case Shelf life: 6 months | Adult ECG Monitoring |
| M4606A CE | Skin prep paper | 100 per bag | use to prep skin for electrodes |



M4612A / M4613A



13941E / 13942E



M4606A

Disposable Adult Metallic ECG Lead Wire Electrodes - AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|---|---|----------------------|
| 989803156201 | Preattached Metallic 3-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2 inches x 1.8 inches | 1 set per pack, 100 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |
| 989803156241 | Preattached Metallic 5-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2 inches x 1.8 inches | 1 set per pack, 60 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |



989803156201



989803156241

Disposable Adult Metallic ECG Lead Wire Electrodes - IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|--|---|----------------------|
| 989803156211 CE | Preattached Metallic 3-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm | 1 set per pack, 100 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |
| 989803156251 CE | Preattached Metallic 5-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm | 1 set per pack, 60 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |

Disposable Adult Radiolucent ECG Lead Wire Electrode - AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|--|---|----------------------|
| 989803156221 | Preattached radiolucent 3-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2" x 1.8" | 1 set per pack, 100 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |
| 989803156261 | Preattached radiolucent 5-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2" x 1.8" | 1 set per pack, 60 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |

Disposable Adult Radiolucent ECG Lead Wire Electrode - IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|---|---|----------------------|
| 989803156231 CE | Preattached radiolucent 3-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm | 1 set per pack, 100 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |
| 989803156271 CE | Preattached radiolucent 5-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm | 1 set per pack, 60 sets/box Shelf life: 24 months unopened | Adult ECG Monitoring |

Disposable Pediatric ECG Monitoring Electrodes

| Product # | Description | Packaging | Use With: |
|--------------|--|---|-----------------------------|
| 13950B CE | Disposable cloth electrode with KCl foam gel retaining ring, snap. (Individually wrapped) Diameter: 1.25" (32mm) | 30 per box 10 boxes/case 6 month shelf life | Pediatric/Neonatal Monitors |
| 13951C CE | Disposable, repositionable solid gel electrode, snap. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 1-1/4" (22mm x 33mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |
| 13955C CE | Disposable, repositionable solid gel electrode, snap. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 7/8" (22mm x 22mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |



13950B




13951C



13955C

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|---------------------------------------|---|---|-----------------------------|
| 13952A | Preattached metallic (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |
| REVISÉD 13953A CE | Preattached metallic (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |



13952A

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: AAMI Color Coding


| Product # | Description | Packaging | Use With: |
|---------------------|--|---|-----------------------------|
| 13952B CE | Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |
| 13953B CE | Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |



13952B

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|---|---|-----------------------------|
| 13952D CE | Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |
| 13953D CE | Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |



13953D

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------|--|---|-----------------------------|
| 13952E CE | Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |
| 13953E CE | Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm) | 3 per pack, 300/case 6 month shelf life | Pediatric/Neonatal Monitors |

EEG Electrode (electroencephalogram)

| Product # | Description | Packaging | Use With: |
|--------------|-------------------------------|---------------------|----------------------|
| M1935A CE | Disposable EEG snap electrode | 25/pack 100/case | Adult EEG Monitoring |




M1935A


Invasive Blood Pressure

Disposable Blood Pressure Kits

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M1567A CE | Single line blood pressure kit. Contains one disposable pressure transducer, infusion set, tube, stopcocks, flushing device and continuous flow system. | 20 kits | IBP Modules: M1006A, M1006B, M3012A, M3016A, M8105A, M8102A, M3002A |
| M1568A CE | Dual line blood pressure kits for measuring CVP, ABP and other vascular pressure measurements. Contains one disposable pressure transducer, infusion set, red and blue color coded tubing, stopcocks, flushing device, and continuous flow system. | 20 kits | IBP Modules: M1006A, M1006B, M3012A, M3016A, M8105A, M8102A, M3002A |
| M1634A CE | Reusable adapter cable for M1567A and M1568A kits. | 1 cable | Multi Measurement Modules: M3002A, MP2: M8102A, MP5: M8105A |
| M2271A | Holder for disposable pressure transducer. Holds up to four pressure transducers. Requires M2272A pole mount. | 1 holder | N/A |
| M2272A | IV Pole Mount. Supports one M2271A holder and clamps to IV pole. | 1 mount | N/A |



M1567A



M1568A

Reusable Invasive Blood Pressure Solution

| Product # | Description | Packaging | Use With: |
|----------------|--|--------------|---|
| CPJ840J6 CE | Reusable Pressure Transducer. Used to measure arterial and venous blood pressure. Includes one transducer holder. Cable length: 3m. | 1 transducer | M8102A, M8105A, M1006A, M1006B, M3000A, M3001A, M3002A, M1312A, M3014A, M3016A, 863065, 863066, 863068 |
| CPJ84022 CE | Single use sterile domes | 50 domes | CPJ840J5, CPJ840J6, M8102A, M1006A, M1006B, M3000A, M3001A, M3002A, M1312A, M3014A, M3016A, 863065, 863066, 863068, M8105A |
| CPJ84026 CE | Sterile line set for CPJ840J6 transducer. Not available in all countries. Contact your local sales office. | 10 sets | N/A |
| CPJ84046 CE | Transducer holder | 4 holders | N/A |
| CPJ84447 CE | IV pole mount | 1 mount | CPJ840J5 or CPJ840J6 transducers with instruments: M8102A, M8105A, M1006A, M1006B, M3000A, M3001A, M3002A, M1312A, M3014A, M3016A, 863065, 863066, 863068 |



CPJ840J6

Non-Invasive Blood Pressure

Reusable NIBP Cuffs

Multi-Patient Comfort Cuffs - Single Hose

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1571A CE | Infant, color coded: orange tubing Bladder width: 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10cm to 15cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| M1572A CE | Pediatric, color coded: green tubing Bladder width: 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14cm to 21.5cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, TeleMon: M2636C |
| M1573A CE | Small Adult, color coded: light blue tubing Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5cm to 28cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, TeleMon: M2636C |

Multi-Patient Comfort Cuffs - Single Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1574A CE | Adult, color coded: navy blue tubing Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, TeleMon: M2636C |
| M1575A CE | Large Adult, color coded: burgundy tubing Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34cm to 43cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, TeleMon: M2636C |
| M1576A CE | Thigh, color coded: grey tubing Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42cm to 54cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; Monitors: MP2: M8102A, MP5: M8105A, M8105AT; C3: 862474, 862478, SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |



Multi-Patient Easy Care Cuffs - 1 Hose

| Product # | Description | Packaging | Use With: |
|---------------|--|-----------|---|
| M4552B CE | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A |
| M4552B5 CE | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A |
| M4553B CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4553B5 CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |

Multi-Patient Easy Care Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|---------------|--|-----------|---|
| M4554B CE | Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference = 20.5cm to 28.5cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4554B5 CE | Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference = 20.5cm to 28.5cm | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4555B CE | Adult, navy blue. Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference = 27.5cm to 36cm | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4555B5 CE | Adult, navy blue. Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference = 27.5cm to 36cm | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, Telemon: M2636C |

Multi-Patient Easy Care Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|---------------|---|-----------|---|
| M4556B CE | Adult Extra-Long, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm). Limb circumference = 27.5cm to 36cm. | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4556B5 CE | Adult Extra-Long, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm). Limb circumference = 27.5cm to 36cm. | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, Telemon: M2636C |
| M4557B CE | Large Adult, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm. | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4557B5 CE | Large Adult, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm. | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; ICG Monitor: 862146, C3: 862474, 862478, SureSigns Monitors: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |

Multi-Patient Easy Care Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|---------------|--|-----------|---|
| M4558B CE | Large Adult Extra-Long, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm. | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; C1: 863051, 863052, 863053, 863054, SureSigns Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4558B5 CE | Large Adult Extra-Long, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm. | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; C1: 863051, 863052, 863053, 863054, SureSigns Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4559B CE | Thigh, grey. Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm). Limb circumference = 44cm to 56cm. | 1 cuff | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; C3: 862474, 862478, C1: 863051, 863052, 863053, 863054, SureSigns Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M4559B5 CE | Thigh, grey. Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm). Limb circumference = 44cm to 56cm. | 5 cuffs | Module: M1008B, Multi-Measurement Modules: M3000A, M3001A, M3002A; C3: 862474, 862478, C1: 863051, 863052, 863053, 863054, SureSigns Monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088, MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |



Multi-Patient Easy Care Cuffs - 2 Hose

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|------------------------------------|
| M4562B CE | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm | 1 cuff | For use with Double Hose Monitors. |
| M4563B CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm | 1 cuff | For use with Double Hose Monitors. |
| M4564B CE | Small Adult, light blue. Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference = 20.5cm to 28.5cm | 1 cuff | For use with Double Hose Monitors. |
| M4565B CE | Adult, navy blue. Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference = 27.5cm to 36cm | 1 cuff | For use with Double Hose Monitors. |
| M4566B CE | Adult Extra-Long, navy blue. Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm). Limb circumference = 27.5cm to 36cm. | 1 cuff | For use with Double Hose Monitors. |
| M4567B CE | Large Adult, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm. | 1 cuff | For use with Double Hose Monitors. |
| M4568B CE | Large Adult Extra-Long, burgundy. Bladder width: 6.7" (17cm) Bladder length: 15" (38cm). Limb circumference = 35cm to 45cm.) | 1 cuff | For use with Double Hose Monitors. |
| M4569B CE | Thigh, grey. Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm). Limb circumference = 44cm to 56cm. | 1 cuff | For use with Double Hose Monitors. |



Multi-Patient Traditional Cuffs - 1 Hose

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| 40401A CE | Infant Bladder width = 2.4" (6cm) Bladder length = 4.7" (12cm) Limb circumference = 10cm to 19cm | 1 cuff | SureSigns: 863051, 863052, 863053, 863054, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; C3 Monitors: 862474, 862478; C1 Monitors: 863051, 863052, 863053, 863054; A1 Monitors: M3921A to M3924A; A3 Monitors: M3926A to M3929A |
| 40401B CE | Pediatric Bladder width = 3.5" (8.9cm) Bladder length = 7.1" (18cm) Limb circumference = 18cm to 26cm | 1 cuff | SureSigns: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; C3 Monitors: 862474, 862478; C1 Monitors: 863051, 863052, 863054; A1 Monitors: M3921A to M3924A; A3 Monitors: M3926A to M3929A; TeleMon: M2636C |
| 40401C CE | Adult Bladder width = 4.8" (12.1cm) Bladder length = 9.1" (23cm) Limb circumference = 25cm to 35cm | 1 cuff | SureSigns: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; C3 Monitors: 862474, 862478; C1 Monitors: 863051, 863052, 863053, 863054; A1 Monitors: M3921A to M3924A; A3 Monitors: M3926A to M3929A; TeleMon: M2636C |

Multi-Patient Traditional Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 40401D CE | Large Adult Bladder width = 6" (15.2cm) Bladder length = 13" (33cm) Limb circumference = 33cm to 47cm | 1 cuff | SureSigns: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; C3 Monitors: 862474, 862478; C1 Monitors: 863051, 863052, 863053, 863054; A1 Monitors: M3921A to M3924A; A3 Monitors: M3926A to M3929A; TeleMon: M2636C |
| 40401E CE | Thigh Bladder width = 7" (17.8cm) Bladder length = 14.2" (36cm) Limb circumference = 46cm to 66cm | 1 cuff | SureSigns: 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; C3 Monitors: 862474, 862478; C1 Monitors: 863051, 863052, 863053, 863054; A1 Monitors: M3921A to M3924A; A3 Monitors: M3926A to M3929A |



Multi-Patient NIBP Cuff Kits - 1 Hose

| Product # | Description | Packaging | Use With: |
|--------------|---|--|---|
| 864288 CE | Easy Care Adult Kit - 4 sizes Small Adult, Adult, Large Adult, Thigh | consists of 4 cuffs: M4554B, M4555B, M4557B, M4559B | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| 864289 CE | Easy Care Pediatric Kit - 4 sizes Infant, Pediatric, Small Adult, Adult | consists of 4 cuffs: M4552B, M4553B, M4554B, M4555B | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| 864290 CE | Easy Care Assortment Kit - 6 sizes Infant, Pediatric, Small Adult, Adult, Large Adult, Thigh | consists of 6 cuffs: M4552B, M4553B, M4554B, M4555B, M4557B, M4559B | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| 864291 CE | Easy Care Adult Long Kit - 6 sizes Small Adult, Adult, Adult Long, Large Adult, Large Adult Long, Thigh | consists of 6 cuffs: M4554B, M4555B, M4556B, M4557B, M4558B, M4559B | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| M1577A CE | Reusable NIBP Comfort Cuff Kit Adult, Infant, Pediatric, Small Adult | consists of 4 cuffs: M1571A, M1572A, M1573A, M1574A | M1008B, MMS: M3000A, M3001A, M3002A; MP2: M8105A; MP5: M8105A; M8105AT; C3: 862474, 862478; SureSigns Monitors: 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082 |
| M1578A CE | Reusable NIBP Comfort Cuff Kit Adult, Large Adult, Small Adult, Thigh | consists of 4 cuffs: M1573A, M1574A, M1575A, M1576A | M1008B, MMS: M3000A, M3001A, M3002A; MP2: M8105A; MP5: M8105A; M8105AT; C3: 862474, 862478; SureSigns Monitors: 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082 |

Multi-Patient NIBP Cuff Kits - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|--|---|
| M1579A CE | Reusable NIBP Comfort Cuff Kit Small Adult, Large Adult, Adult, Infant, Pediatric, Thigh | consists of 6 cuffs: M1571A, M1572A, M1573A, M1574A, M1575A, M1576A | M11008B, MMS: M3000A, M3001A, M3002A; MP2: M8105A; MP5: M8105A; M8105AT; C3: 862474, 862478; SureSigns Monitors: 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082 |
| 40400A CE | Traditional Reusable NIBP Cuff Kit Pediatric, Adult, Large Adult | consists of 3 cuffs: 40401B, 40401C, 40401D | SureSigns: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| 40400B CE | Traditional Reusable NIBP Cuff Kit Infant, Pediatric, Adult, Large Adult, Thigh | consists of 5 cuffs: 40401A, 40401B, 40401C, 40401D, 40401E | SureSigns: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |



M1577A



M1578A

Disposable NIBP Cuffs

Single Patient Disposable Cuffs

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M1874A CE | Infant (orange) Bladder width = 2.2" (5.5cm) Bladder length = 6.1" (15.4cm) Limb circumference = 10cm to 15cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A, MP5 M8105A, M8105AT; SureSigns Monitors: 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |
| M1875A CE | Pediatric (green) Bladder width = 3.1" (8cm) Bladder length = 8.5" (21.5cm) Limb circumference = 14cm to 21.5cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A, MP5 M8105A, M8105AT; SureSigns Monitors: 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |
| M1876A CE | Small Adult (light blue) Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference = 20.5cm to 28cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A, MP5 M8105A, M8105AT; SureSigns Monitors: 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |
| M1877A CE | Adult (navy blue) Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27cm to 35cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A; MP5 M8105A, M8105AT; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |

Single Patient Disposable Cuffs (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M1878A CE | Large Adult (burgundy) Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference = 34cm to 43cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A, MP5 M8105A, M8105AT; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |
| M1879A CE | Thigh (grey) Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm). Limb circumference: 42cm to 54cm | 10 cuffs | Modules: M1008B, Multi-Measurement Module: M3000A, M3001A, M3002A; Patient Monitors: C3: 862474, 862478; MP2 M8102A, MP5 M8105A, M8105AT; SureSigns Monitors: 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082 |



Soft Single Patient Disposable Cuffs - 1 Hose

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M4572B CE | Infant (orange) Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm | 10 cuffs | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863061, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A |
| M4573B CE | Pediatric (green) Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm | 10 cuffs | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |
| M4574B CE | Small Adult (light blue) Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference = 20.5cm to 28.5cm | 10 cuffs | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |
| M4575B CE | Adult (navy blue) Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference = 27.5cm to 36cm | 10 cuffs | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |

Soft Single Patient Disposable Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M4576B CE | Adult Extra-Long (navy blue) Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm. | 10 cuffs | M11008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |



Soft Single Patient Disposable Cuffs - 1 Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M4577B CE | Large Adult (burgundy) Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35cm to 45cm. | 10 cuffs | M11008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |
| M4578B CE | Large Adult Extra-Long (burgundy) Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35cm to 45cm. | 10 cuffs | M11008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A; TeleMon: M2636C |
| M4579B CE | Thigh (grey) Bladder width: 8.3" (21cm) Bladder length 18.9" (48cm) Limb circumference = 44cm to 56cm | 10 cuffs | M11008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863051, 863052, 863053, 863055, 853056, 863057, 863059, 863060, 863061, 863062, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A |

Soft Single Patient Disposable Cuffs - 2 Hose

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|------------------------------------|
| M4582B CE | Infant (orange) Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10-15cm | 10 cuffs | For use with Double Hose Monitors. |
| M4583B CE | Pediatric (green) Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14-21.5cm | 10 cuffs | For use with Double Hose Monitors. |
| M4584B CE | Small Adult (light blue) Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm). Limb circumference = 20.5cm to 28.5cm | 10 cuffs | For use with Double Hose Monitors. |
| M4585B CE | Adult (navy blue) Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm). Limb circumference = 27.5cm to 36cm | 10 cuffs | For use with Double Hose Monitors. |



Soft Single Patient Disposable Cuffs - 2 Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|------------------------------------|
| M4586B CE | Adult Extra-Long (navy blue) Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5cm to 36cm. | 10 cuffs | For use with Double Hose Monitors. |
| M4587B CE | Large Adult (burgundy) Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35cm to 45cm. | 10 cuffs | For use with Double Hose Monitors. |
| M4588B CE | Large Adult Extra-Long (burgundy) Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35cm to 45cm. | 10 cuffs | For use with Double Hose Monitors. |
| M4589B CE | Thigh (grey) Bladder width: 8.3" (21cm) Bladder length 18.9" (48cm) Limb circumference = 44cm to 56cm | 10 cuffs | For use with Double Hose Monitors. |

NIBP Accessories and Hoses

Blood Pressure Air Hoses

| Product # | Description | Packaging | Use With: |
|--------------------|---|------------|---|
| 989803136931 CE | Adult/Pediatric NIBP air hose Length: 11.71' (3.57m) | 1 air hose | VS1 monitors: 863055, 863056, 863057, 863058, 863059, 863060, 863061 |
| M1598B CE | Adult NIBP air hose Length: 4.92' (1.5m) | 1 air hose | M1008B, MMS: M3000A, M3001A, M3002A; SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M1599B CE | Adult NIBP air hose length: 10' (3.0m) | 1 air hose | M1008B, MMS: M3000A, M3001A, M3002A; C3: 862474, 862478, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 863085, 863086, 863087, 863088; MP5: M8105A, M8105AT, MP2: M8102A, TeleMon: M2636C |
| M3918A CE | Adult/Pediatric air hose Length: 10' (3.0m) | 1 air hose | A1: M3921A, M3922A, M3923A, M3924A, A3: M3926A, M3927A, M3928A, M3929A; ICG: 862146; C1 Monitors: 863051, 863052, 863053, 863054 |



989803136931



M1599B



M3918A

NIBP Cuff Caps

| Product # | Description | Packaging | Use With: |
|--------------------|---------------------------------|-----------|-----------|
| 989803151601 CE | Cuff Caps for Single Tube Cuffs | 10 caps | n/a |
| 989803151611 CE | Cuff Caps for Double Tube Cuffs | 20 caps | n/a |



989803151601



989803151611

Bulb and Valve Kit for Manual Sphygmomanometer

| Product # | Description | Packaging | Use With: |
|--------------|---|------------------|---|
| M4580A CE | Bulb and Valve Kit Manual Sphygmomanometer. Replacement coil tubing (M4581A) is available for this product. | 1 bulb and valve | 2-hose reusable cuff models M456xB and soft single patient disposable models M458xB only. |
| M4581A CE | Coiled Tubing Bulb and Valve Kit, Replacement, length: 8' | 1 tubing | M4580A Bulb and Valve Kit |



M4580A



M4581A

NIBP Cuff Adapters

| Product # | Description | Packaging | Use With: |
|----------------------|---|-----------|---|
| REVISÉD M4550A CE | 1 Hose adapter (Welch-Allyn) Threaded female style connector | 1 adapter | Philips 1-hose reusable and disposable cuffs excluding 40401x series. |
| REVISÉD M4551A CE | 2 Hose Y-adapter (Welch-Allyn) Y style cuff adapter | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| REVISÉD M4560A CE | 1 Hose adapter Female quick-connect style connector | 1 adapter | Philips 1-hose reusable and disposable cuffs excluding 40401x series. |
| REVISÉD M4561A CE | 2 Hose Y-adapter Female quick-connect style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| M4570A | 1 Hose adapter (Draeger) Luer-lock style cuff adapter | 1 adapter | Philips 1-hose reusable and disposable cuffs excluding 40401x series. |
| M4571A | 2 Hose Y-Adapter (Draeger) Y Style | 1 adapter | Philips 2-hose reusable and disposable cuffs. |



M4550A



M4551A



M4560A



M4561A



M4570A



M4571A

NIBP Cuff Adapters (continued)

| Product # | Description | Packaging | Use With: |
|----------------------|--|-----------|---|
| M4590A | 1 Hose Adapter (Colin, Datascope-Passport, Spacelabs, Welch-Allyn, Mindray) Locking luer style connection | 1 adapter | Philips 1-hose reusable and disposable cuffs excluding 40401x series. |
| REVISÉD M4591A CE | 1 Hose Adapter (Welch-Allyn Protocol) Male quick-connect style connector | 1 adapter | Philips 1-hose reusable and disposable cuffs excluding 40401x series. |
| M4592A | 2 Hose Adapter (GE/Marquette) Male quick connect style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| M4593A | 2 Hose Adapter (Datex/Ohmeda) One male and one female quick-connect style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| REVISÉD M4594A CE | Y Adapter (Spacelabs, Siemens, Mindray) Bayonet style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| REVISÉD M4595A CE | Y Adapter (Colin, Datascope-Passport, Spacelabs, Welch-Allyn, Mindray) Locking luer style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |
| REVISÉD M4596A CE | Y Adapter (Welch-Allyn Protocol) Male quick-connect style connector | 1 adapter | Philips 2-hose reusable and disposable cuffs. |



M4590A



M4591A



M4592A



M4593A



M4594A



M4595A



M4596A

Paper

Recording Paper

| Product # | Description | Packaging | Use With: |
|--------------|--|---|--|
| 40477A CE | Thermal array recorder paper | 10 rolls | M1116A/M1116B Recorder, M3176B Recorder Module: MP 20 Junior, MP20, MP30, MP40, MP50, MP60, MP70, MP90, M3001A |
| 40477B CE | Thermal array recorder paper | 80 rolls | M1116A/M1116B Recorder, M3176B Recorder Module: MP 20 Junior, MP20, MP30, MP40, MP50, MP60, MP70, MP90, M3001A |
| M4816A CE | Thermal array recorder paper | 10 rolls | M3925A Recorder for A1/C1 Compact Patient Monitors; M3176C USB 2-Channel Recorder; A3/C3 Recorder |
| M4817A CE | Thermal array recorder paper | 80 rolls | M3925A Recorder for A1/C1 Compact Patient Monitors; M3176C USB 2-Channel Recorder; A3/C3 Recorder |
| PSE11268 | Seiko recorder paper, 4-channel | 24 rolls | M3160A Recorder |
| 40469A | Multi-channel thermal array recorder, z-fold | 1 case = 20 packs (300 sheets per pack) | M1117A Recorder |
| M1885A CE | Thermal array recorder paper, A4 size | 1 case = 6 packs (500 sheets per pack) | M1118A Recorder |
| 989803138131 | ICG Paper | 20 rolls | 862146 ICG integrated recorder |
| 989803136891 | Thermal array recorder paper | 5 rolls | VS1, VS3, VM4, VM6, VM8 integrated recorder |



40477A



989803138131



989803136891

Pulse Oximetry

Philips Reusable SpO₂ Sensors

| Product # | Description | Packaging | Use With: |
|---------------|--|-----------|--|
| M1191B CE | Adult Finger Glove Sensor. Length = 2m Patient size: >110 lb (>50kg) 8 pin connector | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1191BL CE | Adult Finger Glove Sensor. Length = 3m Patient size: >110 lb (>50kg) 8 pin connector | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |



M1191B

Philips Reusable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1191T CE | Adult Finger Glove Sensor: Length = 0.45m Patient size: >110 lb (>50kg) 9 pin connector | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1192A CE | Pediatric/Small-Adult Glove Sensor: Length = 1.5m Patient size: 33-110 lb (15-50kg) 8 pin connector | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |



M1191T



M1192A

Philips Reusable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1192T CE | Pediatric/Small-Adult Glove Sensor. Length = 0.45m Patient size: 33-110 lb (15-50kg). 9 pin connector. | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1193A CE | Reusable Neonatal Wrap Sensor. Length = 1.5m Neonatal foot/hand. Patient size: 2.2 lbs. to 8.8 lbs. (1-4kg) 8 pin connector | 1 sensor | Monitors: SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |



M1192T



M1193A

Philips Reusable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1193T CE | Reusable Neonatal Wrap Sensor. Length = 0.90m Neonatal foot/hand. Patient size: 2.2 lbs. to 8.8 lbs. (1-4kg) 9 pin connector | 1 sensor | Monitors: SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1194A CE | Pediatric/Adult Ear Clip Sensor. Length = 1.5m Patient size: >88 lb (>40kg). 8 pin connector. | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |



M1193T



M1194A

Philips Reusable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1195A CE | Reusable Infant Finger Glove Sensor. Length = 1.5m Infant finger. Patient size: 8.8 lbs. to 33 lbs. (4-15kg) 8 pin connector | 1 sensor | Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1196A CE | Pediatric/Adult Finger Clip Sensor. Length = 3m Patient size: >88 lb (>40kg). 8 pin connector. | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1196T CE | Pediatric/Adult Finger Clip Sensor. Length = 90cm Patient size: >88 lb (>40kg). 9 pin connector. | 1 sensor | Telemetry: M2600B telemetry system, M2601B transmitter, TRx4841A, TRx4851A transceivers; MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A; Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1627A CE | Wristband Velcro straps for M1191A/B SpO ₂ sensors | | n/a |



M1195A



M1196A



M1196T



M1627A

Philips Disposable SpO₂ Sensors

| Product # | Description | Packaging | Use With: |
|--------------|--|------------|--|
| M1131A CE | Adult/Pediatric Finger Clip Disposable Sensor (single patient use). Length = 45cm. Patient size: >44 lb. (>20kg) 9 pin connector | 20 sensors | Telemetry: M2601A transmitter, MP5: M8105AT. Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A with SpO ₂ Option A01 |
| M1132A CE | Infant Disposable Wrap Sensor (single patient use). Length = 0.90m. Infant finger. Patient size: 6.6 lbs – 22 lbs (3-10kg) 9 pin connector | 20 sensors | Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A Telemetry: M8105AT with SpO ₂ Opt. A01 |
| M1133A CE | Disposable Wrap Sensor (single patient use). Length = 90cm. Patient sizes: Adult (finger) >88 lbs (>40kg) Neonate: (foot/hand) <6.6 lbs (<3kg) Infant: Infant: (big toe/thumb) 22-44 lbs (10-20kg) | 20 sensors | Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A Telemetry: M8105AT with SpO ₂ Opt. A01 |

Philips Disposable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--|------------|---|
| M1134A CE | Adhesive-free Neonatal/Infant/Adult Disposable Wrap Sensor (single patient use). Length = 90cm. Patient sizes: Adult: >88 lb (>40kg) Neonate: < 6.6 lbs (<3kg) Infant: 22-44 lbs (10-20kg) 9 pin connector | 20 sensors | Monitors: C3: 862474, 862478; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863264, 863265, 863266; M3921A, M3922A, M3923A, M3924A, M3926A, M3927A, M3928A, M3929A, Multi-Measurement Modules: M3000A, M3001A, M3002A; SpO ₂ Modules: M1020A, M1020B; MP2: M8102A, MP5: M8105A Telemetry: M8105AT with SpO ₂ Opt. A01 |



M1131A



M1132A



M1133A



M1134A

SpO₂ Adapter Cables for Philips Sensors

| Product # | Description | Packaging | Use With: |
|---------------|---|-----------|--|
| M1900B CE | SpO ₂ Adapter Cable (12-pin) Adapts 9-pin sensors to 12-pin sockets Length: 9.8 feet (3m) | 1 cable | SpO ₂ Module: M1020A |
| M1940A CE | SpO ₂ Adapter Cable (12-pin) Adapts 8-pin sensors to 12-pin sockets Length: 6.6 feet (2m) | 1 cable | SpO ₂ Modules: M1020A |
| M1941A CE | SpO ₂ Extension Cable (8-pin) Adds additional length to 8-pin sensors to 8-pin sockets Length: 6.6 feet (2m) | 1 cable | SpO ₂ Modules: M1020B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A, M8105AT with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266 |
| M1943A CE | SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 3.3 feet (1.1m) | 1 cable | SpO ₂ Modules: M1020B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A, M8105AT with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266 |
| M1943AL CE | SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 9.8 feet (3m) | 1 cable | SpO ₂ Modules: M1020B, Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A, M8105AT with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266 |
| M1943NL CE | OxiMax® SpO ₂ Adapter Cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 9.8 feet (3m) | 1 cable | M1020B-A02, M3001A-A02, M3001A-A04, M3002A-A04, M8102A-SP4, M8105A-SP4, M8105AT-SP4 |
| M4787A CE | SpO ₂ Adapter Cable (9-pin) Adds additional length to 9-pin sensors to 9-pin sockets Length: 9.8 feet (3m) | 1 cable | A1: M3921A, M3922A, M3923A, M3924A, M3925A; A3: M3926A, M3927A, M3928A, M3929A; C1: 863051, 863053, 863054 |



M1900B



M1940A



M1941A



M1943A/AL



M1943NL

Adhesive SpO₂ Sensors (available from Philips in Europe and Japan only)

| Product # | Description | Packaging | Use With: |
|--------------|---|--|--|
| M1901B CE | Neonatal Disposable Sensor, sterile Patient size: <7 lbs. or >88 lbs. (<3kg or >40kg) 9 pin connector | 24 sensors and tapes Shelf life: 1 year | SpO ₂ Modules: M1020B Opt A01 or A02; Multi-Measurement Modules: M3000A Opt. A01, M3001A Opt. A01, A02 or A04; M3002A Opt. A01 or A04; MP5 Monitor: M8105A/AT Opt. A01 or Opt. SP4; M8102A Opt. A01 or Opt. SP4; Suresigns Monitors: 863051, 863052, 863053, 863054, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088 |
| M1902B CE | Infant Disposable Sensor, sterile Patient size: 7-44 lbs (3-20kg) 9 pin connector | 24 sensors and tapes Shelf life: 1 year | |
| M1903B CE | Pediatric Disposable Sensor, sterile Patient size: 22-110 lbs. (10-50kg) 9 pin connector | 24 sensors and tapes Shelf life: 1 year | |
| M1904B CE | Adult Disposable Sensor, sterile Patient size: >66 lbs (>30kg) 9 pin connector | 24 sensors and tapes Shelf life: 1 year | |



M1901B



M1902B



M1903B



M1904B

Masimo® LNOP® SpO₂ Adapter Cables

| Product # | Description | Packaging | Use With: |
|-------------------|---|-----------|---|
| M1020-61100 CE | Masimo LNOP Adapter Cable Adapts LNOP sensors to Philips monitors. Length: 11.8 feet (3.6m) | 1 cable | SpO ₂ Module: M1020B; Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A with SpO ₂ Opt. A01; SureSigns: 863063, 863064, 863065, 863066, 863068, 863080, 863081, 863082, 863071, 863072, 863073, 863074, 863077, 863266, 863265, 863266 SpO ₂ SET® Module: M1020B Opt. A03; Multi-Measurement SET® Module: M3001A Opt. A03, M3002A Opt. A03 SpO ₂ Module: M1020A |
| M1020-61102 CE | Masimo LNOP Adapter Cable Adapts LNOP sensors to Philips Viridia monitors. Length: 11.8 feet (3.6m) | 1 cable | SpO ₂ Module: M1020A |



M1020-61100



M1020-61102

Masimo® LNCS® SpO₂ Adapter Cable

| Product # | Description | Packaging | Use With: |
|--------------------|---|-----------|---|
| 989803148221 CE | Masimo LNCS Adapter Cable Adapts LNCS sensors to Philips monitors. Length: 10 feet (3m) | 1 cable | SpO ₂ Module: M1020A, M1020B; Multi-Measurement Modules: M3000A, M3001A, M3002A; MP2: M8102A; MP5: M8105A with SpO ₂ Opt. A01; SpO ₂ SET® Module: M1020B Opt. A03; Multi-Measurement SET® Module: M3001A Opt. A03, M3002A Opt. A03 |



989803148221

Temperature Monitoring

Disposable Skin Surface Temperature Probe

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| 21091A CE | Disposable sterile skin surface temperature probe Length: 30 inches (76cm) Shelf life: at least 12 months | 10 probes | Temperature Module: M1029A Multi-Measurement Modules: M3000A, M3001A, M3002A SureSigns: 863064, 863065, 863066, 863068, 863087, 863088 MP2 Monitor: M8102A, MP5 Monitor: M8105A |



21091A

Disposable Esophageal/Rectal Temperature Probes (meets EN12470 Standards)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M1837A CE | Esophageal/Rectal Temperature Probe, 9 FR Length: 36 inches (91cm), sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21090A CE | Esophageal/Rectal Temperature Probe, 12FR Length: 36 inches (91cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |



M1837A



21090A

Disposable Esophageal/Stethoscope Temperature Probes (meets EN12470 Standards)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| 21093A CE | Esophageal/Stethoscope Temperature Probe, 12 FR Length: 20 inches (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21094A CE | Esophageal/Stethoscope Temperature Probe, 18FR Length: 20 inches (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21095A CE | Esophageal/Stethoscope Temperature Probe, 24 FR Length: 20 inches (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |



21093A



21094A

Disposable Foley Catheter Temperature Probes (meets EN12470 Standards)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| M2255A CE | Foley Catheter Temperature Probe, 14FR Length: 17 inches (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21096A CE | Foley Catheter Temperature Probe, 16FR Length: 17 inches (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21097A CE | Foley Catheter Temperature Probe, 18FR Length: 17 inches (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |



M2255A



21096A



21097A

Disposable Temperature Probes Extension and Adapter Cables

| Product # | Description | Packaging | Use With: |
|--------------------|--|-----------|---|
| 21082A CE | Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 10 feet (3m) | 1 cable | Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21082B CE | Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 5 feet (1.5m) | 1 cable | Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 989803162591 CE | Disposable temperature probe extension cable for Series 400 monitors Length: 10 feet (3m) | 1 cable | Series 400 Temperature Monitors with 1/4-inch phone plug |
| 989803162601 CE | Disposable temperature probe extension cable for Series 400 monitors Length: 5 feet (1.5m) | 1 cable | Series 400 Temperature Monitors with 1/4-inch phone plug |



21082A



989803162591

Reusable Autoclavable Temperature Probes for Continuous Monitoring (meets EN12470 Standard)

| Product # | Description | Packaging | Use With: |
|--------------------|--|-----------|--|
| 21075A CE | Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 2-prong plug, reusable, continuous monitoring Length: 10 feet (3m) | 1 probe | Temperature Module: M1029A, C3: 862474, 862478, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21076A CE | Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 2-prong plug, reusable, continuous monitoring Length: 10 feet (3m) | 1 probe | Temperature Module: M1029A, MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 21078A CE | Skin Surface Temperature Probe (autoclavable) 0.37 inches (9.5mm) diameter 2-prong plug, reusable tape-on, continuous monitoring Length: 10 feet (3m) | 1 probe | Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A |
| 989803162621 CE | Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 1/4-inch phone plug, reusable, continuous monitoring Length: 10 feet (3m) | 1 probe | Series 400 Temperature Monitors with 1/4-inch phone plug. |



21075A



21076A



989803162621

Reusable Autoclavable Temperature Probes for Continuous Monitoring (meets EN12470 Standard) (continued)

| Product # | Description | Packaging | Use With: |
|--------------------|--|-----------|--|
| 989803162631 CE | Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 1/4-inch phone plug, reusable, continuous monitoring Length: 10 feet (3m) | 1 probe | Series 400 Temperature Monitors with 1/4-inch phone plug. |
| 989803162641 CE | Skin Surface Temperature Probe (autoclavable) 0.37 inches (9.5mm) diameter 1/4-inch phone plug, reusable tape-on, continuous monitoring Length: 10 feet (3m) | 1 probe | Series 400 Temperature Monitors with 1/4-inch phone plug. |
| 989803163281 CE | Adapter cable required to connect reusable temp probes with 1/4 inch phone plug to legacy HP/Agilent and Philip 2-pin input connectors. | 1 probe | Use with temperature probes 989803162621, 989803162631, 989803162641 |
| 989803163271 CE | Adapter cable required to connect reusable temp probes with 2-pin connector to Series 400 temperature monitors with 1/4-inch phone plug input connector. | 1 probe | Use with temperature probes 21075A, 21076A, 21078A |



989803162631



989803162641



989803163281



989803163271

Reusable Temperature Probes for SureSigns VS1 Monitors

| Product # | Description | Packaging | Use With: |
|--|---|------------|--|
| 989803136921 CE | Temperature probe covers, disposable | 200 covers | 989803136901, 989803136911 temperature probes |
| 989803136901 CE | Temperature Probe - Oral reusable, predictive temp monitoring | 1 probe | VS1, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062 |
| REVISÉD 989803136911 CE LATEX | Temperatue Probe - Rectal reusable, predictive temp monitoring | 1 probe | VS1, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863062 |



989803136921



989803136901



989803136911

Reusable Temperature Probes for SureSigns VM4/VS3/MP5 Monitors

| Product # | Description | Packaging | Use With: |
|--------------------|---|--------------|---|
| 989803143381 CE | Reusable Oral Temp Probe & Well Kit predictive temp monitoring | 1 probe | SureSigns VS3: 863070, 863073, 863074; VM4: 863063; IntelliVue MP5: M8105A, M8105AT |
| 989803143391 CE | Reusable Rectal Temp Probe & Well Kit predictive temp monitoring | 1 probe | SureSigns VS3: 863070, 863073, 863074; VM4: 863063; IntelliVue MP5: M8105A, M8105AT |
| M4823A CE | Temperature Probe Covers | 1,000 covers | 989803143381, 989803143391 temperature probes |



989803143381



989803143391

Reusable Temperature Probes for SureSigns A1/A3/C1/VM4-8 Series Monitors

| Product # | Description | Packaging | Use With: |
|--------------|-----------------------------------|-----------|---|
| M4821A CE | Reusable Oral Temperature Probe | 1 probe | A1/C1 monitors, M3922A, M3924A, 863052, 863054, 863051, 863053; VM4/VM6/VM8 monitors 863063, 863064, 863065, 863066, 863068 |
| M4822A CE | Reusable Rectal Temperature Probe | 1 probe | A1/C1 series monitors, M3922A, M3924A, 863052, 863054, 863051, 863053 |



M4821A



M4822A

Batteries, Cases, Pouches, Sethoscope Parts

Batteries

Patient Monitoring Batteries

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| 989803144631 | Lithium Ion battery | 1 battery | SureSigns: 863063, 863064, 863065, 863066, 863067, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863077, 863079, 863080, 863081, 863082, 863085, 863086, 863088 |
| M4605A | Rechargeable Lithium Ion Battery, 10.8 volts | 1 battery | IntelliVue MP20, MP30, MP40, MP50, MP5, MP5T |
| M4607A | Rechargeable Lithium Ion Battery, 10.8 volts | 1 battery | IntelliVue MP2: M8102A, Multi-Measurement Module X2: M3002A |
| M4607-60001 | Battery Charger Adapter | 1 adapter | IntelliVue Multi-Measurement Module X2: M3002A, MP2: M8102A |



989803144631



M4605A



M4607A

Telemetry Batteries

| Product # | Description | Packaging | Use With: |
|-----------|---|--------------|------------------------|
| 40455A | Zinc air battery, 8.4 volts Shelf life: 9 months | 12 batteries | 78100A, M1400A, M2600A |
| ULBU9VLJ | Lithium battery, 9 volts Shelf life: 10 month | 10 batteries | M2600A, M2601A |
| M4790A | Rechargeable battery pack | 1 combiner | M2636A |

Cases, Pouches, Bags

Patient Monitoring Cases, Bags, Pads

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 989803163331 | MP2/X2 Carry case standard, red | 1 case | M3002A, M8102A |
| 989803163341 | MP2/X2 Carry case standard, black | 1 case | M3002A, M8102A |
| 989803163351 | MP2/X2 Carry case (mini), red | 1 case | M3002A, M8102A |
| 989803163361 | MP2/X2 Carry case (mini), black | 1 case | M3002A, M8102A |
| 989803163371 | MP2/X2 Carry case (large), black | 1 case | M3002A, M8102A |
| 989803163631 | MP2/X2 Carry case replacement kit Kit contains replacement parts for MP2/X2 carry cases: 2 replacement velcro straps, 1 pouch for SpO ₂ sensors, 1 shoulder strap, 1 handle strap, and 1 cable sheath. | 1 kit | M3002A, M8102A |
| 989803153441 | MP2/X2 Carrying strap | 1 strap | M8102A MP2, M3002A Multi-Measurement X2 |
| 989803157171 | X2/MP2 Anti-slip pad | 5 pads | M8102A MP2, M3002A Multi-Measurement X2 |
| 989803163381 | NGV display filter for the MP2/X2 display for night time use. | | M3002A, M8102A |
| M3919A | Canvas tote bag | 1 bag | A3 Monitor: M3923A, M3924A, M3925A, M3926A, M3927A, M3928A, M3929A |
| M3920A | Snap-on accessory bag for the A3 canvas tote bag | 1 bag | M3919A Bag |



989803163331



989803163341



989803163351



989803163361



989803163371

Telemetry Pouches and Protective Covers

| Product # | Description | Packaging | Use With: |
|---------------|--|-----------------------|--|
| 9300-0768-050 | Disposable water and chemical resistant telemetry pouch with snaps | 1 box of 50 pouches | 78100A, M1400A, M2600A |
| 9300-0768-200 | Disposable water and chemical resistant telemetry pouch with snaps | 4 boxes of 50 pouches | 78100A, M1400A, M2600A |
| 989803137831 | Transceiver pouch shoulder strap and viewing window. | 1 box of 50 pouches | IntelliVue Telemetry: 862231 (TRx4851A) 862439 (TRx4841A) 862108 (M2601B) |
| 989803140371 | Transceiver pouch shoulder strap and viewing window. | 4 boxes of 50 pouches | |
| 989803140431 | TeleMon service port protective cover for TRx (ECG only) | 10 covers | |
| 989803140451 | TeleMon service port protective cover for TRx+ (ECG and SpO ₂) | 10 covers | |
| 989803140441 | SpO ₂ connector protective cover | 10 covers | |



989803137831



40431



40441



40451

Index

| | | | | | |
|-------------|--------|--------------------|-----|-------------------|-----|
| 13901A..... | 11, 18 | 40401E..... | 65 | 989803148801..... | 47 |
| 13902A..... | 11, 18 | 40455A..... | 100 | 989803148811..... | 47 |
| 13904A..... | 11, 18 | 40469A..... | 79 | 989803148821..... | 47 |
| 13905A..... | 11, 18 | 40477A..... | 79 | 989803148841..... | 39 |
| 13907A..... | 11 | 40477B..... | 79 | 989803148861..... | 39 |
| 13941E..... | 48 | 40489E..... | 47 | 989803151601..... | 76 |
| 13942E..... | 48 | 40493D..... | 47 | 989803151611..... | 76 |
| 13950B..... | 50 | 40493E..... | 47 | 989803151971..... | 40 |
| 13951C..... | 50 | | | 989803151981..... | 40 |
| 13952A..... | 51 | 864288..... | 66 | 989803151991..... | 40 |
| 13952B..... | 51 | 864289..... | 66 | 989803152001..... | 40 |
| 13952D..... | 52 | 864290..... | 66 | 989803152051..... | 41 |
| 13952E..... | 52 | 864291..... | 66 | 989803152061..... | 41 |
| 13953A..... | 51 | | | 989803152071..... | 41 |
| 13953B..... | 51 | 9300-0768-050..... | 102 | 989803152081..... | 41 |
| 13953D..... | 52 | 9300-0768-200..... | 102 | 989803152131..... | 42 |
| 13953E..... | 52 | | | 989803152141..... | 42 |
| 13955C..... | 50 | 989803136891..... | 79 | 989803152151..... | 42 |
| 14363A..... | 20 | 989803136901..... | 98 | 989803152161..... | 42 |
| 15244A..... | 22 | 989803136911..... | 98 | 989803153031..... | 44 |
| | | 989803136921..... | 98 | 989803153041..... | 44 |
| 21075A..... | 96 | 989803136931..... | 75 | 989803153051..... | 44 |
| 21076A..... | 96 | 989803137831..... | 102 | 989803153441..... | 101 |
| 21078A..... | 96 | 989803138131..... | 79 | 989803156201..... | 49 |
| 21082A..... | 95 | 989803140371..... | 102 | 989803156211..... | 49 |
| 21082B..... | 95 | 989803140401..... | 44 | 989803156221..... | 50 |
| 21090A..... | 92 | 989803140411..... | 44 | 989803156231..... | 50 |
| 21091A..... | 92 | 989803140421..... | 44 | 989803156241..... | 49 |
| 21093A..... | 93 | 989803140431..... | 102 | 989803156251..... | 49 |
| 21094A..... | 93 | 989803140441..... | 102 | 989803156261..... | 50 |
| 21095A..... | 93 | 989803140451..... | 102 | 989803156271..... | 50 |
| 21096A..... | 94 | 989803143171..... | 34 | 989803157171..... | 101 |
| 21097A..... | 94 | 989803143181..... | 34 | 989803159571..... | 13 |
| 23001A..... | 22 | 989803143191..... | 34 | 989803159581..... | 13 |
| 23001B..... | 22 | 989803143201..... | 34 | 989803160241..... | 13 |
| 23002A..... | 22 | 989803143381..... | 98 | 989803160251..... | 13 |
| | | 989803143391..... | 98 | 989803160261..... | 13 |
| 40400A..... | 67 | 989803143481..... | 43 | 989803160271..... | 13 |
| 40400B..... | 67 | 989803143491..... | 43 | 989803160281..... | 13 |
| 40401A..... | 64 | 989803144471..... | 16 | 989803160301..... | 13 |
| 40401B..... | 64 | 989803144531..... | 17 | 989803162591..... | 95 |
| 40401C..... | 64 | 989803144631..... | 100 | 989803162601..... | 95 |
| 40401D..... | 65 | 989803146911..... | 43 | 989803162621..... | 96 |

| | | | | | |
|-------------------|-----|-------------|--------|-------------|--------|
| 989803162631..... | 97 | M1504A..... | 39 | M1624A..... | 32 |
| 989803162641..... | 97 | M1509A..... | 39 | M1625A..... | 37 |
| 989803163271..... | 97 | M1510A..... | 35 | M1626A..... | 33 |
| 989803163281..... | 97 | M1520A..... | 36 | M1627A..... | 84 |
| 989803163331..... | 101 | M1530A..... | 36 | M1629A..... | 37 |
| 989803163341..... | 101 | M1532A..... | 26 | M1631A..... | 37 |
| 989803163351..... | 101 | M1533A..... | 27 | M1633A..... | 37 |
| 989803163361..... | 101 | M1537A..... | 26 | M1634A..... | 53 |
| 989803163371..... | 101 | M1538A..... | 27 | M1635A..... | 37 |
| 989803163381..... | 101 | M1540C..... | 35 | M1636A..... | 39 |
| 989803163631..... | 101 | M1550C..... | 35 | M1638A..... | 39 |
| CPJ840J6..... | 54 | M1557A..... | 26 | M1639A..... | 37 |
| CPJ84022..... | 54 | M1558A..... | 27 | M1642A..... | 22 |
| CPJ84026..... | 54 | M1560C..... | 36 | M1643A..... | 22 |
| CPJ84046..... | 54 | M1567A..... | 53 | M1644A..... | 24, 30 |
| CPJ84447..... | 54 | M1568A..... | 53 | M1645A..... | 25, 31 |
| M1020-61100..... | 91 | M1570C..... | 36 | M1646A..... | 22 |
| M1020-61102..... | 91 | M1571A..... | 55 | M1647A..... | 24, 30 |
| M1131A..... | 86 | M1572A..... | 55 | M1648A..... | 25, 31 |
| M1132A..... | 86 | M1573A..... | 55 | M1649A..... | 34 |
| M1133A..... | 86 | M1574A..... | 56 | M1655A..... | 11 |
| M1134A..... | 87 | M1575A..... | 56 | M1655B..... | 11 |
| M1191B..... | 80 | M1576A..... | 56 | M1656B..... | 12 |
| M1191BL..... | 80 | M1577A..... | 66 | M1657B..... | 12 |
| M1191T..... | 81 | M1578A..... | 66 | M1658A..... | 12 |
| M1192A..... | 81 | M1579A..... | 67 | M1659A..... | 12 |
| M1192T..... | 82 | M1580A..... | 36 | M1662A..... | 12 |
| M1193A..... | 82 | M1590A..... | 36 | M1663A..... | 24 |
| M1193T..... | 83 | M1598B..... | 75 | M1664A..... | 39 |
| M1194A..... | 83 | M1599B..... | 75 | M1665A..... | 26 |
| M1195A..... | 84 | M1600A..... | 36 | M1667A..... | 28 |
| M1196A..... | 84 | M1601A..... | 36 | M1668A..... | 30 |
| M1196T..... | 84 | M1602A..... | 24 | M1669A..... | 32 |
| M1421A..... | 45 | M1603A..... | 36 | M1671A..... | 32 |
| M1422A..... | 45 | M1604A..... | 25 | M1672A..... | 33 |
| M1423A..... | 45 | M1605A..... | 36 | M1673A..... | 32 |
| M1425A..... | 45 | M1608A..... | 36 | M1674A..... | 33 |
| M1430A..... | 45 | M1609A..... | 36 | M1675A..... | 32 |
| M1431A..... | 45 | M1610A..... | 36 | M1678A..... | 33 |
| M1433A..... | 45 | M1611A..... | 37 | M1679A..... | 39 |
| M1465A..... | 20 | M1612A..... | 11, 18 | M1680A..... | 26, 28 |
| M1500A..... | 35 | M1613A..... | 37 | M1681A..... | 27, 29 |
| M1501A..... | 39 | M1615A..... | 37 | M1682A..... | 26, 28 |
| M1502A..... | 39 | M1619A..... | 37 | M1683A..... | 27, 29 |
| M1503A..... | 39 | M1621A..... | 37 | M1684A..... | 26, 28 |
| | | M1622A..... | 32 | M1685A..... | 27, 29 |
| | | M1623A..... | 37 | M1837A..... | 92 |

| | | | | | |
|--------------|--------|--------------|-----|--------------|----|
| M1874A..... | 68 | M2524A..... | 13 | M4554B5..... | 58 |
| M1875A..... | 68 | M2526A..... | 13 | M4555B..... | 58 |
| M1876A..... | 68 | M2533A..... | 19 | M4555B5..... | 58 |
| M1877A..... | 68 | M2536A..... | 19 | M4556B..... | 59 |
| M1878A..... | 69 | M2590A..... | 45 | M4556B5..... | 59 |
| M1879A..... | 69 | M2591A..... | 45 | M4557B..... | 59 |
| M1885A..... | 79 | M2592A..... | 45 | M4557B5..... | 59 |
| M1900B..... | 88 | M2593A..... | 45 | M4558B..... | 60 |
| M1901B..... | 90 | M2594A..... | 45 | M4558B5..... | 60 |
| M1902B..... | 90 | M2595A..... | 45 | M4559B..... | 60 |
| M1903B..... | 90 | M2596A..... | 45 | M4559B5..... | 60 |
| M1904B..... | 90 | M2597A..... | 45 | M4560A..... | 77 |
| M1920A..... | 13 | M2598A..... | 45 | M4561A..... | 77 |
| M1921A..... | 13 | M2599A..... | 45 | M4562B..... | 62 |
| M1923A..... | 13 | M2741A..... | 15 | M4563B..... | 62 |
| M1931A..... | 38 | M2744A..... | 15 | M4564B..... | 62 |
| M1932A..... | 38 | M2745A..... | 15 | M4565B..... | 62 |
| M1934A..... | 38 | M2746A..... | 15 | M4566B..... | 62 |
| M1935A..... | 38, 52 | M2750A..... | 15 | M4567B..... | 62 |
| M1937A..... | 38 | M2751A..... | 15 | M4568B..... | 62 |
| M1940A..... | 88 | M2756A..... | 16 | M4569B..... | 62 |
| M1941A..... | 88 | M2757A..... | 16 | M4570A..... | 77 |
| M1943A..... | 88 | M2760A..... | 16 | M4571A..... | 77 |
| M1943AL..... | 88 | M2761A..... | 16 | M4572B..... | 70 |
| M1943NL..... | 88 | M2768A..... | 17 | M4573B..... | 70 |
| M1949A..... | 24 | M2772A..... | 17 | M4574B..... | 70 |
| M1968A..... | 24, 30 | M2773A..... | 17 | M4575B..... | 70 |
| M1971A..... | 25, 31 | M2776A..... | 17 | M4576B..... | 71 |
| M1973A..... | 24, 30 | M2777A..... | 17 | M4577B..... | 72 |
| M1974A..... | 25, 31 | M2781A..... | 21 | M4578B..... | 72 |
| M1976A..... | 24 | M2782A..... | 21 | M4579B..... | 72 |
| M1978A..... | 25 | M2783A..... | 21 | M4580A..... | 76 |
| M1979A..... | 24 | M2785A..... | 21 | M4581A..... | 76 |
| M1984A..... | 25 | M2786A..... | 21 | M4582B..... | 73 |
| M2255A..... | 94 | M3913A..... | 35 | M4583B..... | 73 |
| M2267A..... | 14 | M3914A..... | 35 | M4584B..... | 73 |
| M2268A..... | 38 | M3915A..... | 35 | M4585B..... | 73 |
| M2271A..... | 53 | M3918A..... | 75 | M4586B..... | 74 |
| M2272A..... | 53 | M3919A..... | 101 | M4587B..... | 74 |
| M2281A..... | 39 | M3920A..... | 101 | M4588B..... | 74 |
| M2501A..... | 19 | M4550A..... | 77 | M4589B..... | 74 |
| M2505A..... | 19 | M4551A..... | 77 | M4590A..... | 78 |
| M2506A..... | 19 | M4552B..... | 57 | M4591A..... | 78 |
| M2513A..... | 19 | M4552B5..... | 57 | M4592A..... | 78 |
| M2516A..... | 19 | M4553B..... | 57 | M4593A..... | 78 |
| M2520A..... | 13 | M4553B5..... | 57 | M4594A..... | 78 |
| M2522A..... | 13 | M4554B..... | 58 | M4595A..... | 78 |

| | |
|-------------------|-----|
| M4596A..... | 78 |
| M4605A | 100 |
| M4606A | 48 |
| M4607-60001 | 100 |
| M4607A..... | 100 |
| M4608A | 23 |
| M4612A..... | 48 |
| M4613A..... | 48 |
| M4620A..... | 23 |
| M4680A | 13 |
| M4681A..... | 13 |
| M4686A | 13 |
| M4687A..... | 13 |
| M4689A..... | 14 |
| M4690A..... | 14 |
| M4691A..... | 14 |
| M4787A..... | 88 |
| M4790A..... | 100 |
| M4793A..... | 46 |
| M4794A..... | 46 |
| M4795A..... | 46 |
| M4796A..... | 46 |
| M4816A..... | 79 |
| M4817A..... | 79 |
| M4821A..... | 99 |
| M4822A | 99 |
| M4823A..... | 98 |
| | |
| PSE11268..... | 79 |
| | |
| ULBU9VLJ..... | 100 |



**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare
healthcare@philips.com

Asia

+49 7031 463 2254

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

Please visit www.philips.com/healthcarestore



© 2010 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.
4522 962 54721 * JAN 2010